

SUPPLIES & ACCESSORIES
US PRODUCT CATALOG

SUPPLIES & ACCESSORIES

NON-INVASIVE BLOOD PRESSURE CUFFS



TRULINK® NYLON REUSABLE BLOOD PRESSURE CUFFS



TruLink Nylon Reusable Blood Pressure Cuffs are available in seven adult/pediatric sizes for accurate measurements.

- Color coded for accurate sizing
- Manufactured without PVC/DEHP
- Provides a fluid barrier
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

TruLink Nylon Reusable BP Cuffs

Single Tube

For use with Command Modules 90496/91496, DM3 Monitor 91330 and mCare Monitor 91220



715-1121-10	Single tube, HP	Child	13 - 20 cm	1 Each
715-1131-10	Single tube, HP	Small Adult	18 - 26 cm	1 Each
715-1141-10	Single tube, HP	Adult	26 - 35 cm	1 Each
715-1151-10	Single tube, HP	Adult Long	29 - 38 cm	1 Each
715-1161-10	Single tube, HP	Large Adult	32 - 42 cm	1 Each
715-1171-10	Single tube, HP	Large Adult Long	35 - 44 cm	1 Each
715-1181-10	Single tube, HP	Thigh Cuff	42 - 50 cm	1 Each

Single Tube Hose



714-0018-00	90496/91496/91330	Adult	275 cm/9 ft	1 Each
714-0060-00	90496/91496/91330	Adult	366 cm/12 ft	1 Each
714-0061-00	91220	Adult	275 cm/9 ft	1 Each

TRULINK® SOFT DISPOSABLE BLOOD PRESSURE CUFFS



TruLink Soft Disposable Blood Pressure Cuffs are available in seven adult/pediatric sizes for accurate measurements.

- Soft fabric helps prevent perspiration under cuff
- Color coded for accurate sizing
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

TruLink Soft Disposable BP Cuffs

Single Tube

For use with Command Modules 90496/91496, DM3 Monitor 91330 and mCare Monitor 91220



Bayonet (HP)

715-1123-10	Single tube, HP	Child	13 - 20 cm	5/Bag
715-1133-10	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
715-1143-10	Single tube, HP	Adult	26 - 35 cm	5/Bag
715-1153-10	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
715-1163-10	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
715-1173-10	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
715-1183-10	Single tube, HP	Thigh Cuff	42 - 50 cm	5/Bag

Single Tube Hose



714-0018-00	90496/91496/91330	Adult	275 cm/9 ft	1 Each
714-0060-00	90496/91496/91330	Adult	366 cm/12 ft	1 Each
714-0061-00	91220	Adult	275 cm/9 ft	1 Each

TRULINK® VINYL DISPOSABLE BLOOD PRESSURE CUFFS



TruLink Vinyl Disposable Blood Pressure Cuffs are available in seven adult/pediatric sizes for accurate measurements.

- Comfortable soft vinyl
- Easy to clean with PDI cloths
- Color coded for accurate sizing
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

TruLink Vinyl Disposable BP Cuffs

Single Tube

For use with Command Modules 90496/91496, DM3 Monitor 91330 and mCare Monitor 91220



Bayonet (HP)

715-1122-10	Single Tube, HP	Child	13 - 20 cm	5/Bag
715-1132-10	Single Tube, HP	Small Adult	18 - 26 cm	5/Bag
715-1142-10	Single Tube, HP	Adult	26 - 35 cm	5/Bag
715-1152-10	Single Tube, HP	Adult Long	29 - 38 cm	5/Bag
715-1162-10	Single Tube, HP	Large Adult	32 - 42 cm	5/Bag
715-1172-10	Single Tube, HP	Large Adult Long	35 - 44 cm	5/Bag
715-1182-10	Single Tube, HP	Thigh Cuff	42 - 50 cm	5/Bag

Single Tube Hose



714-0018-00	90496/91496/91330	Adult	275 cm/9 ft	1 Each
714-0060-00	90496/91496/91330	Adult	366 cm/12 ft	1 Each
714-0061-00	91220	Adult	275 cm/9 ft	1 Each

TRULINK® YELLOW VINYL DISPOSABLE BLOOD PRESSURE CUFFS



TruLink Yellow Vinyl Disposable Blood Pressure Cuffs intended for use in isolation areas are available in three adult/pediatric sizes.

- Color coded for accurate sizing
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

TruLink Yellow Vinyl Disposable BP Cuffs

Single Tube

For use with Command Modules 90496/91496, DM3 Monitor 91330 and mCare Monitor 91220



Bayonet (HP)

715-1135-10	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
715-1145-10	Single tube, HP	Adult	26 - 35 cm	5/Bag
715-1165-10	Single tube, HP	Large Adult	32 - 42 cm	5/Bag

Single Tube Hose



714-0018-00	90496/91496/91330	Adult	275 cm/9 ft	1 Each
714-0060-00	90496/91496/91330	Adult	366 cm/12 ft	1 Each
714-0061-00	91220	Adult	275 cm/9 ft	1 Each

TRULINK® NEONATAL NON-PVC BLOOD PRESSURE CUFFS



TruLink Neonatal Non-PVC Blood Pressure Cuffs are available in five neonatal sizes for accurate measurements.

- Color coded for accurate sizing
- Genuine low profile Velcro® closures
- Manufactured without natural latex
- Fabric fold is more comfortable than a weld



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

TruLink Neonatal Non-PVC BP Cuffs

Single Tube

For use with Command Modules 90496/91496 and mCare Monitor 91220



Mini Male Luer
(ML316)

715-2114-10	Single tube, 3/16" ML	Size 1	3-6 cm	10/Box
715-2124-10	Single tube, 3/16" ML	Size 2	4-8 cm	10/Box
715-2134-10	Single tube, 3/16" ML	Size 3	6-11 cm	10/Box
715-2144-10	Single tube, 3/16" ML	Size 4	7-13 cm	10/Box
715-2154-10	Single tube, 3/16" ML	Size 5	8-15 cm	10/Box

Single Tube Hose



714-0019-01	Single tube, 90496/91496	Neonatal	275 cm/9 ft	1 Each
714-0062-00	Single tube, 91220	Neonatal	275 cm/9 ft	1 Each

SPACELABS NYLON REUSABLE BLOOD PRESSURE CUFFS



Spacelabs brand cuffs provide a comprehensive solution for your NIBP needs. The Nylon Reusable Blood Pressure Cuffs are available in both regular and long sizes.

- Long size cuffs for children, small adults, adults and large adults
- Long cuffs feature additional length allowing ¼ extra turn for greater security and elimination of “pop-offs”



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

Spacelabs Nylon Reusable BP Cuffs

Single Tube

For use with Command Modules 90496/91496, DM3 Monitor 91330 and mCare Monitor 91220



Metal Bayonet (MHP)

714-0020-00	Single tube, MHP	Child	12 - 19 cm	1 Each
714-0070-00	Single tube, MHP	Child Long	12 - 19 cm	1 Each
714-0021-00	Single tube, MHP	Small Adult	17 - 25 cm	1 Each
714-0071-00	Single tube, MHP	Small Adult Long	17 - 25 cm	1 Each
714-0023-00	Single tube, MHP	Adult	23 - 33 cm	1 Each
714-0073-00	Single tube, MHP	Adult Long	23 - 33 cm	1 Each
714-0025-00	Single tube, MHP	Large Adult	31 - 40 cm	1 Each
714-0075-00	Single tube, MHP	Large Adult Long	31 - 40 cm	1 Each
714-0027-00	Single tube, MHP	Thigh	38 - 50 cm	1 Each

Single Tube Hose



714-0018-00	90496/91496/91330	Adult	275 cm/9 ft	1 Each
714-0060-00	90496/91496/91330	Adult	366 cm/12 ft	1 Each
714-0061-00	91220	Adult	275 cm/9 ft	1 Each

SPACELABS VINYL DISPOSABLE BLOOD PRESSURE CUFFS



Spacelabs brand cuffs provide a comprehensive solution for your NIBP needs. The Vinyl Disposable Blood Pressure Cuffs are available in seven sizes for accurate measurement.

- Long size cuffs for adults and large adults
- Long cuffs feature additional length allowing ¼ extra turn for greater security and elimination of “pop-offs”



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

Spacelabs Vinyl Disposable BP Cuffs

Single Tube

For use with Command Modules 90496/91496, DM3 Monitor 91330 and mCare Monitor 91220



Bayonet (HP)

714-1020-00	Single tube, HP	Child	12 - 19 cm	5/Pack
714-1021-00	Single tube, HP	Small Adult	17 - 25 cm	5/Pack
714-1023-00	Single tube, HP	Adult	23 - 33 cm	5/Pack
714-1073-00	Single tube, HP	Adult Long	23 - 33 cm	5/Pack
714-1025-00	Single tube, HP	Large Adult	31 - 40 cm	5/Pack
714-1075-00	Single tube, HP	Large Adult Long	31 - 40 cm	5/Pack
714-1027-00	Single tube, HP	Thigh	38 - 50 cm	5/Pack

Single Tube Hose



714-0018-00	90496/91496/91330	Adult	275 cm/9 ft	1 Each
714-0060-00	90496/91496/91330	Adult	366 cm/12 ft	1 Each
714-0061-00	91220	Adult	275 cm/9 ft	1 Each

SPACELABS SOFT CUFF NEONATAL DISPOSABLE BLOOD PRESSURE CUFFS



Spacelabs brand cuffs provide a comprehensive solution for your NIBP needs.



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

Spacelabs Soft Cuff Neonatal BP Cuffs

Single Tube

For use with Command Modules 90496/91496 and mCare Monitor 91220



Mini Male Luer
(ML316)

714-1028-01	Single tube, 3/16" ML	Size 1	3 - 6 cm	10/Pack
714-1029-01	Single tube, 3/16" ML	Size 2	4 - 8 cm	10/Pack
714-1030-01	Single tube, 3/16" ML	Size 3	6 - 11 cm	10/Pack
714-1032-01	Single tube, 3/16" ML	Size 4	7 - 13 cm	10/Pack
714-1033-01	Single tube, 3/16" ML	Size 5	8 - 15 cm	10/Pack

Single Tube Hose



714-0019-01	Single tube, 90496/91496	Neonatal	275 cm/9 ft	1 Each
714-0062-00	Single tube hose, 91220	Neonatal	275 cm/9 ft	1 Each

NIBP accessories for use with Spacelabs and TruLink branded cuffs.

- Hose
- Fittings/Adapters
- Conversion Accessories



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

Non-Invasive Blood Pressure Hose

Single Tube Hose

For use with Spacelabs Command Modules 90496/91496, DM3 Monitor 91330 and mCare Monitor 91220



714-0018-00	90496/91496/91330	Adult	275 cm/9 ft	1 Each
714-0060-00	90496/91496/91330	Adult	366 cm/12 ft	1 Each
714-0061-00	91220	Adult	275 cm/9 ft	1 Each

Double Tube to Single Tube Hose

For use with Command Modules 90496/91496



714-0018-12	90496/91496	Adult	275 cm/9 ft	1 Each
-------------	-------------	-------	-------------	--------

Hose Connectors

Slip Fit Connector

Connectors fit Spacelabs standard NIBP tubing. Adapts hoses for use with quick connect fittings.



714-0038-00	3mm ID Adult NIBP hoses			1 Each
714-0047-00	4mm ID Adult NIBP hoses			1 Each
714-0046-00	5mm ID Adult NIBP hoses			1 Each
714-0045-00	6mm ID Adult NIBP hoses			1 Each

NIBP accessories for use with Spacelabs and TruLink branded cuffs.

- Hose
- Fittings/Adapters
- Conversion Accessories



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

NIBP Cuff Conversion Accessories

Hose Conversion Kit



714-0043-00	Dual hose conversion kit Use with TruLink and Spacelabs cuffs and dual hose blood pressure monitors incorporating thread type fittings	Adult		1 Each
-------------	---	-------	--	--------

Y-Adapters



714-0056-00	NIBP dual to single hose Y-adapter Connects dual-tube NIBP cuffs to single-tube 3mm ID hose and monitor/modules	Adult		5/Pack
-------------	--	-------	--	--------



714-0057-00	NIBP dual to single hose Y-adapter Connects dual-tube NIBP cuffs to single tube NIBP hoses and module/monitors	Adult		5/Pack
-------------	---	-------	--	--------

NIBP accessories for use with Spacelabs and TruLink branded neonatal cuffs.

- Hose
- Fittings/Adapters
- Conversion Accessories



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

Neonatal BP Cuff Accessories

Single Tube Hose

For use with Command Modules 90496/91496 and mCare Monitor 91220



714-0019-01	90496/91496	Neonatal	275 cm/9 ft	1 Each
714-0062-00	91220	Neonatal	275 cm/9 ft	1 Each

Double Tube Cuff to Single Tube Hose

For use with Command Modules 90496/91496



714-0019-12	90496/91496	Neonatal	275 cm/9 ft	1 Each
-------------	-------------	----------	-------------	--------

Metal Cuff Fitting

For use with NIBP hose 714-0019-01



714-0053-10	Hose fitting, metal	Neonatal		10/Pack
-------------	---------------------	----------	--	---------

Plastic Cuff Fittings

For use with NIBP Hoses 714-0019-01 & 714-0019-12



714-0052-50	Cuff fitting, plastic	Neonatal		50/Pack
714-0054-50	Cuff fitting, plastic	Infant		50/Pack

Hose Connectors



714-0057-10	NIBP non-luer hose metal connector Connects neonatal cuffs to non-Spacelabs monitors with a 3mm ID hose	Neonatal		10/Pack
-------------	--	----------	--	---------

ULTRACHECK® NYLON REUSABLE BLOOD PRESSURE CUFFS



UltraCheck Nylon Reusable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Fluid barrier
- Meets AHA / AAMI standards
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

UltraCheck Nylon Reusable BP Cuffs

Single Tube				
US0814HP	Single tube, HP	Infant	8 - 14 cm	Each
US0814HP-05	Single tube, HP	Infant	8 - 14 cm	5/Bag
US1320HP	Single tube, HP	Child	13 - 20 cm	Each
US1320HP-05	Single tube, HP	Child	13 - 20 cm	5/Bag
US1826HP	Single tube, HP	Small Adult	18 - 26 cm	Each
US1826HP-05	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
US2635HP	Single tube, HP	Adult	26 - 35 cm	Each
US2635HP-05	Single tube, HP	Adult	26 - 35 cm	5/Bag
US2938HP	Single tube, HP	Adult Long	29 - 38 cm	Each
US2938HP-05	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
US3242HP	Single tube, HP	Large Adult	32 - 42 cm	Each
US3242HP-05	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
US3544HP	Single tube, HP	Large Adult Long	35 - 44 cm	Each
US3544HP-05	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
US4250HP	Single tube, HP	Thigh	42 - 50 cm	Each
US4250HP-05	Single tube, HP	Thigh	42 - 50 cm	5/Bag



ULTRACHECK® NYLON REUSABLE BLOOD PRESSURE CUFFS



UltraCheck Nylon Reusable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Fluid barrier
- Meets AHA / AAMI standards
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

UltraCheck Nylon Reusable BP Cuffs

Single Tube



Metal Bayonet (MHP)

US0814MHP-05	Single tube, MHP	Infant	8 - 14 cm	5/Bag
US1320MHP-05	Single tube, MHP	Child	13 - 20 cm	5/Bag
US1826MHP-05	Single tube, MHP	Small Adult	18 - 26 cm	5/Bag
US2635MHP-05	Single tube, MHP	Adult	26 - 35 cm	5/Bag
US2938MHP-05	Single tube, MHP	Adult Long	29 - 38 cm	5/Bag
US3242MHP-05	Single tube, MHP	Large Adult	32 - 42 cm	5/Bag
US3544MHP-05	Single tube, MHP	Large Adult Long	35 - 44 cm	5/Bag
US4250MHP-05	Single tube, MHP	Thigh	42 - 50 cm	5/Bag



Male Luer (ML)

US0814ML-05	Single tube, ML	Infant	8 - 14 cm	5/Bag
US1320ML-05	Single tube, ML	Child	13 - 20 cm	5/Bag
US1826ML-05	Single tube, ML	Small Adult	18 - 26 cm	5/Bag
US2635ML-05	Single tube, ML	Adult	26 - 35 cm	5/Bag
US2938ML-05	Single tube, ML	Adult Long	29 - 38 cm	5/Bag
US3242ML-05	Single tube, ML	Large Adult	32 - 42 cm	5/Bag
US3544ML-05	Single tube, ML	Large Adult Long	35 - 44 cm	5/Bag
US4250ML-05	Single tube, ML	Thigh	42 - 50 cm	5/Bag



Screw (DM)

US0814DM-05	Single tube, DM	Infant	8 - 14 cm	5/Bag
US1320DM-05	Single tube, DM	Child	13 - 20 cm	5/Bag
US1826DM-05	Single tube, DM	Small Adult	18 - 26 cm	5/Bag
US2635DM-05	Single tube, DM	Adult	26 - 35 cm	5/Bag
US2938DM-05	Single tube, DM	Adult Long	29 - 38 cm	5/Bag
US3544DM-05	Single tube, DM	Large Adult	32 - 42 cm	5/Bag
US3242DM-05	Single tube, DM	Large Adult	32 - 42 cm	5/Bag
US4250DM-05	Single tube, DM	Thigh	42 - 50 cm	5/Bag

ULTRACHECK® NYLON REUSABLE BLOOD PRESSURE CUFFS



UltraCheck Nylon Reusable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Fluid barrier
- Meets AHA / AAMI standards
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

UltraCheck Nylon Reusable BP Cuffs

Single Tube



Twist Lock (MQ)

US0814MQ-05	Single tube, MQ	Infant	8 - 14 cm	5/Bag
US1320MQ-05	Single tube, MQ	Child	13 - 20 cm	5/Bag
US1826MQ-05	Single tube, MQ	Small Adult	18 - 26 cm	5/Bag
US2635MQ-05	Single tube, MQ	Adult	26 - 35 cm	5/Bag
US2938MQ-05	Single tube, MQ	Adult Long	29 - 38 cm	5/Bag
US3242MQ-05	Single tube, MQ	Large Adult	32 - 42 cm	5/Bag
US3544MQ-05	Single tube, MQ	Large Adult Long	35 - 44 cm	5/Bag
US4250MQ-05	Single tube, MQ	Thigh	42 - 50 cm	5/Bag



Female Luer (FL)

US0814FL-05	Single tube, FL	Infant	8 - 14 cm	5/Bag
US1320FL-05	Single tube, FL	Child	13 - 20 cm	5/Bag
US1826FL-05	Single tube, FL	Small Adult	18 - 26 cm	5/Bag
US2635FL-05	Single tube, FL	Adult	26 - 35 cm	5/Bag
US2938FL-05	Single tube, FL	Adult Long	29 - 38 cm	5/Bag
US3242FL-05	Single tube, FL	Large Adult	32 - 42 cm	5/Bag
US3544FL-05	Single tube, FL	Large Adult Long	35 - 44 cm	5/Bag
US4250FL-05	Single tube, FL	Thigh	42 - 50 cm	5/Bag



Male Twist Lock (MQML)

US0814MQML-05	Single tube, MQML	Infant	8 - 14 cm	5/Bag
US1320MQML-05	Single tube, MQML	Child	13 - 20 cm	5/Bag
US1826MQML-05	Single tube, MQML	Small Adult	18 - 26 cm	5/Bag
US2635MQML-05	Single tube, MQML	Adult	26 - 35 cm	5/Bag
US2938MQML-05	Single tube, MQML	Adult Long	29 - 38 cm	5/Bag
US3242MQML-05	Single tube, MQML	Large Adult	32 - 42 cm	5/Bag
US3544MQML-05	Single tube, MQML	Large Adult Long	35 - 44 cm	5/Bag
US4250MQML-05	Single tube, MQML	Thigh	42 - 50 cm	5/Bag

ULTRACHECK® NYLON REUSABLE BLOOD PRESSURE CUFFS



UltraCheck Nylon Reusable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Fluid barrier
- Meets AHA / AAMI standards
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

UltraCheck Nylon Reusable BP Cuffs

Single Tube



US0814-05	Single tube, no fitting	Infant	8 - 14 cm	5/Bag
US1320-05	Single tube, no fitting	Child	13 - 20 cm	5/Bag
US1826-05	Single tube, no fitting	Small Adult	18 - 26 cm	5/Bag
US2635-05	Single tube, no fitting	Adult	26 - 35 cm	5/Bag
US2938-05	Single tube, no fitting	Adult Long	29 - 38 cm	5/Bag
US3242-05	Single tube, no fitting	Large Adult	32 - 42 cm	5/Bag
US3544-05	Single tube, no fitting	Large Adult Long	35 - 44 cm	5/Bag
US4250-05	Single tube, no fitting	Thigh	42 - 50 cm	5/Bag



US0814MB	Single tube, MB w/gauge	Infant	8 - 14 cm	1 Each
US1320MB	Single tube, MB w/gauge	Child	13 - 20 cm	1 Each
US1826MB	Single tube, MB w/gauge	Small Adult	18 - 26 cm	1 Each
US2635MB	Single tube, MB w/gauge	Adult	26 - 35 cm	1 Each
US2938MB	Single tube, MB w/gauge	Adult Long	29 - 38 cm	1 Each
US3242MB	Single tube, MB w/gauge	Large Adult	32 - 42 cm	1 Each
US3544MB	Single tube, MB w/gauge	Large Adult Long	35 - 44 cm	1 Each
US4250MB	Single tube, MB w/gauge	Thigh	42 - 50 cm	1 Each

ULTRACHECK® NYLON REUSABLE BLOOD PRESSURE CUFFS



UltraCheck Nylon Reusable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Fluid barrier
- Meets AHA / AAMI standards
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

UltraCheck Nylon Reusable BP Cuffs

Double Tube



Screw (DM)

UD0814DM-05	Double tube, DM	Infant	8 - 14 cm	5/Bag
UD1320DM-05	Double tube, DM	Child	13 - 20 cm	5/Bag
UD1826DM-05	Double tube, DM	Small Adult	18 - 26 cm	5/Bag
UD2635DM-05	Double tube, DM	Adult	26 - 35 cm	5/Bag
UD2938DM-05	Double tube, DM	Adult Long	29 - 38 cm	5/Bag
UD3242DM-05	Double tube, DM	Large Adult	32 - 42 cm	5/Bag
UD3544DM-05	Double tube, DM	Large Adult Long	35 - 44 cm	5/Bag
UD4250DM-05	Double tube, DM	Thigh	42 - 50 cm	5/Bag



Twist Lock (MQ)

UD0814MQ-05	Double tube, MQ	Infant	8 - 14 cm	5/Bag
UD1320MQ-05	Double tube, MQ	Child	13 - 20 cm	5/Bag
UD1826MQ-05	Double tube, MQ	Small Adult	18 - 26 cm	5/Bag
UD2635MQ-05	Double tube, MQ	Adult	26 - 35 cm	5/Bag
UD2938MQ-05	Double tube, MQ	Adult Long	29 - 38 cm	5/Bag
UD3242MQ-05	Double tube, MQ	Large Adult	32 - 42 cm	5/Bag
UD3544MQ-05	Double tube, MQ	Large Adult Long	35 - 44 cm	5/Bag
UD4250MQ-05	Double tube, MQ	Thigh	42 - 50 cm	5/Bag



Male & Female MQ (MQMF)

UD0814MQMF-05	Double tube, MQMF	Infant	8 - 14 cm	5/Bag
UD1320MQMF-05	Double tube, MQMF	Child	13 - 20 cm	5/Bag
UD1826MQMF-05	Double tube, MQMF	Small Adult	18 - 26 cm	5/Bag
UD2635MQMF-05	Double tube, MQMF	Adult	26 - 35 cm	5/Bag
UD2938MQMF-05	Double tube, MQMF	Adult Long	29 - 38 cm	5/Bag
UD3242MQMF-05	Double tube, MQMF	Large Adult	32 - 42 cm	5/Bag
UD3544MQMF-05	Double tube, MQMF	Large Adult Long	35 - 44 cm	5/Bag
UD4250MQMF-05	Double tube, MQMF	Thigh	42 - 50 cm	5/Bag

ULTRACHECK® NYLON REUSABLE BLOOD PRESSURE CUFFS



UltraCheck Nylon Reusable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Fluid barrier
- Meets AHA / AAMI standards
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

UltraCheck Nylon Reusable BP Cuffs

Double Tube



UD0814-05	Double tube, no fitting	Infant	8 - 14 cm	5/Bag
UD1320-05	Double tube, no fitting	Child	13 - 20 cm	5/Bag
UD1826-05	Double tube, no fitting	Small Adult	18 - 26 cm	5/Bag
UD2635-05	Double tube, no fitting	Adult	26 - 35 cm	5/Bag
UD2938-05	Double tube, no fitting	Adult Long	29 - 38 cm	5/Bag
UD3242-05	Double tube, no fitting	Large Adult	32 - 42 cm	5/Bag
UD3544-05	Double tube, no fitting	Large Adult Long	35 - 44 cm	5/Bag
UD4250-05	Double tube, no fitting	Thigh	42 - 50 cm	5/Bag

Double Tube



UD0814MB-05	Double tube, MB	Infant	8 - 14 cm	5/Bag
UD1320MB-05	Double tube, MB	Child	13 - 20 cm	5/Bag
UD1826MB-05	Double tube, MB	Small Adult	18 - 26 cm	5/Bag
UD2635MB-05	Double tube, MB	Adult	26 - 35 cm	5/Bag
UD2938MB-05	Double tube, MB	Adult Long	29 - 38 cm	5/Bag
UD3242MB-05	Double tube, MB	Large Adult	32 - 42 cm	5/Bag
UD3544MB-05	Double tube, MB	Large Adult Long	35 - 44 cm	5/Bag
UD4250MB-05	Double tube, MB	Thigh	42 - 50 cm	5/Bag

ULTRACHECK® CURVE™ REUSABLE BLOOD PRESSURE CUFFS



Worldwide obesity rates are climbing. This disease is strongly associated with increased mortality and morbidity, cardiovascular disease, diabetes, and many other illnesses. Blood pressure is the number one vital sign in detecting these issues and accurate blood pressure is critical. The Curve Cuffs are designed to fit a bariatric patient's arm correctly. Their conical shape fits snugly against the arm and properly compresses the artery.

- Manufactured without DEHP
- Will not slide during inflation
- Genuine Velcro® closures
- Manufactured without natural latex
- Adapters allow for compatibility with monitors worldwide



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

UltraCheck Curve Reusable BP Cuffs

Single Tube



BRUS3854	Single tube, no fitting	Extra Large Adult	38 - 54 cm	1 Each
BRUS3854DM	Single tube, DM	Extra Large Adult	38 - 54 cm	1 Each
BRUS3854FL	Single tube, FL	Extra Large Adult	38 - 54 cm	1 Each
BRUS3854HP	Single tube, HP	Extra Large Adult	38 - 54 cm	1 Each
BRUS3854MHP	Single tube, metal HP	Extra Large Adult	38 - 54 cm	1 Each
BRUS3854ML	Single tube, ML	Extra Large Adult	38 - 54 cm	1 Each
BRUS3854MQ	Single tube, MQ	Extra Large Adult	38 - 54 cm	1 Each

Double Tube



BRUD3854	Double tube, no fitting	Extra Large Adult	38 - 54 cm	1 Each
BRUD3854DM	Double tube, DM	Extra Large Adult	38 - 54 cm	1 Each
BRUD3854HP	Double tube, HP	Extra Large Adult	38 - 54 cm	1 Each
BRUD3854ML	Double tube, ML	Extra Large Adult	38 - 54 cm	1 Each
BRUD3854MQ	Double tube, MQ	Extra Large Adult	38 - 54 cm	1 Each
BRUD3854MQMF	Double tube, MQMF	Extra Large Adult	38 - 54 cm	1 Each
BRUD3854MQML	Double tube, MQML	Extra Large Adult	38 - 54 cm	1 Each

Manual Cuffs



BRUS3854MB	Single tube, manual w/ gauge	Extra Large Adult	38 - 54 cm	1 Each
------------	------------------------------	-------------------	------------	--------



BRUD3854MB	Double tube, manual	Extra Large Adult	38 - 54 cm	1 Each
------------	---------------------	-------------------	------------	--------

SOFTCHECK® FABRIC DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Fabric Disposable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Soft fabric helps prevent perspiration under cuff
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Fabric Disposable BP Cuffs

Single Tube



Bayonet (HP)

ST0814HP-05	Single tube, HP	Infant	8 - 14 cm	5/Bag
ST0814HP-20	Single tube, HP	Infant	8 - 14 cm	20/Box
ST1320HP-05	Single tube, HP	Child	13 - 20 cm	5/Bag
ST1320HP-20	Single tube, HP	Child	13 - 20 cm	20/Box
ST1826HP-05	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
ST1826HP-20	Single tube, HP	Small Adult	18 - 26 cm	20/Box
ST2635HP-05	Single tube, HP	Adult	26 - 35 cm	5/Bag
ST2635HP-20	Single tube, HP	Adult	26 - 35 cm	20/Box
ST2938HP-05	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
ST2938HP-20	Single tube, HP	Adult Long	29 - 38 cm	20/Box
ST3242HP-05	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
ST3242HP-20	Single tube, HP	Large Adult	32 - 42 cm	20/Box
ST3544HP-05	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
ST3544HP-20	Single tube, HP	Large Adult Long	35 - 44 cm	20/Box
ST4250HP-05	Single tube, HP	Thigh	42 - 50 cm	5/Bag
ST4250HP-20	Single tube, HP	Thigh	42 - 50 cm	20/Box



Male Luer (ML)

ST2635ML-05	Single tube, ML	Adult	26 - 35 cm	5/Bag
ST2938ML-05	Single tube, ML	Adult Long	29 - 38 cm	5/Bag
ST2938ML-20	Single tube, ML	Adult Long	29 - 38 cm	20/Box
ST3242ML-05	Single tube, ML	Large Adult	32 - 42 cm	5/Bag
ST3242ML-20	Single tube, ML	Large Adult	32 - 42 cm	20/Box
ST3544ML-05	Single tube, ML	Large Adult Long	35 - 44 cm	5/Bag
ST3544ML-20	Single tube, ML	Large Adult Long	35 - 44 cm	20/Box
ST4250ML-05	Single tube, ML	Thigh	42 - 50 cm	5/Bag
ST4250ML-20	Single tube, ML	Thigh	42 - 50 cm	20/Box

SOFTCHECK® FABRIC DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Fabric Disposable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Soft fabric helps prevent perspiration under cuff
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Fabric Disposable BP Cuffs

Single Tube



Twist Lock (MQ)

ST0814MQ-05	Single tube, MQ	Infant	8 - 14 cm	5/Bag
ST0814MQ-20	Single tube, MQ	Infant	8 - 14 cm	20/Box
ST1320MQ-05	Single tube, MQ	Child	13 - 20 cm	5/Bag
ST1320MQ-20	Single tube, MQ	Child	13 - 20 cm	20/Box
ST1826MQ-05	Single tube, MQ	Small Adult	18 - 26 cm	5/Bag
ST1826MQ-20	Single tube, MQ	Small Adult	18 - 26 cm	20/Box
ST2635MQ-05	Single tube, MQ	Adult	26 - 35 cm	5/Bag
ST2635MQ-20	Single tube, MQ	Adult	26 - 35 cm	20/Box
ST2938MQ-05	Single tube, MQ	Adult Long	29 - 38 cm	5/Bag
ST2938MQ-20	Single tube, MQ	Adult Long	29 - 38 cm	20/Box
ST3242MQ-05	Single tube, MQ	Large Adult	32 - 42 cm	5/Bag
ST3242MQ-20	Single tube, MQ	Large Adult	32 - 42 cm	20/Box
ST3544MQ-05	Single tube, MQ	Large Adult Long	35 - 44 cm	5/Bag
ST3544MQ-20	Single tube, MQ	Large Adult Long	35 - 44 cm	20/Box
ST4250MQ-05	Single tube, MQ	Thigh	42 - 50 cm	5/Bag
ST4250MQ-20	Single tube, MQ	Thigh	42 - 50 cm	20/Box



Screw (DM)

ST0814DM-05	Single tube, DM	Infant	8 - 14 cm	5/Bag
ST1320DM-05	Single tube, DM	Child	13 - 20 cm	5/Bag
ST1826DM-05	Single tube, DM	Small Adult	18 - 26 cm	5/Bag
ST2635DM-05	Single tube, DM	Adult	26 - 35 cm	5/Bag
ST2938DM-05	Single tube, DM	Adult Long	29 - 38 cm	5/Bag
ST3242DM-05	Single tube, DM	Large Adult	32 - 42 cm	5/Bag
ST3544DM-05	Single tube, DM	Large Adult Long	35 - 44 cm	5/Bag
ST4250DM-05	Single tube, DM	Thigh	42 - 50 cm	5/Bag

SOFTCHECK® FABRIC DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Fabric Disposable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Soft fabric helps prevent perspiration under cuff
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Fabric Disposable BP Cuffs

Single Tube



Female Luer (FL)

ST0814FL-05	Single tube, FL	Infant	8 - 14 cm	5/Bag
ST1320FL-05	Single tube, FL	Child	13 - 20 cm	5/Bag
ST1826FL-05	Single tube, FL	Small Adult	18 - 26 cm	5/Bag
ST2635FL-05	Single tube, FL	Adult	26 - 35 cm	5/Bag
ST2938FL-05	Single tube, FL	Adult Long	29 - 38 cm	5/Bag
ST3242FL-05	Single tube, FL	Large Adult	32 - 42 cm	5/Bag
ST3544FL-05	Single tube, FL	Large Adult Long	35 - 44 cm	5/Bag
ST4250FL-05	Single tube, FL	Thigh	42 - 50 cm	5/Bag



Male Twist Lock (MQML)

ST0814MQML-05	Single tube, MQML	Infant	8 - 14 cm	5/Bag
ST1320MQML-05	Single tube, MQML	Child	13 - 20 cm	5/Bag
ST1826MQML-05	Single tube, MQML	Small Adult	18 - 26 cm	5/Bag
ST2635MQML-05	Single tube, MQML	Adult	26 - 35 cm	5/Bag
ST2938MQML-05	Single tube, MQML	Adult Long	29 - 38 cm	5/Bag
ST3242MQML-05	Single tube, MQML	Large Adult	32 - 42 cm	5/Bag
ST3544MQML-05	Single tube, MQML	Large Adult Long	35 - 44 cm	5/Bag
ST4250MQML-05	Single tube, MQML	Thigh	42 - 50 cm	5/Bag



No Fitting

ST0814-05	Single tube, no fitting	Infant	8 - 14 cm	5/Bag
ST1320-05	Single tube, no fitting	Child	13 - 20 cm	5/Bag
ST1826-05	Single tube, no fitting	Small Adult	18 - 26 cm	5/Bag
ST2635-05	Single tube, no fitting	Adult	26 - 35 cm	5/Bag
ST2938-05	Single tube, no fitting	Adult Long	29 - 38 cm	5/Bag
ST3242-05	Single tube, no fitting	Large Adult	32 - 42 cm	5/Bag
ST3544-05	Single tube, no fitting	Large Adult Long	35 - 44 cm	5/Bag
ST4250-05	Single tube, no fitting	Thigh	42 - 50 cm	5/Bag

SOFTCHECK® FABRIC DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Fabric Disposable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Soft fabric helps prevent perspiration under cuff
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Fabric Disposable BP Cuffs



Screw (DM)

Double Tube					
DT0814DM-05	Double tube, DM	Infant	8 - 14 cm	5/Bag	
DT0814DM-20	Double tube, DM	Infant	8 - 14 cm	20/Box	
DT1320DM-05	Double tube, DM	Child	13 - 20 cm	5/Bag	
DT1320DM-20	Double tube, DM	Child	13 - 20 cm	20/Box	
DT1826DM-05	Double tube, DM	Small Adult	18 - 26 cm	5/Bag	
DT1826DM-20	Double tube, DM	Small Adult	18 - 26 cm	20/Box	
DT2635DM-05	Double tube, DM	Adult	26 - 35 cm	5/Bag	
DT2635DM-20	Double tube, DM	Adult	26 - 35 cm	20/Box	
DT2938DM-05	Double tube, DM	Adult Long	29 - 38 cm	5/Bag	
DT2938DM-20	Double tube, DM	Adult Long	29 - 38 cm	20/Box	
DT3242DM-05	Double tube, DM	Large Adult	32 - 42 cm	5/Bag	
DT3242DM-20	Double tube, DM	Large Adult	32 - 42 cm	20/Box	
DT3544DM-05	Double tube, DM	Large Adult Long	35 - 44 cm	5/Bag	
DT3544DM-20	Double tube, DM	Large Adult Long	35 - 44 cm	20/Box	
DT4250DM-05	Double tube, DM	Thigh	42 - 50 cm	5/Bag	
DT4250DM-20	Double tube, DM	Thigh	42 - 50 cm	20/Box	

SOFTCHECK® FABRIC DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Fabric Disposable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Soft fabric helps prevent perspiration under cuff
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Fabric Disposable BP Cuffs



Twist Lock (MQ)

Double Tube					
DT0814MQ-05	Double tube, MQ	Infant	8 - 14 cm	5/Bag	
DT0814MQ-20	Double tube, MQ	Infant	8 - 14 cm	20/Box	
DT1320MQ-05	Double tube, MQ	Child	13 - 20 cm	5/Bag	
DT1320MQ-20	Double tube, MQ	Child	13 - 20 cm	20/Box	
DT1826MQ-05	Double tube, MQ	Small Adult	18 - 26 cm	5/Bag	
DT1826MQ-20	Double tube, MQ	Small Adult	18 - 26 cm	20/Box	
DT2635MQ-05	Double tube, MQ	Adult	26 - 35 cm	5/Bag	
DT2635MQ-20	Double tube, MQ	Adult	26 - 35 cm	20/Box	
DT2938MQ-05	Double tube, MQ	Adult Long	29 - 38 cm	5/Bag	
DT2938MQ-20	Double tube, MQ	Adult Long	29 - 38 cm	20/Box	
DT3242MQ-05	Double tube, MQ	Large Adult	32 - 42 cm	5/Bag	
DT3242MQ-20	Double tube, MQ	Large Adult	32 - 42 cm	20/Box	
DT3544MQ-05	Double tube, MQ	Large Adult Long	35 - 44 cm	5/Bag	
DT3544MQ-20	Double tube, MQ	Large Adult Long	35 - 44 cm	20/Box	
DT4250MQ-05	Double tube, MQ	Thigh	42 - 50 cm	5/Bag	
DT4250MQ-20	Double tube, MQ	Thigh	42 - 50 cm	20/Box	

SOFTCHECK® FABRIC DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Fabric Disposable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Soft fabric helps prevent perspiration under cuff
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Fabric Disposable BP Cuffs

Double Tube



DT0814MQMF-05	Double tube, MQMF	Infant	8 - 14 cm	5/Bag
DT1320MQMF-05	Double tube, MQMF	Child	13 - 20 cm	5/Bag
DT1826MQMF-05	Double tube, MQMF	Small Adult	18 - 26 cm	5/Bag
DT2635MQMF-05	Double tube, MQMF	Adult	26 - 35 cm	5/Bag
DT2938MQMF-05	Double tube, MQMF	Adult Long	29 - 38 cm	5/Bag
DT3242MQMF-05	Double tube, MQMF	Large Adult	32 - 42 cm	5/Bag
DT3544MQMF-05	Double tube, MQMF	Large Adult Long	35 - 44 cm	5/Bag
DT4250MQMF-05	Double tube, MQMF	Thigh	42 - 50 cm	5/Bag



DT0814-05	Double tube, no fitting	Infant	8 - 14 cm	5/Bag
DT1320-05	Double tube, no fitting	Child	13 - 20 cm	5/Bag
DT1826-05	Double tube, no fitting	Small Adult	18 - 26 cm	5/Bag
DT2635-05	Double tube, no fitting	Adult	26 - 35 cm	5/Bag
DT2938-05	Double tube, no fitting	Adult Long	29 - 38 cm	5/Bag
DT3242-05	Double tube, no fitting	Large Adult	32 - 42 cm	5/Bag
DT3544-05	Double tube, no fitting	Large Adult Long	35 - 44 cm	5/Bag
DT4250-05	Double tube, no fitting	Thigh	42 - 50 cm	5/Bag



DT0814MB-05	Double tube, MB	Infant	8 - 14 cm	5/Bag
DT1320MB-05	Double tube, MB	Child	13 - 20 cm	5/Bag
DT1826MB-05	Double tube, MB	Small Adult	18 - 26 cm	5/Bag
DT2635MB-05	Double tube, MB	Adult	26 - 35 cm	5/Bag
DT2938MB-05	Double tube, MB	Adult Long	29 - 38 cm	5/Bag
DT3242MB-05	Double tube, MB	Large Adult	32 - 42 cm	5/Bag
DT3544MB-05	Double tube, MB	Large Adult Long	35 - 44 cm	5/Bag
DT4250MB-05	Double tube, MB	Thigh	42 - 50 cm	5/Bag

SOFTCHECK® CURVE™ DISPOSABLE BLOOD PRESSURE CUFFS



Worldwide obesity rates are climbing. This disease is strongly associated with increased mortality and morbidity, cardiovascular disease, diabetes, and many other illnesses. Blood pressure is the number one vital sign in detecting these issues and accurate blood pressure is critical. The Curve Cuffs are designed to fit a bariatric patient's arm correctly. Their conical shape fits snugly against the arm and properly compresses the artery.

- Will not slide during inflation
- Genuine Velcro® closures
- Manufactured without natural latex
- Adapters allow for compatibility with monitors worldwide



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Curve Disposable BP Cuffs

Single Tube



BRST3854HP-05	Single tube, HP	Extra Large Adult	38 - 54 cm	5/Bag
BRST3854DM-05	Single tube, DM	Extra Large Adult	38 - 54 cm	5/Bag
BRST3854FL-05	Single tube, FL	Extra Large Adult	38 - 54 cm	5/Bag
BRST3854ML-05	Single tube, ML	Extra Large Adult	38 - 54 cm	5/Bag
BRST3854MQ-05	Single tube, MQ	Extra Large Adult	38 - 54 cm	5/Bag
BRST3854MQML-05	Single tube, MQML	Extra Large Adult	38 - 54 cm	5/Bag
BRST3854-05	Single tube, no fitting	Extra Large Adult	38 - 54 cm	5/Bag

Double Tube



BRDT3854HP-05	Double tube, HP	Extra Large Adult	38 - 54 cm	5/Bag
BRDT3854DM-05	Double tube, DM	Extra Large Adult	38 - 54 cm	5/Bag
BRDT3854FL-05	Double tube, FL	Extra Large Adult	38 - 54 cm	5/Bag
BRDT3854ML-05	Double tube, ML	Extra Large Adult	38 - 54 cm	5/Bag
BRDT3854MQ-05	Double tube, MQ	Extra Large Adult	38 - 54 cm	5/Bag
BRDT3854MQMF-05	Double tube, MQMF	Extra Large Adult	38 - 54 cm	5/Bag
BRDT3854MQML-05	Double tube, MQML	Extra Large Adult	38 - 54 cm	5/Bag
BRDT3854-05	Double tube, no fitting	Extra Large Adult	38 - 54 cm	5/Bag

Manual Cuff



BRDT3854MB-05	Double tube, MB	Extra Large Adult	38 - 54 cm	5/Bag
---------------	-----------------	-------------------	------------	-------

SOFTCHECK® VINYL DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Vinyl Disposable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Comfortable soft vinyl
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Vinyl Disposable BP Cuffs



Bayonet (HP)

Single Tube					
VST0814HP-05	Single tube, HP	Infant	8 - 14 cm	5/Bag	
VST0814HP-20	Single tube, HP	Infant	8 - 14 cm	20/Box	
VST1320HP-05	Single tube, HP	Child	13 - 20 cm	5/Bag	
VST1320HP-20	Single tube, HP	Child	13 - 20 cm	20/Box	
VST1826HP-05	Single tube, HP	Small Adult	18 - 26 cm	5/Bag	
VST1826HP-20	Single tube, HP	Small Adult	18 - 26 cm	20/Box	
VST2635HP-05	Single tube, HP	Adult	26 - 35 cm	5/Bag	
VST2635HP-20	Single tube, HP	Adult	26 - 35 cm	20/Box	
VST2938HP-05	Single tube, HP	Adult Long	29 - 38 cm	5/Bag	
VST2938HP-20	Single tube, HP	Adult Long	29 - 38 cm	20/Box	
VST3242HP-05	Single tube, HP	Large Adult	32 - 42 cm	5/Bag	
VST3242HP-20	Single tube, HP	Large Adult	32 - 42 cm	20/Box	
VST3544HP-05	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag	
VST3544HP-20	Single tube, HP	Large Adult Long	35 - 44 cm	20/Box	
VST4250HP-05	Single tube, HP	Thigh	42 - 50 cm	5/Bag	
VST4250HP-20	Single tube, HP	Thigh	42 - 50 cm	20/Box	

SOFTCHECK® VINYL DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Vinyl Disposable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Comfortable soft vinyl
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Vinyl Disposable BP Cuffs

Single Tube				
VST0814ML-05	Single tube, ML	Infant	8 - 14 cm	5/Bag
VST0814ML-20	Single tube, ML	Infant	8 - 14 cm	20/Box
VST1320ML-05	Single tube, ML	Child	13 - 20 cm	5/Bag
VST1320ML-20	Single tube, ML	Child	13 - 20 cm	20/Box
VST1826ML-05	Single tube, ML	Small Adult	18 - 26 cm	5/Bag
VST1826ML-20	Single tube, ML	Small Adult	18 - 26 cm	20/Box
VST2635ML-05	Single tube, ML	Adult	26 - 35 cm	5/Bag
VST2635ML-20	Single tube, ML	Adult	26 - 35 cm	20/Box
VST2938ML-05	Single tube, ML	Adult Long	29 - 38 cm	20/Box
VST2938ML-20	Single tube, ML	Adult Long	29 - 38 cm	20/Box
VST3242ML-05	Single tube, ML	Large Adult	32 - 42 cm	5/Bag
VST3242ML-20	Single tube, ML	Large Adult	32 - 42 cm	20/Box
VST3544ML-05	Single tube, ML	Large Adult Long	35 - 44 cm	5/Bag
VST3544ML-20	Single tube, ML	Large Adult Long	35 - 44 cm	20/Box
VST4250ML-05	Single tube, ML	Thigh	42 - 50 cm	5/Bag
VST4250ML-20	Single tube, ML	Thigh	42 - 50 cm	20/Box



Male Luer (ML)

SOFTCHECK® VINYL DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Vinyl Disposable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Comfortable soft vinyl
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Vinyl Disposable BP Cuffs

Single Tube



Twist Lock (MQ)

VST0814MQ-05	Single tube, MQ	Infant	8 - 14 cm	5/Bag
VST0814MQ-20	Single tube, MQ	Infant	8 - 14 cm	20/Box
VST1320MQ-05	Single tube, MQ	Child	13 - 20 cm	5/Bag
VST1320MQ-20	Single tube, MQ	Child	13 - 20 cm	20/Box
VST1826MQ-05	Single tube, MQ	Small Adult	18 - 26 cm	5/Bag
VST2635MQ-05	Single tube, MQ	Adult	26 - 35 cm	5/Bag
VST2635MQ-20	Single tube, MQ	Adult	26 - 35 cm	20/Box
VST2938MQ-05	Single tube, MQ	Adult Long	29 - 38 cm	20/Box
VST2938MQ-20	Single tube, MQ	Adult Long	29 - 38 cm	20/Box
VST3242MQ-05	Single tube, MQ	Large Adult	32 - 42 cm	5/Bag
VST3242MQ-20	Single tube, MQ	Large Adult	32 - 42 cm	20/Box
VST4250MQ-05	Single tube, MQ	Thigh	42 - 50 cm	5/Bag
VST4250MQ-20	Single tube, MQ	Thigh	42 - 50 cm	20/Box



Male Twist Lock (MQML)

VST0814MQML-20	Single tube, MQ	Infant	8 - 14 cm	20/Box
VST1320MQML-20	Single tube, MQML	Child	13 - 20 cm	20/Box
VST1826MQML-20	Single tube, MQML	Small Adult	18 - 26 cm	20/Box
VST2635MQML-20	Single tube, MQML	Adult	26 - 35 cm	20/Box
VST2938MQML-20	Single tube, MQML	Adult Long	29 - 38 cm	20/Box
VST3242MQML-20	Single tube, MQML	Large Adult	32 - 42 cm	20/Box
VST3544MQML-20	Single tube, MQML	Thigh	35 - 44 cm	20/Box
VST4250MQML-20	Single tube, MQML	Thigh	42 - 50 cm	20/Box

SOFTCHECK® VINYL DISPOSABLE BLOOD PRESSURE CUFFS





SoftCheck Vinyl Disposable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Comfortable soft vinyl
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Vinyl Disposable BP Cuffs

Double Tube					
 Screw (DM)	VDT0814DM-05	Double tube, DM	Infant	8 - 14 cm	5/Bag
	VDT1320DM-05	Double tube, DM	Child	13 - 20 cm	5/Bag
	VDT1826DM-05	Double tube, DM	Small Adult	18 - 26 cm	5/Bag
	VDT2635DM-05	Double tube, DM	Adult	26 - 35 cm	5/Bag
	VDT2938DM-05	Double tube, DM	Adult Long	29 - 38 cm	5/Bag
	VDT3242DM-05	Double tube, DM	Large Adult	32 - 42 cm	5/Bag
	VDT4250DM-05	Double tube, DM	Thigh	42 - 50 cm	5/Bag
 No Fitting	VDT0814-20	Double tube, no fitting	Infant	8 - 14 cm	20/Box
	VDT1320-20	Double tube, no fitting	Child	13 - 20 cm	20/Box
	VDT1826-05	Double tube, no fitting	Small Adult	18 - 26 cm	5/Bag
	VDT2635-20	Double tube, no fitting	Adult	26 - 35 cm	20/Box
	VDT2938-20	Double tube, no fitting	Adult Long	29 - 38 cm	20/Box
	VDT3242-20	Double tube, no fitting	Large Adult	32 - 42 cm	20/Box
	VDT4250-20	Double tube, no fitting	Thigh	42 - 50 cm	20/Box

SOFTCHECK® YELLOW VINYL DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Yellow Vinyl Disposable Blood Pressure Cuffs intended for use in isolation areas are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Comfortable soft vinyl
- Color coded for accurate sizing
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Yellow Vinyl Disposable BP Cuffs

Single Tube



Bayonet (HP)

YVST1826HP-05	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
YVST2635HP-05	Single tube, HP	Adult	26 - 35 cm	5/Bag
YVST3242HP-05	Single tube, HP	Large Adult	32 - 42 cm	5/Bag

SOFTCHECK® FABRIC NEONATAL DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Neonatal Disposable Blood Pressure Cuffs are some of the world's finest quality blood pressure cuffs. In a NICU environment, this means better care for babies as well as greater effectiveness for practitioners.



- Meets AHA / AAMI standards
- Color coded for accurate sizing
- 5 sizes for accurate measurements
- Genuine low profile Velcro® closures

Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Fabric Neonatal Disposable BP Cuffs

Single Tube



Mini Male Luer
(ML316)

NN1ST-M316-10	Single tube, 3/16" ML	Size 1	3-6 cm	10/Box
NN2ST-M316-10	Single tube, 3/16" ML	Size 2	4-8 cm	10/Box
NN3ST-M316-10	Single tube, 3/16" ML	Size 3	6-11 cm	10/Box
NN4ST-M316-10	Single tube, 3/16" ML	Size 4	7-13 cm	10/Box
NN5ST-M316-10	Single tube, 3/16" ML	Size 5	8-15 cm	10/Box



Male Luer (ML)

NN1ST-10	Single tube, ML	Size 1	3-6 cm	10/Box
NN2ST-10	Single tube, ML	Size 2	4-8 cm	10/Box
NN3ST-10	Single tube, ML	Size 3	6-11 cm	10/Box
NN4ST-10	Single tube, ML	Size 4	7-13 cm	10/Box
NN5ST-10	Single tube, ML	Size 5	8-15 cm	10/Box



Bayonet (HP)

NN1ST-HP-10	Single tube, HP	Size 1	3-6 cm	10/Box
NN2ST-HP-10	Single tube, HP	Size 2	4-8 cm	10/Box
NN3ST-HP-10	Single tube, HP	Size 3	6-11 cm	10/Box
NN4ST-HP-10	Single tube, HP	Size 4	7-13 cm	10/Box
NN5ST-HP-10	Single tube, HP	Size 5	8-15 cm	10/Box

Double Tube



Male Luer (ML)

NN1DT-10	Double tube, ML	Size 1	3-6 cm	10/Box
NN2DT-10	Double tube, ML	Size 2	4-8 cm	10/Box
NN3DT-10	Double tube, ML	Size 3	6-11 cm	10/Box
NN4DT-10	Double tube, ML	Size 4	7-13 cm	10/Box
NN5DT-10	Double tube, ML	Size 5	8-15 cm	10/Box

SOFTCHECK® VINYL NEONATAL DISPOSABLE BLOOD PRESSURE CUFFS








SoftCheck Neonatal Disposable Blood Pressure Cuffs are some of the world's finest quality blood pressure cuffs. In a NICU environment, this means better care for babies as well as greater effectiveness for practitioners.



- Comfortable soft vinyl, created with babies' skin integrity in mind
- Easy to clean
- Color coded for accurate sizing
- 5 sizes for accurate measurements
- Genuine low profile Velcro® closures

Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Vinyl Neonatal Disposable BP Cuffs

Single Tube					
 (ML316)	VNN1ST-M316-10	Single tube, 3/16" ML	Size 1	3-6 cm	10/Box
	VNN2ST-M316-10	Single tube, 3/16" ML	Size 2	4-8 cm	10/Box
	VNN3ST-M316-10	Single tube, 3/16" ML	Size 3	6-11 cm	10/Box
	VNN4ST-M316-10	Single tube, 3/16" ML	Size 4	7-13 cm	10/Box
	VNN5ST-M316-10	Single tube, 3/16" ML	Size 5	8-15 cm	10/Box
 SnapQuik (SQ)	VNN1ST-SQ-10	Single tube, SQ	Size 1	3-6 cm	10/Box
	VNN2ST-SQ-10	Single tube, SQ	Size 2	4-8 cm	10/Box
	VNN3ST-SQ-10	Single tube, SQ	Size 3	6-11 cm	10/Box
	VNN4ST-SQ-10	Single tube, SQ	Size 4	7-13 cm	10/Box
	VNN5ST-SQ-10	Single tube, SQ	Size 5	8-15 cm	10/Box
 Male Luer (ML)	VNN1ST-10	Single tube, ML	Size 1	3-6 cm	10/Box
	VNN2ST-10	Single tube, ML	Size 2	4-8 cm	10/Box
	VNN3ST-10	Single tube, ML	Size 3	6-11 cm	10/Box
	VNN4ST-10	Single tube, ML	Size 4	7-13 cm	10/Box
	VNN5ST-10	Single tube, ML	Size 5	8-15 cm	10/Box
 Bayonet (HP)	VNN1ST-HP-10	Single tube, HP	Size 1	3-6 cm	10/Box
	VNN2ST-HP-10	Single tube, HP	Size 2	4-8 cm	10/Box
	VNN3ST-HP-10	Single tube, HP	Size 3	6-11 cm	10/Box
	VNN4ST-HP-10	Single tube, HP	Size 4	7-13 cm	10/Box
	VNN5ST-HP-10	Single tube, HP	Size 5	8-15 cm	10/Box
 Twist Lock (MQ)	VNN1ST-MQ-10	Single tube, MQ	Size 1	3-6 cm	10/Box
	VNN2ST-MQ-10	Single tube, MQ	Size 2	4-8 cm	10/Box
	VNN3ST-MQ-10	Single tube, MQ	Size 3	6-11 cm	10/Box
	VNN4ST-MQ-10	Single tube, MQ	Size 4	7-13 cm	10/Box
	VNN5ST-MQ-10	Single tube, MQ	Size 5	8-15 cm	10/Box

SOFTCHECK® VINYL NEONATAL DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Neonatal Disposable Blood Pressure Cuffs are some of the world's finest quality blood pressure cuffs. In a NICU environment, this means better care for babies as well as greater effectiveness for practitioners.



- Comfortable soft vinyl, created with babies' skin integrity in mind
- Easy to clean
- Color coded for accurate sizing
- 5 sizes for accurate measurements
- Genuine low profile Velcro® closures

Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Vinyl Neonatal Disposable BP Cuffs



Male Luer (ML)

Double Tube					
VNN1DT-10	Double tube, ML	Size 1	3-6 cm	10/Box	
VNN2DT-10	Double tube, ML	Size 2	4-8 cm	10/Box	
VNN3DT-10	Double tube, ML	Size 3	6-11 cm	10/Box	
VNN4DT-10	Double tube, ML	Size 4	7-13 cm	10/Box	
VNN5DT-10	Double tube, ML	Size 5	8-15 cm	10/Box	

SOFTCHECK FREE® NEONATAL DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Free Neonatal Disposable Blood Pressure Cuffs are some of the world's finest quality blood pressure cuffs. In a NICU environment, this means better care for babies as well as greater effectiveness for practitioners.



- Created with extremely soft material
- Color coded for accurate sizing
- 5 sizes for accurate measurements
- Genuine low profile Velcro® closures

Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Free Non-PVC Neonatal Disposable BP Cuffs

Single tube



PNN1ST-M316-10	Single tube, M316	Size 1	3-6 cm	10/Box
PNN2ST-M316-10	Single tube, M316	Size 2	4-8 cm	10/Box
PNN3ST-M316-10	Single tube, M316	Size 3	6-11 cm	10/Box
PNN4ST-M316-10	Single tube, M316	Size 4	7-13 cm	10/Box
PNN5ST-M316-10	Single tube, M316	Size 5	8-15 cm	10/Box



PNN1ST-10	Single tube, ML	Size 1	3-6 cm	10/Box
PNN2ST-10	Single tube, ML	Size 2	4-8 cm	10/Box
PNN3ST-10	Single tube, ML	Size 3	6-11 cm	10/Box
PNN4ST-10	Single tube, ML	Size 4	7-13 cm	10/Box
PNN5ST-10	Single tube, ML	Size 5	8-15 cm	10/Box



PNN1ST-HP-10	Single tube, HP	Size 1	3-6 cm	10/Box
PNN2ST-HP-10	Single tube, HP	Size 2	4-8 cm	10/Box
PNN3ST-HP-10	Single tube, HP	Size 3	6-11 cm	10/Box
PNN4ST-HP-10	Single tube, HP	Size 4	7-13 cm	10/Box
PNN5ST-HP-10	Single tube, HP	Size 5	8-15 cm	10/Box

Double Tube



PNN1DT-10	Double tube, ML	Size 1	3-6 cm	10/Box
PNN2DT-10	Double tube, ML	Size 2	4-8 cm	10/Box
PNN3DT-10	Double tube, ML	Size 3	6-11 cm	10/Box
PNN4DT-10	Double tube, ML	Size 4	7-13 cm	10/Box
PNN5DT-10	Double tube, ML	Size 5	8-15 cm	10/Box

Eighty-five percent of all blood pressure cuffs are made with the bayonet (HP) fitting. Spacelabs Healthcare can convert your monitors and sphygmomanometers to this fitting to save you time, energy and costs. Plus, your patient can move from Emergency, to the OR, to PACU, to the general floor with one cuff that attaches to every monitor in the hospital.



Order Code	Description	Qty
Blood Pressure Cuff Fittings & Adapters		
Used in conversions to HP Cuffs		
ZFT-HP-20	Male Bayonet Connector for Disposable Cuffs	20/Pack
ZFT-HP532-20	Male Bayonet Connector for Reusable Cuffs	20/Pack
ZFT-QC8220	Female Metal Bayonet Connector, 1/8" barb	1 Each
ZFT-QC8230	Female Metal Bayonet Connector, 5/32" barb	1 Each
ZFT-BPF230	Female Plastic Bayonet Connector, 1/8" barb	1 Each
ZFT-BPF240	Female Plastic Bayonet Connector, 5/32" barb	1 Each
ZHB-BPCOILADPRX1	Coiled Tube Manual Bulb Adapter w/metal female bayonet fitting	1 Each
ADP-182YBPF	Y Adapter w/1/8" barb fittings & plastic female bayonet connector	1 Each
ADP-182YRX4	Y Adapter w/1/8" barb fittings & metal female bayonet connector	1 Each
ADP-5322YBPF	Y Adapter w/5/32" barb fittings & plastic female bayonet connector	1 Each
ADP-5322YRX4	Y Adapter w/5/32" barb fittings & metal female bayonet connector	1 Each
ZFT-ML2YFL	Y Adapter w/Female Luer to 2 Male Luers	1 Each
Used in conversions to HP Cuffs-Straight		
ZFT-RX1DM	5" Tube w/Metal Female Bayonet and Screw Fitting	1 Each
ZFT-RX1FL	5" Tube w/Metal Female Bayonet and Female Luer	1 Each
ZFT-RX1ML	5" Tube w/Metal Female Bayonet and Male Luer	1 Each
ZFT-RX1MQ	5" Tube w/Metal Female Bayonet and Twist-Lock Fitting	1 Each
ZFT-RX1MQML	5" Tube w/Metal Female Bayonet and Male Twist-Lock Fitting	1 Each
Used in conversions to HP Cuffs-Ys		
ZFT-DM2YRX1	Y Adapter, Metal Female Bayonet to 2 Screw fittings	1 Each
ZFT-MQ2YRX1	Y Adapter w/Metal Female Bayonet to 2 Twist Lock fittings	1 Each
ZFT-MQML2YRX1	Y Adapter w/Metal Female Bayonet to 2 Male Twist Lock fittings	1 Each
ZFT-MQMF2YRX1	Y Adapter w/Metal Female Bayonet to one female Twist Lock and one male Twist Lock fitting	1 Each

Eighty-five percent of all blood pressure cuffs are made with the bayonet (HP) fitting. Spacelabs Healthcare can convert your monitors and sphygmomanometers to this fitting to save you time, energy and costs. Plus, your patient can move from Emergency, to the OR, to PACU, to the general floor with one cuff that attaches to every monitor in the hospital.



Order Code	Description	Qty
------------	-------------	-----

Blood Pressure Cuff Fittings & Adapters

Miscellaneous		
ZFT-CS901	Brass Male Screw Fitting	1 Each
ZFT-DM-20	Female Screw Connector for Disposable Cuffs	20/Pack
ZFT-DM532-20	Female Screw Connector for Reusable Cuffs	20/Pack
ZFT-DMML-20	Nylon Male Screw Fitting	20/Pack
ZFT-FL532-20	Female Luer for Disposable Cuffs	20/Pack
ZFT-FLL240-20	Female Luer for Reusable Cuffs	20/Pack
ZFT-HPDMAD-05	Y Adapter - Male Bayonet to 2 Male Screw	5/Pack
ZFT-ML532-20	Male Slip Luer for Disposable and Reusable Cuffs	20/Pack
ZFT-MLL532-20	Male Locking Luer	20/Pack
ZFT-MQ-20	Female Twist-Lock Fitting for Disposable and Reusable Cuffs	20/Pack
ZFT-MQMF2YHP	Y Adapter - Twist Lock and Male Twist Lock to Male Bayonet	1 Each
ZFT-MQML-20	Male Twist-Lock Fitting for Disposable Cuffs	20/Pack
ZFT-MQML5-20	Male Twist-Lock Fitting for Reusable Cuffs	20/Pack
ZFT-QC8030	Male Metal Bayonet, 5/32" barb	1 Each
ZGA-BPGBHB	BP Gauge w/Black 2 oz. Handbulb	1 Each
ZGA-BPHBRX1	Aneroid Gauge BPC Inflator, 8 in. tubing w/Metal Female Bayonet	1 Each
ZHB-4FBPADPRX1	4'. single tube Manual Bulb BP Cuff Adapter w/Metal Female Bayonet	1 Each
ZHB-HBTHM	Black 2oz Handbulb w/Metal Thumbwheel	1 Each

SUPPLIES & ACCESSORIES

PRESSURE INFUSORS







Unifusor is designed for rapid IV infusion, invasive pressure monitoring, and autotransfusion procedures; the Unifusor line is the clear choice in pressure infusors.



3 Models: Unifusor, Unifusor II & Unifusor Plus

- Manufactured without natural latex
- 100% leak tested during manufacturing process
- Reinforced translucent sleeve (see type of fluid and level through sleeve)
- IV pole loop
- Fluid bag hook

Order Code	Description	Range	Qty
Unifusor Disposable Pressure Infusor			
Disposable with aneroid gauge - sold as complete unit			
	803FGA-24 Unifusor with thumbwheel valve & aneroid gauge	500cc/ .5L	24/Box
	804FGA-24 Unifusor with thumbwheel valve & aneroid gauge	1000cc/ 1L	24/Box
	805FGA-10 Unifusor with thumbwheel valve & aneroid gauge	3000cc/ 3L	10/Box
	803SGA-24 Unifusor with stopcock valve & aneroid gauge	500cc/ .5L	24/Box
	804SGA-24 Unifusor with stopcock valve & aneroid gauge	1000cc/ 1L	24/Box
	805SGA-10 Unifusor with stopcock valve & aneroid gauge	3000cc/ 3L	10/Box
Disposable with piston gauge - sold as complete unit			
	903TGA-24 Unifusor Classic with thumbwheel valve & piston gauge	500cc/ .5L	24/Box
	904TGA-24 Unifusor Classic with thumbwheel valve & piston gauge	1000cc/ 1L	24/Box
	905TGA-10 Unifusor Classic with thumbwheel valve & piston gauge	3000cc/ 3L	10/Box
	903SGA-24 Unifusor Classic with stopcock valve & piston gauge	500cc/ .5L	24/Box
	904SGA-24 Unifusor Classic with stopcock valve & piston gauge	1000cc/ 1L	24/Box
	905SGA-10 Unifusor Classic with stopcock valve & piston gauge	3000cc/ 3L	10/Box

Unifusor is designed for rapid IV infusion, invasive pressure monitoring, and autotransfusion procedures; the Unifusor line is the clear choice in pressure infusors.



3 Models: Unifusor, Unifusor II & Unifusor Plus

- Manufactured without natural latex
- 100% leak tested during manufacturing process
- Reinforced translucent sleeve (see type of fluid and level through sleeve)
- IV pole loop
- Fluid bag hook

Order Code	Description	Range	Qty
------------	-------------	-------	-----

Unifusor II Semi-Disposable Pressure Infusor



Semi-Disposable (white bag with white hose - sold as two pieces in one bag)

2803FGA-24	Unifusor II with aneroid gauge & thumbwheel valve	500cc/ .5L	24/Box
2804FGA-24	Unifusor II with aneroid gauge & thumbwheel valve	1000cc/ 1L	24/Box
2805FGA-10	Unifusor II with aneroid gauge & thumbwheel valve	3000cc/ 3L	10/Box



2803SGA-24	Unifusor II with aneroid gauge & stopcock valve	500cc/ .5L	24/Box
2804SGA-24	Unifusor II with aneroid gauge & stopcock valve	1000cc/ 1L	24/Box
2805SGA-10	Unifusor II with aneroid gauge & stopcock valve	3000cc/ 3L	10/Box



2903TGA-24	Unifusor II with piston gauge & thumbwheel valve	500cc/ .5L	24/Box
2904TGA-24	Unifusor II with piston gauge & thumbwheel valve	1000cc/ 1L	24/Box
2905TGA-10	Unifusor II with piston gauge & thumbwheel valve	3000cc/ 3L	10/Box



2903SGA-24	Unifusor II with piston gauge & stopcock valve	500cc/ .5L	24/Box
2904SGA-24	Unifusor II with piston gauge & stopcock valve	1000cc/ 1L	24/Box
2905SGA-10	Unifusor II with piston gauge & stopcock valve	3000cc/ 3L	10/Box


Semi-Disposable replacement parts



2803X-24	Unifusor II Bag only	500cc/ .5L	24/Box
2804X-24	Unifusor II Bag only	1000cc/ 1L	24/Box
2805X-10	Unifusor II Bag only	3000cc/ 3L	10/Box



899F-04	Unifusor II aneroid gauge & thumbwheel valve assembly		4/Pkg
899S-04	Unifusor II aneroid gauge & stopcock valve assembly		4/Pkg



999T-04	Unifusor II piston gauge & thumbwheel valve assembly		4/Pkg
999S-04	Unifusor II piston gauge & stopcock valve assembly		4/Pkg

Unifusor is designed for rapid IV infusion, invasive pressure monitoring, and autotransfusion procedures; the Unifusor line is the clear choice in pressure infusors.



3 Models: Unifusor, Unifusor II & Unifusor Plus

- Manufactured without natural latex
- 100% leak tested during manufacturing process
- Reinforced translucent sleeve (see type of fluid and level through sleeve)
- IV pole loop
- Fluid bag hook

Order Code	Description	Range	Qty
------------	-------------	-------	-----


Unifusor II Semi-Reusable Pressure Infusor



Semi-Reusable (blue bag with white hose - sold as two pieces in one bag)

1103FGA-24	Unifusor II with aneroid gauge & thumbwheel valve	500cc/ .5L	24/Box
1104FGA-24	Unifusor II with aneroid gauge & thumbwheel valve	1000cc/ 1L	24/Box
1105FGA-10	Unifusor II with aneroid gauge & thumbwheel valve	3000cc/ 3L	10/Box
1103SGA-24	Unifusor II with aneroid gauge & stopcock valve	500cc/ .5L	24/Box
1104SGA-24	Unifusor II with aneroid gauge & stopcock valve	1000cc/ 1L	24/Box
1105SGA-10	Unifusor II with aneroid gauge & stopcock valve	3000cc/ 3L	10/Box
1103TGA-24	Unifusor II with piston gauge & thumbwheel valve	500cc/ .5L	24/Box
1104TGA-24	Unifusor II with piston gauge & thumbwheel valve	1000cc/ 1L	24/Box
1105TGA-10	Unifusor II with piston gauge & thumbwheel valve	3000cc/ 3L	10/Box
1113SGA-24	Unifusor II with piston gauge & stopcock valve	500cc/ .5L	24/Box
1114SGA-24	Unifusor II with piston gauge & stopcock valve	1000cc/ 1L	24/Box
1115SGA-10	Unifusor II with piston gauge & stopcock valve	3000cc/ 3L	10/Box

Semi-Reusable replacement parts






1103XB-05	Unifusor II Bag only	500cc/ .5L	5/Box
1104XB-05	Unifusor II Bag only	1000cc/ 1L	5/Box
1105XB-05	Unifusor II Bag only	3000cc/ 3L	5/Box
899F-04	Unifusor II aneroid gauge & thumbwheel valve assembly		4/Pkg
899S-04	Unifusor II aneroid gauge & stopcock valve assembly		4/Pkg
999T-04	Unifusor II piston gauge & thumbwheel valve assembly		4/Pkg
999S-04	Unifusor II piston gauge & stopcock valve assembly		4/Pkg

Unifusor is designed for rapid IV infusion, invasive pressure monitoring, and autotransfusion procedures; the Unifusor line is the clear choice in pressure infusors.





3 Models: Unifusor, Unifusor II & Unifusor Plus

- Manufactured without natural latex
- 100% leak tested during manufacturing process
- Reinforced translucent sleeve (see type of fluid and level through sleeve)
- IV pole loop
- Fluid bag hook


Order Code	Description	Range	Qty
Unifusor Plus Reusable Pressure Infusor			
Reusable (blue bag with black hose)			
	1103 Unifusor Plus with trigger / gauge assembly	500cc/ .5L	1 Each
	1104 Unifusor Plus with trigger / gauge assembly	1000cc/ 1L	1 Each
	1105 Unifusor Plus with trigger / gauge assembly	3000cc/ 3L	1 Each
	1103X Unifusor Plus Bag only	500cc/ .5L	1 Each
	1103X-05 Unifusor Plus Bag only	500cc/ .5L	5/Box
	1104X Unifusor Plus Bag only	1000cc/ 1L	1 Each
	1104X-05 Unifusor Plus Bag only	1000cc/ 1L	5/Box
	1105X Unifusor Plus Bag only	3000cc/ 3L	1 Each
	1105X-05 Unifusor Plus Bag only	3000cc/ 3L	5/Box
	1199 Trigger/ gauge assembly only		1 Each

SUPPLIES & ACCESSORIES ELECTROCARDIOGRAPHY



The ECG Combiner System can standardize Spacelabs monitors and Ultraview & Aria Telemetry transmitters to use one set of lead wires.

Additional cables, adapters and lead wires are also available to standardize most other manufacturer monitors and telemetry transmitters, allowing your hospital to use one set of lead wires per patient throughout the continuum of care.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Combiner Lead Wire Sets - requires ECG Combiner Cable or Adapter

DIN Shielded, Multi-Pinch, Disposable, AAMI



700-0006-30	3 Lead	Multi-Pinch, Disposable	127cm /50 in	25 Sets
700-0006-31	4 Lead (V-Lead)	Multi-Pinch, Disposable	127cm /50 in	25 Sets
700-0006-32	5 Lead	Multi-Pinch, Disposable	127cm /50 in	25 Sets
700-0006-33	5 Lead (V-Lead)	Multi-Pinch, Disposable	127cm /50 in	25 Sets
700-0006-34	6 Lead	Multi-Pinch, Disposable	127cm /50 in	25 Sets

DIN Shielded, Snap, Disposable, AAMI



700-0007-30	3 Lead	Snap, Disposable	127cm /50 in	25 Sets
700-0007-31	4 Lead (V-Lead)	Snap, Disposable	127cm /50 in	25 Sets
700-0007-32	5 Lead	Snap, Disposable	127cm /50 in	25 Sets
700-0007-33	5 Lead (V-Lead)	Snap, Disposable	127cm /50 in	25 Sets
700-0007-34	6 Lead	Snap, Disposable	127cm /50 in	25 Sets

DIN Shielded, Multi-Pinch, Reusable, AAMI



700-0006-35	3 Lead	Multi-Pinch, Reusable	74cm /29 in	1 Set
700-0006-36	4 Lead (V-Lead)	Multi-Pinch, Reusable	74cm /29 in	1 Set
700-0006-37	5 Lead	Multi-Pinch, Reusable	74cm /29 in	1 Set
700-0006-38	5 Lead (V-Lead)	Multi-Pinch, Reusable	74cm /29 in	1 Set
700-0006-39	6 Lead	Multi-Pinch, Reusable	74cm /29 in	1 Set
700-0010-00	3 Lead	Multi-Pinch, Reusable	130cm/51 in	1 Set
700-0010-01	5 Lead	Multi-Pinch, Reusable	130cm/51 in	1 Set
700-0010-02	5 Lead (V-Lead)	Multi-Pinch, Reusable	130cm/51 in	1 Set

DIN Shielded, Snap, Reusable, AAMI




700-0007-35	3 Lead	Snap, Reusable	74cm /29 in	1 Set
700-0007-36	4 Lead (V-Lead)	Snap, Reusable	74cm /29 in	1 Set
700-0007-37	5 Lead	Snap, Reusable	74cm /29 in	1 Set
700-0007-38	5 Lead (V-Lead)	Snap, Reusable	74cm /29 in	1 Set
700-0007-39	6 Lead	Snap, Reusable	74cm /29 in	1 Set

ECG COMBINER LEAD WIRE SYSTEM, AAMI



The ECG Combiner System can standardize Spacelabs monitors and Ultraview & Aria Telemetry transmitters to use one set of lead wires.

Additional cables, adapters and lead wires are also available to standardize most other manufacturer monitors and telemetry transmitters, allowing your hospital to use one set of lead wires per patient throughout the continuum of care.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Combiner Radiolucent Lead Wire Sets - requires ECG Combiner Cable or Adapter

DIN Shielded, Multi-Pinch, Reusable, AAMI



700-0010-03	3 Lead	Radiolucent, Multi-Pinch, Reusable	74cm /29 in	1 Set
700-0010-04	5 Lead	Radiolucent, Multi-Pinch, Reusable	74cm /29 in	1 Set
700-0010-05	3 Lead	Radiolucent, Multi-Pinch, Reusable	122cm/48 in	1 Set
700-0010-06	5 Lead	Radiolucent, Multi-Pinch, Reusable	122cm/48 in	1 Set
700-0010-07	3 Lead	Radiolucent, Multi-Pinch, Reusable	183cm/72 in	1 Set
700-0010-08	5 Lead	Radiolucent, Multi-Pinch, Reusable	183cm/72 in	1 Set

ECG Combiner Lead Wire Sets with Spacelabs Connection - plugs directly into Command Modules 90496/91496

Shielded, Spacelabs Connection, Multi-Pinch, Reusable, AAMI



700-0010-16	5 Lead	Spacelabs Connection, Multi-Pinch	60 cm/24 in	1 Set
700-0010-14	5 Lead	Spacelabs Connection, Multi-Pinch	102 cm/40 in	1 Set



Spacelabs Connection

Shielded, Spacelabs Connection, Snap, Reusable, AAMI



Snap

700-0010-15	5 Lead	Spacelabs Connection, Snap	60 cm/24 in	1 Set
-------------	--------	----------------------------	-------------	-------

ECG COMBINER LEAD WIRE SYSTEM, AAMI



The ECG Combiner System can standardize Spacelabs monitors and Ultraview & Aria Telemetry transmitters to use one set of lead wires.

Additional cables, adapters and lead wires are also available to standardize most other manufacturer monitors and telemetry transmitters, allowing your hospital to use one set of lead wires per patient throughout the continuum of care.



Reusable



Disposable



Multi-Pinch



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Combiner Cables

Spacelabs Connect Hardwired Cables for monitors using Command Module (90496/91496)



700-0008-55	3 Lead	AAMI	305 cm/10 ft	1 Cable
700-0008-57	5 Lead	AAMI	305 cm/10 ft	1 Cable
700-0008-59	6 Lead	AAMI	305 cm/10 ft	1 Cable
700-0008-56	10 Lead (5+5)	AAMI	305 cm/10 ft	1 Cable
700-0008-58	10 Lead (6+4)	AAMI	305 cm/10 ft	1 Cable

ECG Combiner Adapters

Spacelabs Monitors, AAMI



700-0008-78	3 Lead	Command Modules (90496/ 91496)		1 Adapter
700-0008-79	5 Lead	Command Modules (90496/ 91496)		1 Adapter

Spacelabs Telemetry, AAMI



700-0008-80	5 Lead	Ultraview Telemetry (9034X/ 9134X)		1 Adapter
700-0008-88	5 Lead	Aria Telemetry (96281)		1 Adapter

The ECG Combiner System can standardize Spacelabs monitors and Ultraview & Aria Telemetry transmitters to use one set of lead wires.

Additional cables, adapters and lead wires are also available to standardize most other manufacturer monitors and telemetry transmitters, allowing your hospital to use one set of lead wires per patient throughout the continuum of care.



Reusable



Disposable



Multi-Pinch



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Combiner Cables for Standardization


Other manufacturers



700-0008-30	3 Lead	Datascope	305 cm/10 ft	1 Cable
700-0008-31	5 Lead	Datascope	305 cm/10 ft	1 Cable
700-0008-33	3 Lead	Datex-Ohmeda	305 cm/10 ft	1 Cable
700-0008-34	5 Lead	Datex-Ohmeda	305 cm/10 ft	1 Cable
700-0008-35	6 lead	Datex-Ohmeda	305 cm/10 ft	1 Cable
700-0008-36	3 Lead	GE	305 cm/10 ft	1 Cable
700-0008-37	5 Lead	GE	305 cm/10 ft	1 Cable
700-0008-38	10 Lead (5+5)	GE	305 cm/10 ft	1 Cable
700-0008-39	6 Lead	GE	305 cm/10 ft	1 Cable
700-0008-40	10 Lead (6+4)	GE	305 cm/10 ft	1 Cable
700-0008-47	3 Lead	Philips	305 cm/10 ft	1 Cable

The ECG Combiner System can standardize Spacelabs monitors and Ultraview & Aria Telemetry transmitters to use one set of lead wires.

Additional cables, adapters and lead wires are also available to standardize most other manufacturer monitors and telemetry transmitters, allowing your hospital to use one set of lead wires per patient throughout the continuum of care.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Combiner Cables for Standardization



Other manufacturers					
700-0008-48	10 Lead (5+5)	Philips	305 cm/10 ft	1 Cable	
700-0008-49	5 Lead	Philips	305 cm/10 ft	1 Cable	
700-0008-50	10 Lead (6+4)	Philips	305 cm/10 ft	1 Cable	
700-0008-51	6 Lead	Philips	305 cm/10 ft	1 Cable	
700-0008-41	3 Lead	Medtronic LifePak	305 cm/10 ft	1 Cable	
700-0008-42	5 Lead	Medtronic LifePak	305 cm/10 ft	1 Cable	
700-0008-43	6 lead	Medtronic LifePak	305 cm/10 ft	1 Cable	
700-0008-44	3 Lead	Nihon Kohden	305 cm/10 ft	1 Cable	
700-0008-45	5 Lead	Nihon Kohden	305 cm/10 ft	1 Cable	
700-0008-46	6 lead	Nihon Kohden	305 cm/10 ft	1 Cable	
700-0008-52	3 Lead	Zoll, 90° angle	305 cm/10 ft	1 Cable	
700-0008-53	5 Lead	Zoll, 90° angle	305 cm/10 ft	1 Cable	
700-0008-54	6 lead	Zoll, 90° angle	305 cm/10 ft	1 Cable	

The ECG Combiner System can standardize Spacelabs monitors and Ultraview & Aria Telemetry transmitters to use one set of lead wires.

Additional cables, adapters and lead wires are also available to standardize most other manufacturer monitors and telemetry transmitters, allowing your hospital to use one set of lead wires per patient throughout the continuum of care.



Reusable



Disposable



Multi-Pinch



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Combiner Adapters for Standardization

Other manufacturers

	700-0008-60	5 Lead	Datascope Telemetry	1 Adapter
	700-0008-61	3 Lead	Datascope	1 Adapter
	700-0008-62	5 Lead	Datascope	1 Adapter
	700-0008-63	3 Lead	Datex-Ohmeda	1 Adapter
	700-0008-64	5 Lead	Datex-Ohmeda	1 Adapter
	700-0008-65	5 Lead	GE Apex Pro FH Telemetry	1 Adapter
	700-0008-66	6 Lead	GE Apex Pro FH Telemetry	1 Adapter

The ECG Combiner System can standardize Spacelabs monitors and Ultraview & Aria Telemetry transmitters to use one set of lead wires.

Additional cables, adapters and lead wires are also available to standardize most other manufacturer monitors and telemetry transmitters, allowing your hospital to use one set of lead wires per patient throughout the continuum of care.



Reusable



Disposable



Multi-Pinch



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Combiner Adapters for Standardization

Other manufacturers

	700-0008-67	3 Lead	GE	1 Adapter
	700-0008-68	5 Lead	GE	1 Adapter
	700-0008-69	6 Lead	GE Monitors	1 Adapter
	700-0008-70	3 Lead	Nihon Kohden	1 Adapter
	700-0008-71	6 Lead	Nihon Kohden	1 Adapter
	700-0008-72	3/5 Lead	Philips Telemetry	1 Adapter
	700-0008-73	3 Lead	Philips	1 Adapter
	700-0008-74	5 Lead	Philips	1 Adapter
	700-0008-75	3 Lead	Siemens Drager	1 Adapter
	700-0008-76	5 Lead	Siemens Drager	1 Adapter

Shielded ECG lead wires and cables for use with Spacelabs Command Modules 90496 & 91496 and mCare Monitor 91220 (3- and 5- lead wire sets only).



Shielded

Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Lead Wire Sets, Shielded, Reusable, AAMI



Snap

Snap Sets

700-0007-00	3 Lead, Shielded	61 cm/24 in	1 Set
700-0007-02	3 Lead, Shielded	102 cm/40 in	1 Set
700-0007-08	5 Lead, Shielded	61 cm/24 in	1 Set
700-0007-10	5 Lead, Shielded	102 cm/40 in	1 Set
700-0007-16*	10 Lead, limb/chest, Shielded	61 cm/24 in	1 Set
700-0007-18*	10 Lead, limb/chest, Shielded	102 cm/40 in	1 Set

*10 Lead ECG Lead Wire Sets not compatible with mCare Monitor 91220



Pinch

Pinch Sets

700-0006-00	3 Lead, Shielded	61 cm/24 in	1 Set
700-0006-02	3 Lead, Shielded	102 cm/40 in	1 Set
700-0006-08	5 Lead, Shielded	61 cm/24 in	1 Set
700-0006-10	5 Lead, Shielded	102 cm/40 in	1 Set
700-0006-16	10 Lead, limb/chest, Shielded	61 cm/24 in	1 Set
700-0006-18	10 Lead, limb/chest, Shielded	102 cm/40 in	1 Set

*10 Lead ECG Lead Wire Sets not compatible with mCare Monitor 91220

Shielded ECG lead wires and cables for use with Spacelabs Command Modules 90496 & 91496 and mCare Monitor 91220 (3- and 5- lead wire sets only).



Shielded

Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Lead Wires, Individual Colors, Shielded, Reusable, AAMI



Snap

24" Snap Individual Color

012-0498-00	Red	61 cm/24 in	3/Pack
012-0498-04	Black	61 cm/24 in	3/Pack
012-0498-34	Black	61 cm/24 in	5/Pack
012-0498-05	Brown	61 cm/24 in	3/Pack
012-0498-03	White	61 cm/24 in	3/Pack
012-0498-01	Green	61 cm/24 in	3/Pack
012-0498-11	Yellow/Brown	61 cm/24 in	3/Pack
012-0498-12	Green/Brown	61 cm/24 in	3/Pack
012-0498-13	Blue/Brown	61 cm/24 in	3/Pack
012-0498-14	Orange/Brown	61 cm/24 in	3/Pack
012-0498-15	Violet/Brown	61 cm/24 in	3/Pack



Snap

40" Snap Individual Color

012-0498-06	Red	102 cm/40 in	3/Pack
012-0498-09	Black	102 cm/40 in	3/Pack
012-0498-39	Black	102 cm/40 in	5/Pack
012-0498-27	Brown	102 cm/40 in	3/Pack
012-0498-08	White	102 cm/40 in	3/Pack
012-0498-07	Green	102 cm/40 in	3/Pack

Shielded ECG lead wires and cables for use with Spacelabs Command Modules 90496 & 91496 and mCare Monitor 91220 (3- and 5- lead wire sets only).



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Cables, Shielded, AAMI



With 17 PIN Adapter

700-0008-08	3 Lead, Shielded	305 cm/10 ft	1 Cable
700-0008-06	5 Lead, Shielded	305 cm/10 ft	1 Cable
700-0008-00	10 Lead, Shielded	305 cm/10 ft	1 Cable

ECG Accessories



Electrodes

369062-001	ECG electrode, tape, single use	Adult	100/Box
690-0262-00	ECG electrode, transparent tape, general/stress	Adult	30/Box

Shielded ECG lead wires and cables for use with Spacelabs Command Modules 90496 & 91496 and mCare Monitor 91220 (3- and 5- lead wire sets only).



Shielded

Order Code	Description	Length	Qty
------------	-------------	--------	-----

Neonatal ECG Supplies & Accessories



DIN Unshielded, Neonatal ECG Cable

700-0008-12	3 Lead		1 Cable
-------------	--------	--	---------



DIN Unshielded, Neonatal Electrodes


685-0037-00	Pre-wired neonatal electrodes, single use Pre-gelled. Silver/Silver Chloride	3-Lead	100 Sets/Box
-------------	---	--------	--------------



Electrode

369063-001	ECG neonatal electrode, single use, Individually packaged. Pre-gelled. Silver/Silver Chloride.	Ped/Neo	100/Box
------------	--	---------	---------

DIN Unshielded ECG lead wires and cables for use with Spacelabs Multiparameter Module 90470 ; available in 3- and 5-lead configurations.



DIN Unshielded

Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Lead Wire Sets, DIN Unshielded, Reusable, AAMI

Snap Sets



Snap

012-0283-00	3 Lead	46 cm/18 in	1 Set
012-0283-11	3 Lead	61 cm/24 in	1 Set
012-0283-02	3 Lead	102 cm/40 in	1 Set
012-0285-00	5 Lead	46 cm/18 in	1 Set
012-0285-11	5 Lead	61 cm/24 in	1 Set
012-0285-02	5 Lead	102 cm/40 in	1 Set


Pinch Sets



Pinch

012-0284-00	3 Lead	46 cm/18 in	1 Set
012-0284-11	3 Lead	61 cm/24 in	1 Set
012-0284-02	3 Lead	102 cm/40 in	1 Set
012-0286-00	5 Lead	46 cm/18 in	1 Set
012-0286-11	5 Lead	61 cm/24 in	1 Set
012-0286-02	5 Lead	102 cm/40 in	1 Set

DIN Unshielded ECG lead wires and cables for use with Spacelabs Multiparameter Module 90470 ; available in 3- and 5-lead configurations.



DIN Unshielded

Order Code	Description	Length	Qty
------------	-------------	--------	-----


ECG Lead Wires, Individual Colors, DIN Unshielded, Reusable, AAMI



Snap

Snap Individual Color			
012-0295-00	Red	46 cm/18 in	3/Pack
012-0295-10	Black	46 cm/18 in	3/Pack
012-0295-20	Brown	46 cm/18 in	3/Pack
012-0295-30	White	46 cm/18 in	3/Pack
012-0295-40	Green	46 cm/18 in	3/Pack
012-0295-06	Red	61 cm/24 in	3/Pack
012-0295-26	Black	61 cm/24 in	3/Pack
012-0295-46	Brown	61 cm/24 in	3/Pack
012-0295-16	White	61 cm/24 in	3/Pack
012-0295-36	Green	61 cm/24 in	3/Pack
012-0295-02	Red	102 cm/40 in	3/Pack
012-0295-22	Black	102 cm/40 in	3/Pack
012-0295-42	Brown	102 cm/40 in	3/Pack
012-0295-12	White	102 cm/40 in	3/Pack
012-0295-32	Green	102 cm/40 in	3/Pack

DIN Unshielded ECG lead wires and cables for use with Spacelabs Multiparameter Module 90470 ; available in 3- and 5- lead configurations.



DIN Unshielded

Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Lead Wires, Individual Colors, DIN Unshielded, Reusable, AAMI



Pinch

Pinch Individual Color

012-0296-00	Red	46 cm/18 in	3/Pack
012-0296-20	Black	46 cm/18 in	3/Pack
012-0286-00	Brown	46 cm/18 in	3/Pack
012-0296-10	White	46 cm/18 in	3/Pack
012-0296-04	Red	61 cm/24 in	3/Pack
012-0296-24	Black	61 cm/24 in	3/Pack
012-0296-44	Brown	61 cm/24 in	3/Pack
012-0296-14	White	61 cm/24 in	3/Pack
012-0296-34	Green	61 cm/24 in	3/Pack
012-0296-02	Red	102 cm/40 in	3/Pack
012-0296-22	Black	102 cm/40 in	3/Pack
012-0296-42	Brown	102 cm/40 in	3/Pack
012-0296-12	White	102 cm/40 in	3/Pack
012-0296-32	Green	102 cm/40 in	3/Pack

DIN Unshielded ECG lead wires and cables for use with Spacelabs Multiparameter Module 90470 ; available in 3- and 5-lead configurations.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Adapter and Cable, AAMI



Required for use with Multiparameter Module 90470

103-0059-00	17-pin to 6-pin adapter, allows the ECG cables designed for Command Modules 90496/91496 to work on the Multiparameter Module 90470		1 Each
-------------	--	--	--------

SUPPLIES & ACCESSORIES

TELEMETRY



For use with Spacelabs Telemetry Transmitter s (90341, 90343, 90347, 91341, 91343, 91347 and 96218).

Accessories that help protect your equipment investment and extend product life.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Lead Wire Sets, DIN Unshielded, Reusable, AAMI



Snap

Snap Sets

012-0285-11	5 Lead	61 cm/24 in	1 Set
012-0605-10	5 Lead	64 cm/25.2 in	1 Set



Pinch

Pinch Sets

012-0286-11	5 Lead	61 cm/24 in	1 Set
-------------	--------	-------------	-------

ECG Lead Wires, Individual Colors, DIN Unshielded, Reusable, AAMI



Snap

Snap Individual Color

012-0295-06	Red	61 cm/24 in	3/Pack
012-0295-26	Black	61 cm/24 in	3/Pack
012-0295-46	Brown	61 cm/24 in	3/Pack
012-0295-16	White	61 cm/24 in	3/Pack
012-0295-36	Green	61 cm/24 in	3/Pack

Snap, Individual Color

For use with Spacelabs Ultraview Telemetry 9034X/9134X only



Snap

700-0010-12	Red	74 cm/29 in	5/Pack
700-0010-09	Black	74 cm/29 in	5/Pack
700-0010-10	Brown	74 cm/29 in	5/Pack
700-0010-13	White	74 cm/29 in	5/Pack
700-0010-11	Green	74 cm/29 in	5/Pack

For use with Spacelabs Telemetry Transmitters (90341, 90343, 90347, 91341, 91343, 91347 and 96218).

Accessories that help protect your equipment investment and extend product life.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Lead Wires, Individual Colors, DIN Unshielded, Reusable, AAMI



Pinch

Pinch Individual Color			
012-0296-04	Red	61 cm/24 in	3/Pack
012-0296-14	White	61 cm/24 in	3/Pack
012-0296-24	Black	61 cm/24 in	3/Pack
012-0296-34	Green	61 cm/24 in	3/Pack
012-0296-44	Brown	61 cm/24 in	3/Pack

Telemetry Transmitter Adapters

For use with ECG Combiner Lead Wires



700-0008-80	Ultraview Telemetry Combiner Adapter (9034X/ 9134X)	1 Adapter
700-0008-88	Aria Telemetry Combiner Adapter (96281)	1 Adapter



For use with DIN Unshielded Lead Wires



706-0191-00	Aria Telemetry Grouper (96281)	5/Bag
-------------	--------------------------------	-------

For use with Spacelabs Telemetry Transmitter s (90341, 90343, 90347, 91341, 91343, 91347 and 96218).

Accessories that help protect your equipment investment and extend product life.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

Telemetry Accessories

Disposable Pouches and Electrodes



015-0500-02	Ultraview Telemetry pouch (90347)		100/Box
369062-001	ECG electrode	Adult	100/Box
690-0262-00	ECG transparent electrode, general/stress	Adult	30/Box

For use with Spacelabs Telemetry Transmitter s (90341, 90343, 90347, 91341, 91343, 91347 and 96218).

Accessories that help protect your equipment investment and extend product life.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

SpO₂ Cable & Sensors

SpO₂ Cable



700-0014-00	TruLink SpO ₂ Sensor Adapter Cable Compatible with SL multi-parameter telemetry transmitters (90343/91343/96218) and TruLink SpO ₂ sensors	60 cm/24 in	1 Cable
-------------	--	-------------	---------

Reusable SpO₂ Sensors



015-0660-00	TruLink SpO ₂ finger sensor	Adult >30 kg/66 lbs	1 Each
-------------	--	------------------------	--------

Disposable SpO₂ Sensors



015-0662-00	TruLink SpO ₂ foam sensor	Adult > 30 kg/66 lbs	24/Box
-------------	--------------------------------------	-------------------------	--------



015-0663-00	TruLink SpO ₂ vinyl sensor	Adult >30 kg/66 lbs	24/Box
-------------	---------------------------------------	------------------------	--------



015-0664-00	TruLink SpO ₂ foam sensor	Pediatric 10-50 kg/22-110 lbs	24/Box
-------------	--------------------------------------	----------------------------------	--------



015-0665-00	TruLink SpO ₂ vinyl sensor	Pediatric 10-50 kg/22-110 lbs	24/Box
-------------	---------------------------------------	----------------------------------	--------



015-0667-00	TruLink SpO ₂ cloth sensor	Neonatal <3 kg/7 lbs	24/Box
-------------	---------------------------------------	-------------------------	--------

For use with Spacelabs Telemetry Transmitter s (90341, 90343, 90347, 91341, 91343, 91347 and 96218).

Accessories that help protect your equipment investment and extend product life.




Order Code	Description	Length	Qty
------------	-------------	--------	-----

ABP Cable

ABP Cable			
700-0015-00	TruLink ABP Interface Cable Compatible with SL multi-parameter telemetry transmitters (90343/91343/96217)	60 cm/24 in	1 Cable

SUPPLIES & ACCESSORIES

CAPNOGRAPHY & MULTIGAS



The Nomoline™ multi-patient use adapter for sidestream monitoring offers a cost effective solution for both intubated and non-intubated patients. The sophisticated technology removes water from the sample line and reduces potential for water intrusion and cross contamination. Utilizing the Nomoline multi-patient use adapter in conjunction with single use CO₂ cannulas, we are able to offer a capnography solution at a much lower cost..



(92516) Capno Pod (92517) Capnography Module (92518) Multi-Gas Module

Order Code	Description	Qty
------------	-------------	-----

Sidestream Monitoring Supplies

Multi-patient use Nomoline Adapter - required for sidestream monitoring



103-0234-00	Sidestream Nomoline adapter with female luer lock. To be used with single patient use cannulas , sidestream sampling lines and T-adapters. T-adapters require a sampling line for use.	Multi-use	25/Box
-------------	--	-----------	--------



015-0683-00	Sidestream Nomoline adapter with sampling line	Single Use	25/Box
-------------	--	------------	--------

CO2 Nasal Cannula with male luer connector



704-0176-00	Adult	2 m/6.5 ft	25/Box
704-0177-00	Pediatric	2 m/6.5 ft	25/Box
704-0178-00	Infant	2 m/6.5 ft	25/Box
704-0184-00	Adult	3 m/10 ft	25/Box
704-0185-00	Pediatric	3 m/10 ft	25/Box

CO2 Nasal Cannula w/O2 Delivery and male luer connector



704-0181-00	Adult	2 m/7 ft	25/Box
704-0180-00	Pediatric	2 m/7 ft	25/Box
704-0179-00	Infant	2 m/7 ft	25/Box
704-0186-00	Adult	3 m/10 ft	25/Box

CO2 Oral/Nasal Cannula w/O2 Delivery and male luer connector



704-0183-00	Adult	2 m/7 ft	25/Box
704-0182-00	Pediatric	2 m/7 ft	25/Box
704-0187-00	Adult	3 m/10 ft	25/Box
704-0189-00	Pediatric	3 m/10 ft	25/Box

The Nomoline™ multi-patient use adapter for sidestream monitoring offers a cost effective solution for both intubated and non-intubated patients. The sophisticated technology removes water from the sample line and reduces potential for water intrusion and cross contamination. Utilizing the Nomoline multi-patient use adapter in conjunction with single use CO₂ cannulas, we are able to offer a capnography solution at a much lower cost..



(92516) Capno Pod (92517) Capnography Module (92518) Multi-Gas Module

Order Code	Description	Qty
------------	-------------	-----

Sidestream Monitoring Supplies - Non-intubated patients



Emergency CO2 Nasal Cannula

704-0193-00	Adult	2 m/6.5 ft	25/Box
-------------	-------	------------	--------



Sampling Lines & Airway Adapters, single use

704-0191-00	Adult CO2 Sample Line w/ T-Adapter	2 m/6.5 ft	20/Box
-------------	------------------------------------	------------	--------

704-0192-00	Adult CO2 Sample Line w/ T-Adapter	3 m/10 ft	20/Box
-------------	------------------------------------	-----------	--------



166-7085-00	Sidestream sampling line with male luer lock on both ends	2 m/6.5 ft	25/Box
-------------	---	------------	--------

166-7121-00	Sidestream sampling line with male luer lock on both ends	3 m/10 ft	50/Box
-------------	---	-----------	--------



704-0175-00	Sidestream airway adapter/T-Adapter. Requires a sampling line for use.		25/Box
-------------	---	--	--------

Mainstream Monitoring Supplies for use only with the Capnography Module 92517 on intubated patients



Sensor & Airway Adapters - required for mainstream monitoring

010-1980-00	Mainstream sensor, reusable		1 Each
-------------	-----------------------------	--	--------



704-0173-00	Mainstream airway adapter, single use	Adult/Pediatric	25/Box
-------------	---------------------------------------	-----------------	--------







704-0174-00	Mainstream airway adapter, single use	Infant	10/Box
-------------	---------------------------------------	--------	--------

The Dual Capnography Module (91517) provides both sidestream and mainstream EtCO₂ capabilities in a single module. Supplies available for all patient types.









(91517)
Dual Capnography Module

Order Code	Description	Qty
Sidestream, Non-intubated, single patient use		
	CO₂ Nasal Cannula	
704-0010-00	Adult	10/Box
704-0012-00	Pediatric	10/Box
704-0011-00	Infant/Neonatal	10/Box
	CO₂ Nasal Cannula w/O₂ Delivery	
704-0014-00	Adult	10/Box
704-0016-00	Pediatric	10/Box
704-0015-00	Infant/Neonatal	10/Box
	CO₂ Oral/Nasal Cannula	
704-0017-00	Adult	10/Box
704-0019-00	Pediatric	10/Box
	CO₂ Oral/Nasal Cannula w/O₂ Delivery	
704-0020-00	Adult	10/Box
704-0022-00	Pediatric	10/Box

The Dual Capnography Module (91517) provides both sidestream and mainstream EtCO₂ capabilities in a single module. Supplies available for all patient types.



(91517)
Dual Capnography Module

Order Code	Description		Qty
Sidestream, Intubated, single patient use			
Sidestream - short term up to 12 hours (patient non-humidified)			
	704-0023-00 Airway adapter + sampling line, ET tube > 4.0 mm	Adult/Ped	10/Box
	704-0024-00 Airway adapter + sampling line, ET tube ≤ 4.0 mm	Ped/Infant	10/Box
	704-0027-00 CO ₂ sampling line w/male luer <i>Recommended for heat moisture exchange filters</i>		10/Box
Sidestream - long term up to 120 hours (patient humidified)			
	704-0025-00 Airway adapter + sampling line with Nafion®, ET tube > 4.0 mm	Adult/Ped	10/Box
	704-0026-00 Airway adapter + sampling line with Nafion	Ped/Infant	10/Box
	704-0028-00 CO ₂ sampling line w/male luer and Nafion <i>Recommended for heat moisture exchange filters</i>		10/Box

The Dual Capnography Module (91517) provides both sidestream and mainstream EtCO₂ capabilities in a single module. Supplies available for all patient types.



(91517)
Dual Capnography Module

Order Code	Description	Qty
------------	-------------	-----

Mainstream, for Intubated Patients Only

Mainstream Cable with Sensor



704-0001-00	CO ₂ sensor, reusable	1 Each
-------------	----------------------------------	--------



704-0006-00	Sensor cable holding clip	5/Pack
-------------	---------------------------	--------

Mainstream Airway Adapters, for use with CO₂ Sensor 704-0001-00



704-0005-00	CO ₂ airway adapter, single use	Adult	24/Box
-------------	--	-------	--------



704-0004-00	CO ₂ airway adapter, single use	Neonate	24/Box
-------------	--	---------	--------



704-0002-00	CO ₂ airway adapter, reusable	Adult	1 Each
-------------	--	-------	--------



704-0003-00	CO ₂ airway adapter, reusable	Neonate	1 Each
-------------	--	---------	--------

For use with Ultraview SL™ Multigas Analyzers (90518 & 91518).



(90518)




(91518)

Multigas Analyzers

Order Code	Description	Qty
------------	-------------	-----

Multigas Supplies & Accessories

Anesthetic Agent Samples Lines



015-0312-01	Anesthetic agent sample line with gas sampling tee and airway adapter, single use, Sampling tee: 22mm OD/15 mm ID x 15mm OD Sampling line: 1.22 mm/0.05 inches ID tubing (90518/91518)	213 cm/7 ft.	50/Box
-------------	--	--------------	--------



015-0313-01	Anesthetic agent sample line with gas sampling tee and airway adapter, single use, Sampling tee: 22mm OD/15 mm ID x 15mm OD Sampling line: 1.22 mm/0.05 inches ID tubing (90518/91518)	305 cm/10 ft.	50/Box
-------------	--	---------------	--------



162-0046-01	Anesthetic agent sample line with a male luer connector on each end, single use, 1.22 mm/.048 in ID tubing (90518/91518)	305 cm/10 ft.	50/Box
-------------	--	---------------	--------



162-0047-01	Anesthetic agent sample line with a male luer connector on each end, single use, 1.22 mm/.048 in ID tubing (90518/91518)	213 cm/7 ft.	50/Box
-------------	--	--------------	--------

Exhaust Scavenging



717-0408-00	Exhaust scavenging line with female luer connection on both ends, single use (90518)	2 m/7 ft.	5/Pack
-------------	--	-----------	--------

MULTIGAS ANALYZER SUPPLIES & ACCESSORIES



For use with Ultraview SL™ Multigas Analyzers (90518 & 91518).



(90518)



(91518)

Multigas Analyzers

Order Code	Description	Qty
------------	-------------	-----

Multigas Supplies & Accessories

Water Traps, Reusable



717-0447-00	D-FEND water trap cartridge, Black. Replace at least once every two months or when occluded (91518)	10/Pack
-------------	--	---------



717-0446-00	D-FEND water trap container (91518)	5/Pack
-------------	-------------------------------------	--------



717-0448-00	D-FEND+ water trap cartridge, green Replace at least once every 24 hours or when occluded (91518)	10/Pack
-------------	--	---------



Water Traps, Single Use

378-0050-00	TruLink water trap (90518)	1 Each
-------------	----------------------------	--------

For use with Ultraview SL™ Multigas Analyzers (90518 & 91518).



(90518)




(91518)

Multigas Analyzers

Order Code	Description	Qty
------------	-------------	-----

Multigas Supplies & Accessories

Nafion® Tubing



369014-003	Nafion tubing, single use Use it in-line with sampling line to reduce moisture and prevent clogging (90518)	61 cm/24 in	1 Each
------------	--	-------------	--------

Filter Cartridge



016-0428-00	Nafion tubing, single use Use it in-line with sampling line to reduce moisture and prevent clogging (90518)		1 Each
-------------	--	--	--------



050-1120-00	O ₂ sensor kit, micro-fuel cell (90518-B)		1 Each
-------------	--	--	--------

SUPPLIES & ACCESSORIES PULSE OXIMETRY



For use with Spacelabs Command Modules (90496 & 91496-U).
Cables and connectors are color-coded blue.

Full line of industry standard oximetry products that support your needs. TruLink SpO₂ sensors are available in a wide range of sizes and shapes to fit all patients.



Order Code	Description	Size	Qty
------------	-------------	------	-----

Spacelabs TruLink Technology

Cables



700-0030-00	TruLink SpO ₂ cable	305 cm/10 ft.	1 Each
-------------	--------------------------------	---------------	--------

Reusable sensors



015-0660-00	TruLink SpO ₂ finger sensor	Adult >30 kg/66 lbs	1 Each
-------------	--	------------------------	--------



015-0661-00	TruLink SpO ₂ Y multi-site sensor Includes: 1 Sensor 1 ear clip 10 adult finger/neonate feet/hand wraps 10 pediatric finger or infant toe adhesive wraps	Universal >30 kg/66 lbs <30 kg/66 lbs	1 Each
-------------	--	---	--------

For use with Spacelabs Command Modules (90496 & 91496-U).
Cables and connectors are color-coded blue.

Full line of industry standard oximetry products that support your needs. TruLink SpO₂ sensors are available in a wide range of sizes and shapes to fit all patients.



Order Code	Description	Size	Qty
------------	-------------	------	-----

Spacelabs TruLink Technology

Disposable Sensors

	015-0662-00	TruLink SpO ₂ foam sensor	Adult > 30 kg/66 lbs	24/Box
	015-0663-00	TruLink SpO ₂ vinyl sensor	Adult >30 kg/66 lbs	24/Box
	015-0664-00	TruLink SpO ₂ foam sensor	Pediatric 10-50 kg/22-110 lbs	24/Box
	015-0665-00	TruLink SpO ₂ vinyl sensor	Pediatric 10-50 kg/22-110 lbs	24/Box
	015-0666-00	TruLink SpO ₂ finger sensor	Infant 3-20 kg/7-44 lbs	24/Box
	015-0667-00	TruLink SpO ₂ cloth sensor	Neonatal <3 kg/7 lbs	24/Box

For use with Spacelabs Command Modules (90496 & 91496-U).
Cables and connectors are color-coded blue.

Full line of industry standard oximetry products that support your needs. TruLink SpO₂ sensors are available in a wide range of sizes and shapes to fit all patients.



Order Code	Description	Size	Qty
------------	-------------	------	-----

Spacelabs TruLink Technology

Dixtal Cables



700-0029-00	TruLink SpO ₂ cable compatible with Novametrix® sensors	30 cm/12 in.	1 Each
-------------	--	--------------	--------

Dixtal Reusable Sensors



015-0130-01	Dixtal SpO ₂ finger sensor	Adult >30 kg/66 lbs	1 Each
-------------	---------------------------------------	------------------------	--------



015-0133-01	Dixtal SpO ₂ universal sensor	Universal >30 kg/66 lbs or <30 kg/66 lbs	1 Each
-------------	--	--	--------




015-0321-01	Dixtal SpO ₂ universal sensor ear clip	Universal >30 kg/66 lbs or <30 kg/66 lbs	1 Each
-------------	---	--	--------

Accessories

015-0139-00	SpO ₂ universal sensor adhesive taping system finger tape	Adult/Pediatric	100/Box
-------------	--	-----------------	---------

015-0330-00	SpO ₂ foam wrap	Universal	25/Pack
-------------	----------------------------	-----------	---------

For use with Ultraview SL™ Command Module (91496-M).
Cables and connectors are color-coded gray.



Order Code	Description	Size	Qty
LNCS Masimo SET Technology			
LNCS Adapter Cable			
700-0906-01	TruLink SpO ₂ cable compatible with Masimo LNCS sensors (91496-M/91496-U)	305 cm/10 ft	1 Each
LNCS reusable sensors			
690-0247-00	Masimo LNCS DC-I Digit Sensor	Adult > 30 kg	1 Each
690-0248-00	Masimo LNCS DC-IP Pediatric Digit Sensor	Pediatric 10-50 kg	1 Each
690-0249-00	Masimo LNCS TC-I Tip-Clip Ear Sensor	Adult > 30 kg	1 Each
690-0250-00	Masimo LNCS TF-I Transflectance Forehead Sensor	Adult > 30 kg	1 Each

For use with Ultraview SL™ Command Module (91496-M).
Cables and connectors are color-coded gray.



Order Code	Description	Size	Qty
------------	-------------	------	-----

LNCS Masimo SET Technology





LNCS disposable sensors



690-0251-00	Masimo LNCS Actx adhesive sensor, 18 in.	Adult > 30 kg	20/Box
690-0252-00	Masimo LCNS Pctx adhesive sensor, 18 in.	Pediatric 10 - 50 kg	20/Box
690-0258-00	Masimo LNCS Inf-L adhesive sensor, 3 ft.	Infant 3 kg > 20 kg	20/Box
690-0259-00	Masimo LNCS Neo-L adhesive sensor, 3 ft.	Neonatal < 3 kg or > 40 kg	20/Box

For use with Spacelabs Command Module (91496-M) .
Cables and connectors are color-coded gray.



Order Code	Description	Size	Qty
LNOP Masimo SET Technology			
LNOP Adapter Cable			
 700-0789-01	TruLink SpO ₂ cable compatible with LNOP Masimo sensors	244 cm/8 ft.	1 Each
LNOP Reusable Sensors			
 690-0230-00	Masimo LNOP DC-I sensor	Adult > 30 kg	1 Each
 690-0231-00	Masimo LNOP DC-IP sensor	Pediatric 10-50 kg	1 Each
 690-0232-00	Masimo LNOP YI multi-site sensor	Neonatal	1 Each

For use with Spacelabs Command Module (91496-M) .
Cables and connectors are color-coded gray.



Order Code	Description	Size	Qty
------------	-------------	------	-----

LNOP Masimo SET Technology



LNOP Adhesive Disposable Sensors



690-0217-00	Masimo LNOP Adt adhesive sensor	Adult	20/Box
690-0219-00	Masimo LNOP Pdt adhesive sensor	Pediatric	20/Box
690-0221-00	Masimo LNOP Neo-L adhesive sensor	Neonatal	20/Box
690-0218-00	Masimo LNOP Adtx adhesive sensor, transparent tape	Adult	20/Box
690-0226-00	Masimo LNOP Blue thumb/toe sensor	Infant	10/Box

For use with Spacelabs Command Modules (90496,91496-N & 91496-U) and mCare Monitor 91220. Cables and connectors are color-coded purple.



Order Code	Description	Size	Qty
Nellcor OxiMax Technology			
Cable			
	700-0792-00 TruLink SpO ₂ cable compatible with Nellcor OxiMax sensors (91496-N, 91220, 91330)		1 Each
	700-0030-00 TruLink SpO ₂ cable compatible with Nellcor OxiSensor II and OxiMax sensors		1 Each

For use with Spacelabs Command Modules (90496,91496-N & 91496-U) and mCare Monitor 91220. Cables and connectors are color-coded purple.



Order Code	Description	Size	Qty
------------	-------------	------	-----

Nellcor OxiMax Technology

Reusable Sensors

Note: Nellcor sensors are available for purchase through Spacelabs Healthcare only with the purchase of capital equipment. All future purchases should be made directly through Nellcor(Covidien).



690-0003-01	Nellcor OxiMax DS-100A finger sensor	Adult	1 Each
-------------	--------------------------------------	-------	--------



690-0004-01	Nellcor OXI-A/N sensor, reusable	Adult/Neonatal	1 Each
-------------	----------------------------------	----------------	--------



690-0039-01	Nellcor OXI-P/I sensor	Pediatric/Infant	1 Each
-------------	------------------------	------------------	--------



690-0215-01	Nellcor DS-Y multi-site sensor	Universal	1 Each
-------------	--------------------------------	-----------	--------

SUPPLIES & ACCESSORIES TEMPERATURE



For use with Spacelabs Command Modules (90496 & 91496) when monitoring continuous temperature .

Full line of continuous core temperature monitoring products.



Order Code	Description	Size	Qty
------------	-------------	------	-----

Continuous Temperature

Cables

For use with Command Modules 90496/91496, DM3 Monitor 91330 and mCare Monitor 91220



700-0031-00	TruLink Dual Temperature Cable	30 cm/12 in.	1 Each
700-4000-00	Universal temperature probe extension cable, YSI 400 series compatible *Required for use of single use temperature probes/sensors	305 cm/10 ft	1 Each

Reusable Probes/Sensors

Skin Surface



690-0028-01	YSI 700 series compatible with 6.5 mm DIA phone jack type connections	305 cm/10 ft	1 Each
-------------	---	--------------	--------

Esophageal / Rectal



690-0272-00	12 French, YSI 400 series compatible	305 cm/10 ft	1 Each
-------------	--------------------------------------	--------------	--------








118-0256-01	12 French, YSI 700 series compatible	305 cm/10 ft	1 Each
-------------	--------------------------------------	--------------	--------

For use with Spacelabs Command Modules (90496 & 91496) when monitoring continuous temperature .

Full line of continuous core temperature monitoring products.



	Order Code	Description	Size	Qty
Single use Probes/Sensors - Requires YSI 400 series cable (700-4000-00)				
Skin				
	701-4000-00	Skin Temperature Sensor		20/Box
Esophageal/ Stethoscope				
	701-4009-00	9 French, YSI 400 series compatible		20/Box
	701-4012-00	12 French, YSI 400 series		20/Box
	701-4018-00	18 French, YSI 400 series		20/Box
	701-4024-00	24 French, YSI 400 series compatible		20/Box
Esophageal/ Rectal				
	701-4109-00	9 French, YSI 400 series compatible		20/Box
	701-4112-00	12 French, YSI 400 series compatible		20/Box
Foley catheter with temperature probe				
	701-4214-00	14 French, YSI 400 series		20/Box
	701-4216-00	16 French, YSI 400 series		20/Box
	701-4218-00	18 French, YSI 400 series compatible		20/Box
Tympanic temperature sensor				
	701-4300-00	Adult, YSI 400 series compatible		20/Box

SUPPLIES & ACCESSORIES INVASIVE PRESSURE



INVASIVE PRESSURE CABLES



For use with Spacelabs Command Modules (90496 & 91496) configured with the invasive pressure option. Option B supports two pressures; Option C supports four pressures.



(9X496-B) (9X496-C)
Command Modules

Order Code	Description	Size	Qty
------------	-------------	------	-----

Integrated Adapter Cables (connects directly to transducer)

Edward Lifesciences



700-0293-01	Single pressure integrated cable	396 cm/13 ft	1 Each
-------------	----------------------------------	--------------	--------



700-0292-01	Dual pressure integrated cable	396 cm/13 ft	1 Each
-------------	--------------------------------	--------------	--------

Abbott Laboratories/Hospira/ICU Medical



700-0295-01	Single pressure integrated cable	396 cm/13 ft	1 Each
-------------	----------------------------------	--------------	--------



700-0294-01	Dual pressure integrated cable	396 cm/13 ft	1 Each
-------------	--------------------------------	--------------	--------

Interface Cable & Adapters

Interface Cable



690-0021-00	Compatible with Datex-Ohmeda TC-VTK or Becton Dickinson model 684085 pressure transducers. Requires adapter cable 700-0091-00 or 700-0028-00	366 cm/12 ft	1 Each
-------------	--	--------------	--------

Interface Adapter- Transducer interface cable required from transducer company




700-0791-00	Single pressure cable	122 cm/4 ft	1 Each
-------------	-----------------------	-------------	--------



700-0028-00	Dual pressure cable	122 cm/4 ft	1 Each
-------------	---------------------	-------------	--------

Diaphragm Dome






690-0018-00	Diaphragm dome, single use Rotating luer lock fittings Compatible with Datex-Ohmeda TA101OND		12/Box
-------------	--	--	--------

SUPPLIES & ACCESSORIES CARDIAC OUTPUT





For use with Spacelabs Command Modules (90496 & 91496) configured with the cardiac output option C.



Order Code	Description	Size	Qty
Cardiac Output Cables			
700-0027-00	TruLink Cardiac Output Adapter Cable	30 cm/12 in.	1 Each
			
306655-003	Cardiac Output cable compatible with Baxter Edwards/Abbott injectate temperature probes and most models of thermodilution catheters		1 Each
			
306655-004	Cardiac Output cable compatible with Becton Dickinson/Datex-Ohmeda injectate temperature probe		1 Each
			

SUPPLIES & ACCESSORIES RECORDING & PRINTING







High quality Z-fold paper for use with the Bedside Printer Module (90449) and Central Printer (90469).

Roll paper for use with the Standalone printer (91449) and Spacelabs Portable Monitors.



(91449) Printer

Order Code	Description	Size	Qty
Thermal Paper			
Z-fold Paper for use with 90449 and 90469			
 006-0196-00	TruLink Thermal Paper, Z-Fold Pads, 2-channel printing, blank paper, 230 Sheets/pad	27.4 m/90 ft	10 Pads/Box
 307568-001	TruLink Thermal Paper, Z-fold Pads, 2- or 4-channel printing, blank paper, 300 sheets/pad	45.7 m/150 ft	10 Pads/Box
Roll Paper for use with 91449, 91367, 91369, 91370 & 91390			
 305240-001	TruLink Thermal Paper, Rolls, single-channel printing, 40 mm/1.575 inch wide preprinted grid	61 m/200 ft	10 Rolls/Box
 307438-001	TruLink Thermal Paper, Rolls, single- or dual-channel printing, 50 mm/1.97 inch wide blank paper	30.1 m/100 ft	10 Rolls/Box

SUPPLIES & ACCESSORIES BISPECTRAL INDEX



For use with Spacelabs BISx™ Module (91482).



Order Code	Description	Qty
------------	-------------	-----

BISx Cables & Sensors

BISx Cables

719-0004-01	BISx Patient Interface cable	1 Each
719-0006-20	BISx Pod	1 Each
719-0007-01	BISx host monitor cable	1 Each

BISx Sensors

Note: Sales of BISx sensors are limited to one box with the purchase of capital equipment. All future purchases should be made directly through Covidien.

719-0003-10	BISx Extend Sensor, single patient use	25/Box
719-0001-10	BISx Pediatric Sensor, single patient use	25/Box
719-0002-10	BISx Quatro Sensor, single patient use	25/Box

SUPPLIES & ACCESSORIES

ULTRAVIEW DM3 DUAL MODE VITAL SIGNS MONITOR



ULTRAVIEW DM3 VITAL SIGNS MONITOR SUPPLIES & ACCESSORIES



Ultraview DM3 supports measurement of SpO₂, pulse rate, NIBP and temperature, and can be used in any patient care environment where basic vital signs monitoring is needed.



91330

Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

TruLink Nylon Reusable BP Cuffs



Single Tube				
715-1121-10	Single tube, HP	Child	13 - 20 cm	1 Each
715-1131-10	Single tube, HP	Small Adult	18 - 26 cm	1 Each
715-1141-10	Single tube, HP	Adult	26 - 35 cm	1 Each
715-1151-10	Single tube, HP	Adult Long	29 - 38 cm	1 Each
715-1161-10	Single tube, HP	Large Adult	32 - 42 cm	1 Each
715-1171-10	Single tube, HP	Large Adult Long	35 - 44 cm	1 Each
715-1181-10	Single tube, HP	Thigh Cuff	42 - 50 cm	1 Each

TruLink Soft Disposable BP Cuffs



Single Tube				
715-1123-10	Single tube, HP	Child	13 - 20 cm	5/Bag
715-1133-10	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
715-1143-10	Single tube, HP	Adult	26 - 35 cm	5/Bag
715-1153-10	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
715-1163-10	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
715-1173-10	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
715-1183-10	Single tube, HP	Thigh Cuff	42 - 50 cm	5/Bag

TruLink Vinyl Disposable BP Cuffs



Single Tube				
715-1122-10	Single Tube, HP	Child	13 - 20 cm	5/Bag
715-1132-10	Single Tube, HP	Small Adult	18 - 26 cm	5/Bag
715-1142-10	Single Tube, HP	Adult	26 - 35 cm	5/Bag
715-1152-10	Single Tube, HP	Adult Long	29 - 38 cm	5/Bag
715-1162-10	Single Tube, HP	Large Adult	32 - 42 cm	5/Bag
715-1172-10	Single Tube, HP	Large Adult Long	35 - 44 cm	5/Bag
715-1182-10	Single Tube, HP	Thigh Cuff	42 - 50 cm	5/Bag



Single Tube Hose				
714-0018-00	Single tube	Adult	275 cm/9 ft	1 Each
714-0060-00	Single tube	Adult	366 cm/12 ft	1 Each

ULTRAVIEW DM3 VITAL SIGNS MONITOR SUPPLIES & ACCESSORIES



Ultraview DM3 supports measurement of SpO₂, pulse rate, NIBP and temperature, and can be used in any patient care environment where basic vital signs monitoring is needed.



91330

Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

Spacelabs SpO₂ TruLink Technology

Cables



700-0030-00	TruLink SpO ₂ cable		305 cm/10 ft.	1 Each
-------------	--------------------------------	--	---------------	--------

Reusable sensors



015-0660-00	TruLink SpO ₂ finger sensor		Adult >30 kg/66 lbs	1 Each
-------------	--	--	------------------------	--------



015-0661-00	TruLink SpO ₂ Y multi-site sensor Includes: 1 Sensor, 1 ear clip, 10 adult finger/neonate feet/hand wraps and 10 pediatric finger or infant toe adhesive wraps		Universal >30 kg/66 lbs <30 kg/66 lbs	1 Each
-------------	--	--	---	--------

Disposable Sensors



015-0662-00	TruLink SpO ₂ foam sensor		Adult > 30 kg/66 lbs	24/Box
-------------	--------------------------------------	--	-------------------------	--------



015-0663-00	TruLink SpO ₂ vinyl sensor		Adult >30 kg/66 lbs	24/Box
-------------	---------------------------------------	--	------------------------	--------



015-0664-00	TruLink SpO ₂ foam sensor		Pediatric 10-50 kg/22-110 lbs	24/Box
-------------	--------------------------------------	--	----------------------------------	--------



015-0665-00	TruLink SpO ₂ vinyl sensor		Pediatric 10-50 kg/22-110 lbs	24/Box
-------------	---------------------------------------	--	----------------------------------	--------

ULTRAVIEW DM3 VITAL SIGNS MONITOR SUPPLIES & ACCESSORIES



Ultraview DM3 supports measurement of SpO₂, pulse rate, NIBP and temperature, and can be used in any patient care environment where basic vital signs monitoring is needed.



91330

Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

LNCS Masimo SET SpO₂ Technology

LNCS Adapter Cable



700-0906-01	TruLink SpO ₂ cable compatible with Masimo LNCS sensors		305 cm/10 ft	1 Each
-------------	--	--	--------------	--------

LNCS reusable sensors



690-0247-00	Masimo LNCS DC-I Digit Sensor		Adult > 30 kg	1 Each
-------------	-------------------------------	--	------------------	--------



690-0248-00	Masimo LNCS DC-IP Pediatric Digit Sensor		Pediatric 10-50 kg	1 Each
-------------	--	--	-----------------------	--------



690-0249-00	Masimo LNCS TC-I Tip-Clip Ear Sensor		Adult > 30 kg	1 Each
-------------	--------------------------------------	--	------------------	--------



690-0250-00	Masimo LNCS TF-I Transflectance Forehead Sensor		Adult > 30 kg	1 Each
-------------	---	--	------------------	--------

LNCS disposable sensors



690-0251-00	Masimo LNCS Aadx adhesive sensor, 18 in.		Adult > 30 kg	20/Box
-------------	--	--	------------------	--------

ULTRAVIEW DM3 VITAL SIGNS MONITOR SUPPLIES & ACCESSORIES



Ultraview DM3 supports measurement of SpO₂, pulse rate, NIBP and temperature, and can be used in any patient care environment where basic vital signs monitoring is needed.



91330

Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

LNOP Masimo SET SpO₂ Technology

LNOP Adapter Cable



700-0789-01	TruLink SpO ₂ cable compatible with LNOP Masimo sensors		244 cm/8 ft.	1 Each
-------------	--	--	--------------	--------

LNOP Reusable Sensors



690-0230-00	Masimo LNOP DC-I sensor		Adult > 30 kg	1 Each
-------------	-------------------------	--	------------------	--------



690-0231-00	Masimo LNOP DC-IP sensor		Pediatric 10-50 kg	1 Each
-------------	--------------------------	--	-----------------------	--------



690-0232-00	Masimo LNOP YI multi-site sensor		Neonatal	1 Each
-------------	----------------------------------	--	----------	--------

LNOP Adhesive Disposable Sensors



690-0217-00	Masimo LNOP Adt adhesive sensor		Adult	20/Box
-------------	---------------------------------	--	-------	--------



690-0218-00	Masimo LNOP Adtx adhesive sensor, transparent tape		Adult	20/Box
-------------	--	--	-------	--------

ULTRAVIEW DM3 VITAL SIGNS MONITOR SUPPLIES & ACCESSORIES



Ultraview DM3 supports measurement of SpO₂, pulse rate, NIBP and temperature, and can be used in any patient care environment where basic vital signs monitoring is needed.



91330

Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

Nellcor OxiMax SpO₂ Technology

Cable



700-0792-00	TruLink SpO ₂ cable compatible with Nellcor OxiMax sensors			1 Each
-------------	---	--	--	--------

Reusable Sensors

Note: Nellcor sensors are available for purchase through Spacelabs Healthcare only with the purchase of capital equipment. All future purchases should be made directly through Nellcor(Covidien).



690-0003-01	Nellcor OxiMax DS-100A finger sensor		Adult	1 Each
-------------	--------------------------------------	--	-------	--------



690-0004-01	Nellcor OXI-A/N sensor		Adult/Neonatal	1 Each
-------------	------------------------	--	----------------	--------



690-0039-01	Nellcor OXI-P/I sensor		Pediatric/Infant	1 Each
-------------	------------------------	--	------------------	--------



690-0215-01	Nellcor DS-Y multi-site sensor		Universal	1 Each
-------------	--------------------------------	--	-----------	--------

Temperature

Probes

019-0539-00	Replacement Temperature Probe, Oral/Axillary			1 Each
701-0122-00	Disposable Temperature Probe Covers			200/Box

SUPPLIES & ACCESSORIES

AMBULATORY BLOOD PRESSURE MONITOR



Ambulatory blood pressure is rapidly becoming a leading diagnostic tool. It compensates for White Coat Syndrome and allows for a series of measurements rather than relying on one isolated reading.



(90207)



(90217)

Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

Reusable ABP Cuffs



UltraCheck Reusable Nylon Cuffs

016-1320-01	Self-applied, reusable	child	13 - 20 cm	1 Each
016-1826-01	Self-applied, reusable	small adult	18 - 26 cm	1 Each
016-2635-01	Self-applied, reusable	adult	26 - 35 cm	1 Each
016-3242-01	Self-applied, reusable	large adult	32 - 42 cm	1 Each
016-3854-01	Self-applied, reusable	XL adult	38 - 54 cm	1 Each

Spacelabs Reusable Cloth Cuffs

015-0118-04Q	Self-applied with quick connect rectus fitting	child	12 - 20 cm	1 Each
015-0067-04Q	Self-applied with quick connect rectus fitting	small adult	17 - 26 cm	1 Each
015-0068-05Q	Self-applied with quick connect rectus fitting	adult	24 - 32 cm	1 Each
016-0077-05Q	Self-applied with quick connect rectus fitting	large adult	32 - 42 cm	1 Each
016-0109-04Q	Self-applied with quick connect rectus fitting	XL adult	38 - 50 cm	1 Each

Disposable Cuff Covers

For use with Reusable ABP Cuffs

706-0182-00	Universal cuff cover	Standard		50/Box
706-0183-00	Universal cuff cover	Large		50/Box

Ambulatory blood pressure is rapidly becoming a leading diagnostic tool. It compensates for White Coat Syndrome and allows for a series of measurements rather than relying on one isolated reading.



(90207)



(90217)

Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

Disposable ABP Cuffs



Disposable Cuffs

015-1826-00	Single tube, no fitting	Small Adult	18 - 26 cm	20/Box
015-2635-00	Single tube, no fitting	Adult	26 - 35 cm	20/Box
015-3242-00	Single tube, no fitting	Large Adult	32 - 42 cm	20/Box
015-4250-00	Single tube, no fitting	X-large Adult	42 - 50 cm	20/Box
712-1039-00	Quick Release Metal Fittings, female bayonet. For use with above single tube, disposable cuffs with no fitting.			3/Pack

Ambulatory blood pressure is rapidly becoming a leading diagnostic tool. It compensates for White Coat Syndrome and allows for a series of measurements rather than relying on one isolated reading.



(90207)



(90217)

Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

ABP Accessories

Pouches - For use with APB Monitor 90207



016-0263-00	Vinyl carrying pouch with Velcro closure			Each
-------------	--	--	--	------



016-0410-00	Cloth carrying pouch with Velcro closure			Each
-------------	--	--	--	------

Pouches - For use with APB Monitor 90217



016-0340-00	Carrying pouch, cloth			Each
-------------	-----------------------	--	--	------

Pouches - For use with APB Monitor 90227

706-0200-00	Disposable Pouch			50/Pack
-------------	------------------	--	--	---------



015-0715-00	Reusable Pouch			Each
-------------	----------------	--	--	------



015-0716-00	Reusable Shoulder Strap (for use with Pouch 015-0715-00)			Each
-------------	--	--	--	------

Ambulatory blood pressure is rapidly becoming a leading diagnostic tool. It compensates for White Coat Syndrome and allows for a series of measurements rather than relying on one isolated reading.



(90207)



(90217)

Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

ABP Accessories

Support



015-0070-00	Cuff support harness			Each
-------------	----------------------	--	--	------



016-0080-00	Nylon web adjustable waist belt, d-ring adjustable fastener			Each
-------------	---	--	--	------



016-0262-00	Nylon web adjustable shoulder strap, Velcro fasteners For carrying the pouch and monitor over the shoulder			Each
-------------	---	--	--	------



344-0056-00	Garment clip with vinyl strap, 3.5"			Each
-------------	-------------------------------------	--	--	------

Patient Diaries and Paper




000-0027-02	ABP Patient Diary, English, 24 hour Monitoring			Each
-------------	--	--	--	------



006-0210-02	ABP Thermal paper, 57 mm x 32 ft.			6 Rolls/ Pack
-------------	-----------------------------------	--	--	---------------

SUPPLIES & ACCESSORIES

HOLTER



For use with Spacelab's Digital Holter Recorders
Lifecard CF, evo & Aria.



Lifecard CF



evo



Aria

Order Code	Description	For Use With	Qty
Holter Hook-up Kits			
Aria Hook-Up Kits			
23189-106-401	Aria hook-up kit with case, pouch and clear tape	Aria	20/Case
23189-106	4 electrode hook-up kit	Aria, evo	Each
23189-110-401	Hook-up Kit (no pouch)	Aria, evo	20/Case
50-0048	4 electrode hook-up kit	evo	
50-0049	4 electrode hook-up kit with disposable chest pouch	evo	
50-0190	4 electrode hook-up kit	evo	20/Case
50-0191	4 electrode hook-up kit with chest pouch	evo	20/Case
23189-116	LCF hook-up kit	Lifecard CF	Each
23189-116-401	Lifecard hook-up kit	Lifecard CF	20/Case
23189-111-401	Holter hook-up Kit with case, 3-channel	485	20/Case

For use with Spacelab's Digital Holter Recorders
Lifecard CF, evo & Aria.



Lifecard CF



evo



Aria

Order Code	Description	For Use With	Qty
Lead Wires & Electrodes			
Lead Wires			
25452-525	4-Lead wire set, snap, three 14" leads and one 7" white lead	Aria, evo	4LW/set
25452-535	4-Lead wire set, snap, three 14" leads and one 14" white lead	Aria, evo	4LW/set
25452-529	4-Lead wire set, snap, 40" leads, Aria and evo belt clip configuration	Aria, evo	4LW/set
25452-513	5-Lead wire set, snap, 18" (white/black/red/green/brown)	420, 456	5 LW/ set
25452-519	7-Lead wire set, snap, 18" (white/black/red/green/orange/blue/brown)	423, 423A, 459, 461, 461A, 464, 483, 485	7 LW/set
19-7509	6-Lead wire set, snap, 24 in.	Lifecard CF	6LW/Set
19-7510	7-lead wire set, 24 in.	Lifecard CF	7LW/Set
25450-187	Individual color lead wire, snap, white, 14 in.	Aria, evo	



For use with Spacelab's Digital Holter Recorders
Lifecard CF, evo & Aria.



Lifecard CF



evo




Aria

Order Code	Description	For Use With	Qty
------------	-------------	--------------	-----

Lead Wires & Electrodes

Electrodes



27565-9	ECG Electrode, clear tape	Holter	4/Pack
27565-9CASE	ECG Electrode, clear tape	Holter	150/Pack
717-0071-00	ECG Electrode, AMBU	Holter	5/Pack



27565-5	ECG Electrode, foam	Holter	5/Pack
---------	---------------------	--------	--------



27565-7	ECG Electrode, foam	Holter	7/Pack
---------	---------------------	--------	--------

For use with Spacelab's Digital Holter Recorders
Lifecard CF, evo & Aria.



Lifecard CF



evo



Aria

Order Code	Description	For Use With	Qty
Patient Cables, Pouches, Shoulder Straps & Belts			
Patient Cables			
700-0316-00	Patient cable, 3-Channel, 6-lead, 22 in.	Lifecard CF	Each
700-0317-00	Patient cable with Belt Clip, 3-Channel, 3-lead, 41 in.	Lifecard CF	Each
700-0318-00	Patient cable, 3-channel	Lifecard CF	Each
700-0319-00	Patient cable with Neck Pouch, 3-Channel, 3-lead, 26 in.	Lifecard CF	Each
700-0320-00	Patient cable, 2-Channel	Lifecard CF	Each
700-0321-00	Patient cable, 2-Channel, short	Lifecard CF	Each
700-0322-00	Patient cable, 3-Channel, 7-lead, 23 in.	Lifecard CF	Each
700-0323-00	Patient cable with Beth Clip, 2-Channel, 4-lead, long	Lifecard CF	Each
700-0324-00	Patient cable with Neck Pouch, 2-Channel, 4-lead, short	Lifecard CF	Each
26927-503	Patient Cable, 30 in.	485	



For use with Spacelab's Digital Holter Recorders
Lifecard CF, evo & Aria.



Lifecard CF



evo



Aria

Order Code	Description	For Use With	Qty
Patient Cables, Pouches, Shoulder Straps & Belts			
Pouches, Shoulder Straps & Belts			
27249-7-401	Aria disposable pouch	Aria	50/Pack
27380	Aria belt clip (replacement)	Aria	Each
50-0042	evo holder for waist belt	evo	Each
PK0001	Fabric waist belt	evo	Each
50-0050	Disposable chest pouch for evo recorder	evo	50/Pack
706-0189-00	Disposable pouch for evo recorder	evo	100/Pack
46-0408	Fabric pouch set, 2 pouches with laces, cordlock & instructions	Lifecard CF	Each
19-7130	Laces for pouch set	Lifecard CF	Each
19-7131	cordlock for pouch set	Lifecard CF	Each
46-0092	Fabric pouch	Lifecard CF	Each
46-0512	Disposable pouch	Lifecard CF	100/Pack
27685-401	Disposable pouch	LC12	50/Pack
200-1508-001	Pouch & straps	485	Each
24156	Shoulder strap	485	Each



For use with Spacelab's Digital Holter Recorders
Lifecard CF, evo & Aria.



Lifecard CF



evo



Aria

Order Code	Description	For Use With	Qty
------------	-------------	--------------	-----

Recording

Recorder Memory			
CFC90MB	Flashcard Memory, 90MByte	Lifecard CF	Each
CFC256MB	Flashcard LifeCard CF, 256MB	LC12	Each
46-0063	Plastic memory card cases	CFCxxMB	Each
50-0220	LCF / evo pen	LCF / evo	Each
19-7118	Lifecard Compact USB Flashcard Reader	CFCxxMB	Each
27311	Aria download cable	Aria	Each
46-0413	LCF opener tool	LCF	Each
26633	Banana Universal test cable	483	Each
26918-107	Memory card	485	Each

Patient Diaries

Patient Diaries			
073-0267-00	Patient Diary, English CE0123	Aria, evo, Lifecard CF	Each
46-0459	Patient Diary, 7 day, English	Lifecard CF	Each
46-0573	Patient Diary, English, Italian, French, Spanish, Polish	Lifecard CF	Each
46-0435	Patient diary, 7 day, German	Lifecard CF	Each

SUPPLIES & ACCESSORIES GENERIC & SUNTECH





Order Code	Description	For Use With	Qty
Generic Hook-up Kits			
DC Hook-up			
30404	Liquid Medisol adhesive remover, 4 oz.	all electrodes	Each
13-9083	Razor with blade cover	Holter / Event	Each
392003-001	Skin prep pad		100/Box
Stress Test			
91-0013-00	Cable, Tango, RS-232	Suntech Tango	Each
98-0002-50	Microphone pads, 5 PK	Suntech Tango	5/Pack
98-0006-00	K-Sound microphone	Suntech Tango	Each
98-0061-01	Cuff, Orbit-K, small adult	Suntech Tango	Each
Recorder Memory			
41-0080	CD12USB stress tests hook-up kit	CD12-S	Each

SUPPLIES & ACCESSORIES

12-LEAD ECG





Order Code	Description	For Use With	Qty
------------	-------------	--------------	-----

12-Lead ECG Patient Cables

AHA Color Coding

MS1R-107048-A0	ECG cable with banana connector, 4MM	CardioExpress	Each
MS1R-107584-A0	ECG cable with grabber connector, 4MM	CardioExpress	Each
25-0028	10 Electrode patient cable, short	CD12	Each
25-0020	10 Electrode patient cable, 4MM	CD12	Each
010-1885-00	ECG cable, 4MM, banana connector	CD12	Each

12 Lead Hook-up Sets

AHA Color Coding

41-0053US	12-lead hook-up set, AHA	Generic	
MS1R-107371-A0	Adult Precordial Suction Electrodes	CardioExpress	6 Pcs/ Set
MS1R-32461	Child Precordial Suction Electrodes	CardioExpress	6 Pcs/ Set
MS1R-32462	Child Limb Clamp Electrodes	CardioExpress	4 Pcs/ Set
MS1R-32463	Banana Socket Adapters, clip/snap	CardioExpress	10 Pieces

Miscellaneous 12-Lead ECG Accessories

Electrodes

RT-41	Skintact Tab Electrodes	Generic	100/Pack
0315BOX	Nikotab Electrodes	Generic	500/Box
717-0072-00	Adult Foam Electrodes w/snap	Generic	25/Pack

SUPPLIES & ACCESSORIES

EVENT MONITORING





Order Code	Description	For Use With	Qty
Event Monitoring Supplies & Accessories			
Hook-Up Kits			
23-0255	Lanyard clip & assembly	CardioCall	Each
Lead Wires & Electrodes			
23-0059	Patient cable, 3 electrodes, 2	CardioCall	Each
23-0051	Patient cable, 2 electrodes, 1	CardioCall	Each
Programming & Download			
23-0245	Programing card & key	CCVS20	Each
23-0252	Programing card & key	CCST80	Each
23-0026	Download cable	CardioCall	Each

SUPPLIES & ACCESSORIES

ANESTHESIA



ANESTHESIA DELIVERY SYSTEMS SUPPLIES & ACCESSORIES



For use with Spacelabs ARKON Anesthesia Delivery System and BleaseSirius and BleaseFocus Anesthesia Machines.



BleaseFocus

BleaseSirius

ARKON

Order Code	Description	Qty
------------	-------------	-----

Anesthesia Supplies

Nomoline Gas Sampling

For use with ARKON Anesthesia Delivery System only



103-0234-00	Nomoline adapter with female Luer lock, multi-patient use	25/Box
166-7085-00	Sidestream sampling line with male Luer lock on both ends	2m/6.5 ft 25/Box
166-7121-00	Sidestream sampling line with male Luer lock on both ends	3m.10 ft 50/Box

CO₂ Canisters

For use with ARKON Anesthesia Delivery System only



015-0702-00	CO2 Canister w/Spherasorb, white to violet	6/Case
015-0701-00	CO2 Canister w/LoFloSorb, green to violet	6/Case

Patient Flow Sensors

For use with BleaseSirius and BleaseFocus Anesthesia Machines only



10110090	Adult Patient flow sensor	1 Each
54100063	Flow sensor caps	50/Box
376-0561-01	Pediatric Patient flow sensor	1 Each
10110219	1.5M Patient flow sensor tubing assembly	1 Each
10110218	2M Patient flow sensor tubing assembly	1 Each
10110122	2.5M Patient flow sensor tubing assembly	1 Each
10110102	3M Patient flow sensor tubing assembly	1 Each
10110103	6M Patient flow sensor tubing assembly	1 Each
10110104	1.8M Patient flow sensor tubing assembly kit, single ended	1 Each
10110106	1.8M Patient flow sensor & tubing assembly kit, single ended	1 Each
70300001	O2 probe with tee	1 Each

ANESTHESIA DELIVERY SYSTEMS SUPPLIES & ACCESSORIES



For use with Spacelabs ARKON Anesthesia Delivery System and BleaseSirius and BleaseSirius Anesthesia Machines.



BleaseFocus

BleaseSirius

ARKON

Order Code	Description	Qty
------------	-------------	-----

Anesthesia Supplies - not anesthesia machine specific

Reusable Breathing Circuits

M1050-4	Reusable patient "Y" circuit, adult	1 Each
ST7557	Magill breathing circuit w/face mask	1 Each
ST2111	2.0 L rebreathing bag, 22mm neck	1 Each
ST2670	Malleable face mask, size 4	1 Each
2100	Jackson-Rees pediatric circuit	1 Each
2150	Rendell Baker pediatric mask set, 4 sizes	4/Set
42100	Malleable face mask set, 7 sizes	7/Set

Disposable Breathing Circuits

54100035	Pediatric Y piece breathing circuit 2010-022	1 Each
54100036	Infant T piece mapleson f breathing circuit 2120	1 Each
54100031	1.8M Magill-Mapleson, open circuit	1 Each
54100032	Individual deluxe Bain coaxial circuit 2115	1 Each
54100033	1.6M Basic Y piece breathing circuit 200	1 Each
10400001	Basic Y piece breathing circuit w/limb and bag 2010	1 Each

Disposable Face Masks

54100040	Infant disposable face mask	1 Each
54100039	Small disposable face mask	1 Each
54100038	Medium disposable face mask	1 Each
54100037	Large disposable face mask	1 Each

Disposable Reservoir Bags

54100041	0.5 L anesthetic reservoir bag, 15mm female neck	1 Each
54100042	1.0 L anesthetic reservoir bag, 22mm female neck	1 Each
54100043	2.0 L anesthetic reservoir bag, 22mm female neck	1 Each
54100044	3.0 L anesthetic reservoir bag, 22mm female neck	1 Each

Connectors & Filters

54100026	22 mm connector - 30f scavenging 1971	1 Each
54100027	Filta-Guard 1944 standard bacterial/viral breathing filter,	1 Each

Order Code	Description	UOM	Page(s)
1103	Unifusor Plus 500cc w/trigger gauge inflation system	ea	40
1104	Unifusor Plus 1000cc w/trigger gauge inflation system	ea	40
1105	Unifusor Plus 3000cc w/trigger trigger gauge inflation system	ea	40
1199	Unifusor Plus trigger gauge, assembly only	ea	40
2100	Jackson-Rees pediatric circuit	ea	103
2150	Rendell Baker pediatric mask set, 4 sizes	ea	103
24156	Shoulder strap	ea	96
26633	Banana Universal test cable	ea	98
27311	Aria download cable	ea	98
27380	Aria belt replacement clip	ea	97
30404	Liquid Medisol adhesive remover, 4 oz.	ea	99
42100	Malleable face mask set, 7 sizes	ea	103
10110090	Adult Patient flow sensor	ea	102
10110102	3M Patient flow sensor tubing assembly	ea	102
10110103	6M Patient flow sensor tubing assembly	ea	102
10110104	1.8M Patient flow sensor tubing assembly kit, single ended	ea	102
10110106	1.8M Patient flow sensor & tubing assembly kit, single ended	ea	102
10110122	2.5M Patient flow sensor tubing assembly	ea	102
10110218	2M Patient flow sensor tubing assembly	ea	102
10110219	1.5M Patient flow sensor tubing assembly	ea	102
10400001	Basic Y piece breathing circuit w/limb and bag 2010	ea	103
54100026	22 mm connector - 30f scavenging 1971	ea	103
54100027	Filter-Guard 1944 standard bacterial/viral breathing filter, 22mm OD x 22mm ID	ea	103
54100031	1.8M Magill-Mapleson, open circuit	ea	103
54100032	Individual deluxe Bain coaxial circuit 2115	ea	103
54100033	1.6M Basic Y piece breathing circuit 200	ea	103
54100035	Pediatric Y piece breathing circuit 2010-022	ea	103
54100036	Infant T piece mapleson f breathing circuit 2120	ea	103
54100037	Large disposable face mask	ea	103
54100038	Medium disposable face mask	ea	103
54100039	small disposable face mask	ea	103
54100040	Infant disposable face mask	ea	103
54100041	0.5 L Anesthetic Reservoir Bag, 15mm female neck	ea	103
54100042	1.0 L Anesthetic Reservoir Bag, 22mm female neck	ea	103
54100043	2.0 L Anesthetic Reservoir Bag, 22mm female neck	ea	103
54100044	3.0 L Anesthetic Reservoir Bag, 22mm female neck	ea	103
54100063	Flow sensor caps	50/pk	102
70300001	O2 probe with tee	ea	102
000-0027-02	ABP English Patient Diary, 24 hour Monitoring	ea	92
006-0196-00	Z-fold thermal paper, 2- channel printing	10/bx	82
006-0210-02	ABP Thermal paper, 57 mm x 32 ft	6/bx	92
010-1885-00	ECG cable, 4MM, banana connector, AHA	ea	100
010-1980-00	CO2 mainstream sensor, 92517	ea	62
012-0283-00	ECG Lead Wire Set, TruLink, DIN unshielded, 3L, snap, AAMI, 46 cm/18 in	1 Set	52
012-0283-02	ECG Lead Wire Set, TruLink, DIN unshielded, 3L, snap, AAMI, 102 cm/40 in	1 Set	52
012-0283-11	ECG Lead Wire Set, TruLink, DIN unshielded, 3L, snap, AAMI, 61 cm/24 in	1 Set	52
012-0284-00	ECG Lead Wire Set, TruLink, DIN unshielded, 3L pinch, AAMI, 46 cm/18 in	1 Set	52
012-0284-02	ECG Lead Wire Set, TruLink, DIN unshielded, 3L pinch, AAMI, 102 cm/40 in	1 Set	52
012-0284-11	ECG Lead Wire Set, TruLink, DIN unshielded, 3L pinch, AAMI, 61 cm/24 in	1 Set	52
012-0285-00	ECG Lead Wire Set, TruLink, DIN unshielded, 5L, snap, AAMI, 46 cm/18 in	1 Set	52
012-0285-02	ECG Lead Wire Set, TruLink, DIN unshielded, 5L, snap, AAMI, 102 cm/40 in	1 Set	52
012-0285-11	ECG Lead Wire Set, TruLink, DIN unshielded, 5L, snap, AAMI, 61 cm/24 in	1 Set	52, 56
012-0286-00	ECG Lead Wire Set, TruLink, DIN unshielded, 5L, pinch, AAMI, 46 cm/18 in	1 Set	52

Order Code	Description	UOM	Page(s)
012-0286-02	ECG Lead Wire Set, TruLink, DIN unshielded, 5L, pinch, AAMI, 102 cm/40 in	1 Set	52
012-0286-11	ECG Lead Wire Set, TruLink, DIN unshielded, 5L, pinch, AAMI, 61 cm/24 in	1 Set	52, 56
012-0295-00	ECG Lead Wires, TruLink, DIN unshielded, snap, AAMI, 46 cm/18 in, red	3/pk	53
012-0295-02	ECG Lead Wires, TruLink, DIN unshielded, snap, AAMI, 102 cm/40 in, red	3/pk	53
012-0295-06	ECG Lead Wires, TruLink, DIN unshielded, snap, AAMI, 61 cm/24 in, red	3/pk	53, 56
012-0295-10	ECG Lead Wires, TruLink, DIN unshielded, snap, AAMI, 46 cm/18 in, white	3/pk	53
012-0295-12	ECG Lead Wires, TruLink, DIN unshielded, snap, AAMI, 102 cm/40 in, white	3/pk	53
012-0295-16	ECG Lead Wires, TruLink, DIN unshielded, snap, AAMI, 61 cm/24 in, white	3/pk	53, 56
012-0295-20	ECG Lead Wires, TruLink, DIN unshielded, snap, AAMI, 46 cm/18 in, black	3/pk	53
012-0295-22	ECG Lead Wires, TruLink, DIN unshielded, snap, AAMI, 102 cm/40 in, black	3/pk	53
012-0295-26	ECG Lead Wires, TruLink, DIN unshielded, snap, AAMI, 61 cm/24 in, black	3/pk	53, 56
012-0295-30	ECG Lead Wires, TruLink, DIN unshielded, snap, AAMI, 46 cm/18 in, green	3/pk	53
012-0295-32	ECG Lead Wires, TruLink, DIN unshielded, snap, AAMI, 102 cm/40 in, green	3/pk	53
012-0295-36	ECG Lead Wires, TruLink, DIN unshielded, snap, AAMI, 61 cm/24 in, green	3/pk	53, 56
012-0295-40	ECG Lead Wires, TruLink, DIN unshielded, snap, AAMI, 46 cm/18 in, brown	3/pk	53
012-0295-42	ECG Lead Wires, TruLink, DIN unshielded, snap, AAMI, 102 cm/40 in, brown	3/pk	53
012-0295-46	ECG Lead Wires, TruLink, DIN unshielded, snap, AAMI, 61 cm/24 in, brown	3/pk	53, 56
012-0296-00	ECG Lead Wires, TruLink, DIN unshielded, pinch, AAMI, 46 cm/18 in, red	3/pk	54
012-0296-02	ECG Lead Wires, TruLink, DIN unshielded, pinch, AAMI, 102 cm/40 in, red	3/pk	54
012-0296-04	ECG Lead Wires, TruLink, DIN unshielded, pinch, AAMI, 61 cm/24 in, red	3/pk	54, 57
012-0296-10	ECG Lead Wires, TruLink, DIN unshielded, pinch, AAMI, 46 cm/18 in, white	3/pk	54
012-0296-12	ECG Lead Wires, TruLink, DIN unshielded, pinch, AAMI, 102 cm/40 in, white	3/pk	54
012-0296-14	ECG Lead Wires, TruLink, DIN unshielded, pinch, AAMI, 61 cm/24 in, white	3/pk	54, 57
012-0296-20	ECG Lead Wires, TruLink, DIN unshielded, pinch, AAMI, 46 cm/18 in, black	3/pk	54
012-0296-22	ECG Lead Wires, TruLink, DIN unshielded, pinch, AAMI, 102 cm/40 in, black	3/pk	54
012-0296-24	ECG Lead Wires, TruLink, DIN unshielded, pinch, AAMI, 61 cm/24 in, black	3/pk	54, 57
012-0296-32	ECG Lead Wires, TruLink, DIN unshielded, pinch, AAMI, 102 cm/40 in, green	3/pk	54
012-0296-34	ECG Lead Wires, TruLink, DIN unshielded, pinch, AAMI, 61 cm/24 in, green	3/pk	54, 57
012-0296-42	ECG Lead Wires, TruLink, DIN unshielded, pinch, AAMI, 102 cm/40 in, brown	3/pk	55
012-0296-44	ECG Lead Wires, TruLink, DIN unshielded, pinch, AAMI, 61 cm/24 in, brown	3/pk	54, 57
012-0498-00	ECG Lead Wires, TruLink, shielded, snap, AAMI, 61 cm/24 in, red	3/pk	49
012-0498-01	ECG Lead Wires, TruLink, shielded, snap, AAMI, 61 cm/24 in, green	3/pk	49
012-0498-03	ECG Lead Wires, TruLink, shielded, snap, AAMI, 61 cm/24 in, white	3/pk	49
012-0498-04	ECG Lead Wires, TruLink, shielded, snap, AAMI, 61 cm/24 in, black	3/pk	49
012-0498-05	ECG Lead Wires, TruLink, shielded, snap, AAMI, 61 cm/24 in, brown	3/pk	49
012-0498-06	ECG Lead Wires, TruLink, shielded, snap, AAMI, 102 cm/40 in, red	3/pk	49
012-0498-07	ECG Lead Wires, TruLink, shielded, snap, AAMI, 102 cm/40 in, green	3/pk	49
012-0498-08	ECG Lead Wires, TruLink, shielded, snap, AAMI, 102 cm/40 in, white	3/pk	49
012-0498-09	ECG Lead Wires, TruLink, shielded, snap, AAMI, 102 cm/40 in, black	3/pk	49
012-0498-11	ECG Lead Wires, TruLink, shielded, snap, AAMI, 61 cm/24 in, yellow/brown	3/pk	49
012-0498-12	ECG Lead Wires, TruLink, shielded, snap, AAMI, 61 cm/24 in, green/brown	3/pk	49
012-0498-13	ECG Lead Wires, TruLink, shielded, snap, AAMI, 61 cm/24 in, blue/brown	3/pk	49
012-0498-14	ECG Lead Wires, TruLink, shielded, snap, AAMI, 61 cm/24 in, orange/brown	3/pk	49
012-0498-15	ECG Lead Wires, TruLink, shielded, snap, AAMI, 61 cm/24 in, violet/brown	3/pk	49
012-0498-27	ECG Lead Wires, TruLink, shielded, snap, AAMI, 102 cm/40 in, brown	3/pk	49
012-0498-34	ECG Lead Wires, TruLink, shielded, snap, AAMI, 61 cm/24 in, black	5/pk	49
012-0498-39	ECG Lead Wires, TruLink, shielded, snap, AAMI, 102 cm/40 in, black	5/pk	49
012-0605-10	ECG Lead Wire Set, TruLink, DIN unshielded, 5L, snap, 64 cm/25.2 in	1 Set	56
015-0067-04Q	ABP Cuff, Spacelabs Cloth, self-applied, sm adult, 17-26 cm	ea	89
015-0068-05Q	ABP Cuff, Spacelabs Cloth, self-applied, adult, 24-32 cm	ea	89
015-0070-00	Cuff support harness	ea	92
015-0118-04Q	ABP Cuff, Spacelabs Cloth, self-applied, child, 12-20 cm	ea	89
015-0130-01	Dixtal SPO2 finger sensor, reusable, adult	ea	71

Order Code	Description	UOM	Page(s)
015-0133-01	Dixtal SPO2 universal sensor, reusable	ea	71
015-0139-00	SPO2 universal sensor adhesive taping system finger tape, single use	100/bx	71
015-0312-01	Anesthetic agent sampling line + gas sampling tee, 90515/91518	50/bx	66
015-0313-01	Anesthetic agent sampling line + gas sampling tee, 90515/91518	50/bx	66
015-0321-01	Dixtal SPO2 universal sensor ear clip, reusable	10/pk	71
015-0330-00	SPO2 universal foam wrap	25/pk	71
015-0500-02	Telemetry pouch, single use	100/bx	58
015-0660-00	SPO2 finger sensor, adult, reusable	ea	59, 69, 85
015-0661-00	SPO2 Y multi-site sensor, universal, reusable	ea	59, 69, 85
015-0662-00	SPO2 foam sensor, adult, single use	24/bx	59, 70, 85
015-0663-00	SPO2 vinyl sensor, adult, single use	24/bx	59, 70, 85
015-0664-00	SPO2 foam sensor, pediatric, single use	24/bx	59, 70, 85
015-0665-00	SPO2 vinyl sensor, pediatric single use	24/bx	59, 70, 85
015-0666-00	SPO2 finger sensor, adult, single use	24/bx	59, 70, 85
015-0667-00	SPO2 cloth sensor, neonate, single use	24/bx	59, 70, 85
015-0683-00	CO2 sidestream Nomoline adapter w/sampling line, single use	25/bx	61
015-0701-00	CO2 canister w/LoFloSorb, green to violet	6/cs	102
015-0702-00	CO2 canister w/Spherasorb, white to violet	6/cs	102
015-0715-00	Reusable Pouch, 90227	ea	91
015-0716-00	Reusable Shoulder Strap	ea	91
015-1826-00	ABP Cuff, Spacelabs Soft, self-applied, sm adult, 18-26 cm, no fitting	20/bx	90
015-2635-00	ABP Cuff, Spacelabs Soft, self-applied, adult, 26-35 cm, no fitting	20/bx	90
015-3242-00	ABP Cuff, Spacelabs Soft, self-applied, lg adult, 32-42 cm, no fitting	20/bx	90
015-4250-00	ABP Cuff, Spacelabs Soft, self-applied, xl adult, 42-50 cm, no fitting	20/bx	90
016-0077-05Q	ABP Cuff, Spacelabs Cloth, self-applied, lg adult, 24-32 cm	ea	89
016-0080-00	Nylon web adjustable waist belt, D-ring adjustable fastener	ea	92
016-0109-04Q	ABP Cuff, Spacelabs Cloth, self-applied, xl adult, 38-50 cm	ea	89
016-0262-00	Nylon web adjustable shoulder strap w/velcro fasteners	ea	92
016-0263-00	Vinyl carrying pouch wth Velcro closure, 90207	ea	91
016-0340-00	Carrying pouch, cloth, 90217	ea	91
016-0410-00	Cloth carrying pouch with Velcro closure, 90207	ea	91
016-0428-00	Filter cartridge	ea	68
016-1320-01	ABP Cuff, UltraCheck Nylon, self-applied, child, 13 - 20 cm	ea	89
016-1826-01	ABP Cuff, UltraCheck Nylon, self-applied, small adult, 18 - 26 cm	ea	89
016-2635-01	ABP Cuff, UltraCheck Nylon, self-applied, adult, 26 - 35 cm	ea	89
016-3242-01	ABP Cuff, UltraCheck Nylon, self-applied, large adult, 32 - 42 cm	ea	89
016-3854-01	ABP Cuff, UltraCheck Nylon, self-applied, XL adult, 38 - 54 cm	ea	89
019-0539-00	Oral/axillary temperature probe, 91330	ea	88
0315BOX	Nikotab Eelectrode	500/bx	100
050-1120-00	O2 sensor Kit	ea	68
073-0267-00	Patient Diary, English, (Aria, EVO, LCF) CE0123	ea	98
103-0059-00	TruLink ECG cable adapter, 17-pin to 6-pin adapter	ea	55
103-0234-00	CO2 sidestream Nomoline adapter, mulit-use	25/bx	61, 102
1103FGA-24	Unifusor II 500cc w/aneroid gauge and thumbwheel valve, semi-reusable	24/bx	39
1103SGA-24	Unifusor II 500cc w/aneroid gauge and stopcock valve, semi-reusable	24/bx	39
1103TGA-24	Unifusor II 500cc w/piston gauge and thumbwheel valve, semi-reusable	24/bx	39
1103X	Unifusor Plus 500cc, cuff only	ea	40
1103X-05	Unifusor Plus 500cc, cuff only	5/bx	40
1103XB-05	Unifusor II 500cc, cuff only	5/bx	39
1104FGA-24	Unifusor II 1000cc w/aneroid gauge and thumbwheel valve, semi-reusable	24/bx	39
1104SGA-24	Unifusor II 1000cc w/aneroid gauge and stopcock valve, semi-reusable	24/bx	39
1104TGA-24	Unifusor II 1000cc w/piston gauge and thumbwheel valve, semi-reusable	24/bx	39
1104X	Unifusor Plus 1000cc, cuff only	ea	40

Order Code	Description	UOM	Page(s)
1104X-05	Unifusor Plus 1000cc, cuff only	5/bx	40
1104XB-05	Unifusor Plus 1000cc, cuff only	5/bx	39
1105FGA-10	Unifusor II 3000cc w/aneroid gauge and thumbwheel valve, semi-reusable	10/bx	39
1105SGA-10	Unifusor II 3000cc w/aneroid gauge and stopcock valve, semi-reusable	10/bx	39
1105TGA-10	Unifusor II 3000cc w/piston gauge and thumbwheel valve, semi-reusable	10/bx	39
1105X	Unifusor Plus 3000cc, cuff only	ea	40
1105X-05	Unifusor Plus 3000cc, cuff only	5/bx	40
1105XB-05	Unifusor II 3000cc, cuff only	5/bx	39
1113SGA-24	Unifusor II 500cc w/piston gauge and stopcock valve, semi-reusable	24/bx	39
1114SGA-24	Unifusor II 1000cc w/piston gauge and stopcock valve, semi-reusable	24/bx	39
1115SGA-10	Unifusor II 3000cc w/piston gauge and stopcock valve, semi-reusable	10/bx	39
118-0256-01	Esophageal/rectal probe, reusable	ea	78
13-9083	Razor with blade cover	ea	99
162-0046-01	Anesthetic agent sampling line, 90515/91518	50/bx	66
162-0047-01	Anesthetic agent sampling line, 90515/91518	50/bx	66
166-7085-00	CO2 sampling line, sidestream, 2m/6.5 ft	25/bx	62, 102
166-7121-00	CO2 sampling line, sidestream, 3m/10 ft	50/bx	62, 102
19-7118	Lifecard Compact USB Flashcard Reader	ea	98
19-7130	Laces for pouch set	ea	97
19-7131	Cordlock for pouch set	ea	97
19-7509	6-Lead wire set, snap, 61 cm/24 in	ea	94
19-7510	7-lead wire set, 61 cm/24 in	ea	94
200-1508-001	Pouch & straps for 485	ea	97
23-0026	Download cable	ea	101
23-0051	Patient cable, 2 electrodes, 1 ch	ea	101
23-0059	Patient cable, 3 electrodes, 2 ch	ea	101
23-0245	Programing card & key, CCVS20	ea	101
23-0252	Programing card & key, CCST80	ea	101
23-0255	Lanyard clip & assembly	ea	101
23189-106	4 electrode hook-up kit	ea	93
23189-106-401	Aria hook-up kit with case, pouch and clear tape electrodes	20/cs	93
23189-110-401	Hook-up Kit (no pouch)	20/cs	93
23189-111-401	Holter hook-up Kit with case, 3-channel	20/cs	93
23189-116	LCF hook-up kit	ea	93
23189-116-401	Lifecard hook-upn Kit	ea	93
25-0020	10 Electrode patient cable, 4MM, AHA	ea	100
25-0028	10 Electrode patient cable, short, AHA	ea	100
25450-187	indv color lead wire, snap, white, 14 in	ea	94
25452-513	5-Lead wire set, snap, 46 cm/18 in, white/black/red/green/brown	ea	94
25452-519	7-Lead wire set, snap, 46 cm/18 in, white/black/red/green/orange/blue/brown	ea	94
25452-525	4-Lead wire set, snap, three 14 in leads and one 7 in white lead	ea	94
25452-529	4-Lead wire set, snap, 102 cm/40 in, Aria and evo belt clip config	ea	94
25452-535	4-Lead wire set, snap, three 14 in leads and one 14 in white lead	ea	94
26918-107	Memory card for 485	ea	98
26927-503	Patient cable, 30 in	ea	95
27249-7-401	Aria disposable pouch	50/bx	97
27565-5	ECG Electrode, foam	5/pk	95
27565-7	ECG Electrode, foam	7/pk	95
27565-9	ECG Electrode, clear tape	4/pk	95
27565-9CASE	ECG Electrode, clear tape	150/cs	95
27685-401	Disposable pouch, LC12	50/pk	97
2803FGA-24	Unifusor II 500cc w/aneroid gauge and thumbwheel valve, semi-disposable	24/bx	38
2803SGA-24	Unifusor II 500cc w/aneroid gauge and stopcock valve, semi-disposable	24/bx	38

Order Code	Description	UOM	Page(s)
2803X-24	Unifusor II 500cc, cuff only	24/bx	38
2804FGA-24	Unifusor II 1000cc w/aneroid gauge and thumbwheel valve, semi-disposable	24/bx	38
2804SGA-24	Unifusor II 1000cc w/aneroid gauge and stopcock valve, semi-disposable	24/bx	38
2804X-24	Unifusor II 1000cc, cuff only	24/bx	38
2805FGA-10	Unifusor II 3000cc w/aneroid gauge and thumbwheel valve, semi-disposable	10/bx	38
2805SGA-10	Unifusor II 3000cc w/aneroid gauge and stopcock valve, semi-disposable	10/bx	38
2805X-10	Unifusor II 3000cc, cuff only	10/bx	38
2903SGA-24	Unifusor II 500cc w/piston gauge and stopcock valve, semi-disposable	24/bx	38
2903TGA-24	Unifusor II 500cc w/piston gauge and thumbwheel valve, semi-disposable	24/bx	38
2904SGA-24	Unifusor II 1000cc w/piston gauge and stopcock valve, semi-disposable	24/bx	38
2904TGA-24	Unifusor II 1000cc w/piston gauge and thumbwheel valve, semi-disposable	24/bx	38
2905SGA-10	Unifusor II 3000cc w/piston gauge and stopcock valve, semi-disposable	10/bx	38
2905TGA-10	Unifusor II 3000cc w/piston gauge and thumbwheel valve, semi-disposable	10/bx	38
305240-001	Roll thermal paper, single-channel printing	10/bx	82
306655-003	Cardiac output cable, Baxter/Edwards/Abbott	ea	81
306655-004	Cardiac output cable, Becton Dickinson/Datex-Ohmeda	ea	81
307438-001	Roll thermal paper, single- or dual-channel printing	10/pk	82
307568-001	Z-fold thermal paper, 2- or 4- channel printing	10/bx	82
344-0056-00	Garment clip	ea	92
369014-003	Nafion tubing, single use	ea	68
369062-001	ECG electrode, adult, tape, single use	100/bx	50, 58
369063-001	ECG electrode, transparent tape, pediatric/neonatal	100/bx	51
376-0561-01	Pediatric Patient flow sensor	ea	102
378-0050-00	Water trap, single use	ea	66
392003-001	Skin prep pad	100/bx	99
41-0053US	12-lead hook-up set, AHA	ea	100
41-0080	CD12USB stress tets hook-up kit	ea	99
46-0063	Plastic mememory card cases	ea	98
46-0092	Fabric pouch	ea	97
46-0408	Fabric pouch set, 2 pouches with laces, cordlock & instructions	ea	97
46-0413	LCF opener tool	ea	98
46-0435	Patient diary, 7 day, German	ea	98
46-0459	Patient diary, 7 day, English	ea	98
46-0512	Disposable pouch, lifecard cf	100/pk	97
46-0573	Patient diary lcf, eng,ita, fra,spa,pol	ea	98
50-0042	Evo holder for waist belt	ea	97
50-0048	4 electrode hook-up kit	ea	93
50-0049	4 electrode hook-up kit with disposable chest pouch	ea	93
50-0050	Disposable chest pouch for evo recorder	50/pk	97
50-0190	4 electrode hook-up kit	20/cs	93
50-0191	4 electrode hook-up kit with chest pouch	20/cs	93
50-0220	LCF / evo pen	ea	98
685-0037-00	ECG pre-wired neonatal electrode, DIN unshielded, 3-lead, single use	100/bx	51
690-0003-01	Nellcor OxiMax DS-100A finger sensor, adult, reusable	ea	77, 88
690-0004-01	Nellcor OXI-A/N sensor, adult/neonatal, reusable	ea	77, 88
690-0018-00	Diaphragm dome, Datex-Ohmeda TA101OND	12/bx	80
690-0021-00	INVP interface cable for pressure transducers, Datex-Ohmeda/Becton Dickinson	ea	80
690-0028-01	Skin temperature probe, small, YSI729	ea	78
690-0039-01	Nellcor OXI-P/I sensor, pediatric/infant, reusable	ea	77, 88
690-0215-01	Nellcor DS-Y multi-site sensor, universal, reusable	ea	77, 88
690-0217-00	Masimo LNOP Adt adhesive sensor, adult, single use	20/bx	75, 87
690-0218-00	Masimo LNOP Adtx adhesive sensor, transparent tape, adult single use	20/bx	75, 87
690-0219-00	Masimo LNOP Pdt adhesive sensor, pediatric, single use	20/bx	75

Order Code	Description	UOM	Page(s)
690-0221-00	Masimo LNOP Neo-L adhesive sensor, neonatal, single use	20/bx	75
690-0226-00	Masimo LNOP Blue thumb/toe sensor, infant, single use	10/bx	75
690-0230-00	Masimo LNOP DC-I sensor, adult, reusable	ea	74, 87
690-0231-00	Masimo LNOP DC-IP sensor, pediatric, reusable	ea	74, 87
690-0232-00	Masimo LNOP YI multi-site sensor, reusable	ea	74, 87
690-0247-00	Masimo LNCS DC-I sensor, adult, reusable	ea	72, 86
690-0248-00	Masimo LNCS DC-IP sensor, pediatric, reusable	ea	72, 86
690-0249-00	Masimo LNCS TC-I tip clip sensor, reusable	ea	72, 86
690-0250-00	Masimo LNCS TF-I forehead sensor, adult, reusable	ea	72, 86
690-0251-00	Masimo LNCS Amtx adhesive sensor, adult, single use	20/bx	73, 86
690-0252-00	Masimo LNCS Pmtx adhesive sensor, pediatric, single use	20/bx	73
690-0258-00	Masimo LNCS Neo-L adhesive sensor, neonatal, single use	20/bx	73
690-0259-00	Masimo LNCS Inf-L adhesive sensor, infant, single use	20/bx	73
690-0262-00	ECG electrode, transparent tape, general/stress, adult	30/bx	50, 58
690-0272-00	Esophageal/rectal probe, reusable	ea	78
700-0006-00	ECG Lead Wire Set, TruLink, shielded, 3L pinch, AAMI, 61 cm/24 in	1 Set	48
700-0006-02	ECG Lead Wire Set, TruLink, shielded, 3L pinch, AAMI, 102 cm/40 in	1 Set	48
700-0006-08	ECG Lead Wire Set, TruLink, shielded, 5L, pinch, AAMI, 61 cm/24 in	1 Set	48
700-0006-10	ECG Lead Wire Set, TruLink, shielded, 5L, pinch, AAMI, 102 cm/40 in	1 Set	48
700-0006-16	ECG Lead Wire Set, TruLink, shielded, 10L, pinch, AAMI, 61 cm/24 in	1 Set	48
700-0006-18	ECG Lead Wire Set, TruLink, shielded, 10L, pinch, AAMI, 102 cm/40 in	1 Set	48
700-0006-30	ECG Combiner Disp Lead Wire Set, 3L, multi-pinch, AAMI, 127 cm/50 in	25/bx	41
700-0006-31	ECG Combiner Disp Lead Wire Set, 4V, multi-pinch, AAMI, 127 cm/50 in	25/bx	41
700-0006-32	ECG Combiner Disp Lead Wire Set, 5L, multi-pinch, AAMI, 127 cm/50 in	25/bx	41
700-0006-33	ECG Combiner Disp Lead Wire Set, 5V, multi-pinch, AAMI, 127 cm/50 in	25/bx	41
700-0006-34	ECG Combiner Disp Lead Wire Set, 6L, multi-pinch, AAMI, 127 cm/50 in	25/bx	41
700-0006-35	ECG Combiner Lead Wire Set, 3L, multi-pinch, AAMI, 74 cm/29 in	1 Set	41
700-0006-36	ECG Combiner Lead Wire Set, 4V, multi-pinch, AAMI, 74 cm/29 in	1 Set	41
700-0006-37	ECG Combiner Lead Wire Set, 5L, multi-pinch, AAMI, 74 cm/29 in	1 Set	41
700-0006-38	ECG Combiner Lead Wire Set, 5V-Leads, multi-pinch, AAMI, 74 cm/29 in	1 Set	41
700-0006-39	ECG Combiner Lead Wire Set, 6L, multi-pinch, AAMI, 74 cm/29 in	1 Set	41
700-0007-00	ECG Lead Wire Set, TruLink, shielded, 3L, snap, AAMI, 61 cm/24 in	1 Set	48
700-0007-02	ECG Lead Wire Set, TruLink, shielded, 3L, snap, AAMI, 102 cm/40 in	1 Set	48
700-0007-08	ECG Lead Wire Set, TruLink, shielded, 5L, snap, AAMI, 61 cm/24 in	1 Set	48
700-0007-10	ECG Lead Wire Set, TruLink, shielded, 5L, snap, AAMI, 102 cm/40 in	1 Set	48
700-0007-16	ECG Lead Wire Set, TruLink, shielded, 10L, snap, AAMI, 61 cm/24 in	1 Set	48
700-0007-18	ECG Lead Wire Set, TruLink, shielded, 10L, snap, AAMI, 102 cm/40 in	1 Set	48
700-0007-30	ECG Combiner Disp Lead Wire Set, 3L, multi-pinch, AAMI, 127 cm/50 in	25/bx	41
700-0007-31	ECG Combiner Disp Lead Wire Set, 4V, multi-pinch, AAMI, 127 cm/50 in	25/bx	41
700-0007-32	ECG Combiner Disp Lead Wire Set, 5L, multi-pinch, AAMI, 127 cm/50 in	25/bx	41
700-0007-33	ECG Combiner Disp Lead Wire Set, 5V, multi-pinch, AAMI, 127 cm/50 in	25/bx	41
700-0007-34	ECG Combiner Disp Lead Wire Set, 6L, multi-pinch, AAMI, 127 cm/50 in	25/bx	41
700-0007-35	ECG Combiner Lead Wire Set, 3L, multi-pinch, AAMI, 74 cm/29 in	1 Set	41
700-0007-36	ECG Combiner Lead Wire Set, 4V-Leads, multi-pinch, AAMI, 74 cm/29 in	1 Set	41
700-0007-37	ECG Combiner Lead Wire Set, 5L, multi-pinch, AAMI, 74 cm/29 in	1 Set	41
700-0007-38	ECG Combiner Lead Wire Set, 5V-Leads, multi-pinch, AAMI, 74 cm/29 in	1 Set	41
700-0007-39	ECG Combiner Lead Wire Set, 6L, multi-pinch, AAMI, 74 cm/29 in	1 Set	41
700-0008-00	ECG cable, shielded, 10-Lead, AAMI	ea	50
700-0008-06	ECG cable, shielded, 5-Lead, AAMI	ea	50
700-0008-08	ECG cable, shielded, 3-Lead, AAMI	ea	50
700-0008-12	ECG Neonatal cable, DIN unshielded, 3-Lead, AAMI	ea	51
700-0008-30	ECG Combiner Cable, 3-Lead, Datascope, AAMI	ea	44
700-0008-31	ECG Combiner Cable, 5-Lead, Datascope, AAMI	ea	44

Order Code	Description	UOM	Page(s)
700-0008-33	ECG Combiner Cable, 3-Lead, Datex-Ohmeda, AAMI	ea	44
700-0008-34	ECG Combiner Cable, 5-Lead, Datex-Ohmeda, AAMI	ea	44
700-0008-35	ECG Combiner Cable, 6-Lead, Datex-Ohmeda, AAMI	ea	44
700-0008-36	ECG Combiner Cable, 3-Lead, GE, AAMI	ea	44
700-0008-37	ECG Combiner Cable, 5-Lead, GE, AAMI	ea	44
700-0008-38	ECG Combiner Cable, 5+5-Lead, GE, AAMI	ea	44
700-0008-39	ECG Combiner Cable, 6-Lead, GE, AAMI	ea	44
700-0008-40	ECG Combiner Cable, 6+4-Lead, GE, AAMI	ea	44
700-0008-41	ECG Combiner Cable, 3-Lead, Medtronic LifePak 90 angle, AAMI	ea	45
700-0008-42	ECG Combiner Cable, 5-Lead, Medtronic LifePak 90 angle, AAMI	ea	45
700-0008-43	ECG Combiner Cable, 6-Lead, Medtronic LifePak 90 angle, AAMI	ea	45
700-0008-44	ECG Combiner Cable, 3-Lead, Nihon Kohden, AAMI	ea	45
700-0008-45	ECG Combiner Cable, 5-Lead, Nihon Kohden, AAMI	ea	45
700-0008-46	ECG Combiner Cable, 6-Lead, Nihon Kohden, AAMI	ea	45
700-0008-47	ECG Combiner Cable, 3-Lead, Philips, AAMI	ea	44
700-0008-48	ECG Combiner Cable, 5+5-Lead, Philips, AAMI	ea	45
700-0008-49	ECG Combiner Cable, 5-Lead, Philips, AAMI	ea	45
700-0008-50	ECG Combiner Cable, 6+4-Lead, Philips, AAMI	ea	45
700-0008-51	ECG Combiner Cable, 6-Lead, Philips, AAMI	ea	45
700-0008-52	ECG Combiner Cable, 3-Lead, Zoll 90 angle, AAMI	ea	45
700-0008-53	ECG Combiner Cable, 5-Lead, Zoll 90 angle, AAMI	ea	45
700-0008-54	ECG Combiner Cable, 6-Lead, Zoll 90 angle, AAMI	ea	45
700-0008-55	ECG Combiner Cable, 3-Lead Spacelabs monitors	ea	43
700-0008-56	ECG Combiner Cable, 5+5-Lead Spacelabs monitors	ea	43
700-0008-57	ECG Combiner Cable, 5-Lead Spacelabs monitors	ea	43
700-0008-58	ECG Combiner Cable, 6+4-Lead, Spacelabs monitors	ea	43
700-0008-59	ECG Combiner Cable, 6-Lead, Spacelabs monitors	ea	43
700-0008-60	ECG Combiner adapter, 5-Lead, Datascope Telemetry	ea	46
700-0008-61	ECG Combiner adapter, 3-Lead, Datascope	ea	46
700-0008-62	ECG Combiner adapter, 5-Lead, Datascope	ea	46
700-0008-63	ECG Combiner adapter, 3-Lead, Datex-Ohmeda	ea	46
700-0008-64	ECG Combiner adapter, 5-Lead, Datex-Ohmeda	ea	46
700-0008-65	ECG Combiner adapter, 5-Lead, GE Apex Pro FH Telemetry	ea	46
700-0008-66	ECG Combiner adapter, 6-Lead, GE Apex Pro FH Telemetry	ea	46
700-0008-67	ECG Combiner adapter, 3-Lead, GE	ea	47
700-0008-68	ECG Combiner adapter, 5-Lead, GE	ea	47
700-0008-69	ECG Combiner adapter, 6-Lead, GE	ea	47
700-0008-70	ECG Combiner adapter, 3-Lead, Nihon Kohden	ea	47
700-0008-71	ECG Combiner adapter, 6-Lead, Nihon Kohden	ea	47
700-0008-72	ECG Combiner adapter, 3/5-Lead, Philips Telemetry	ea	47
700-0008-73	ECG Combiner adapter, 3-Lead, Philips	ea	47
700-0008-74	ECG Combiner adapter, 5-Lead, Philips	ea	47
700-0008-75	ECG Combiner adapter, 3-Lead, Siemens Draeger	ea	47
700-0008-76	ECG Combiner adapter, 5-Lead, Siemens Draeger	ea	47
700-0008-78	ECG Combiner adapter, 3-Lead Spacelabs Monitor compatible	ea	43
700-0008-79	ECG Combiner adapter, 5-Lead Spacelabs Monitor compatible	ea	43
700-0008-80	ECG Combiner adapter, 5-Lead Spacelabs Telemetry compatible	ea	46,57
700-0008-88	ECG Combine adapter, Aria Telemetry, 96281, AAMI	ea	46,57
700-0010-00	ECG Combiner Lead Wire Set, 3L , multi-pinch, AAMI, 130 cm/51 in	1 Set	41
700-0010-01	ECG Combiner Lead Wire Set, 5L , multi-pinch, AAMI, 130 cm/51 in	1 Set	41
700-0010-02	ECG Combiner Lead Wire Set, 5L , multi-pinch, AAMI, 130 cm/51 in	1 Set	41
700-0010-03	ECG Combiner Radiolucent Lead Wire Set, multi-pinch, AAMI, 74 cm/29 in	1 Set	42
700-0010-04	ECG Combiner Radiolucent Lead Wire Set, 5L , multi-pinch, AAMI, 74 cm/29 in	1 Set	42

Order Code	Description	UOM	Page(s)
700-0010-05	ECG Combiner Radiolucent Lead Wire Set, 3L multi-pinch, AAMI, 122 cm/48 in	1 Set	42
700-0010-06	ECG Combiner Radiolucent Lead Wire Set, 5L multi-pinch, AAMI, 122 cm/48 in	1 Set	42
700-0010-07	ECG Combiner Radiolucent Lead Wire Set, 3L multi-pinch, AAMI, 183 cm/72 in	1 Set	42
700-0010-08	ECG Combiner Radiolucent Lead Wire Set, 5L multi-pinch, AAMI, 183 cm/72 in	1 Set	42
700-0010-09	ECG Lead Wires, DIN unshielded, snap, AAMI, 74 cm/29 in, black	5/pk	56
700-0010-10	ECG Lead Wires, DIN unshielded, snap, AAMI, 74 cm/29 in, brown	5/pk	56
700-0010-11	ECG Lead Wires, DIN unshielded, snap, AAMI, 74 cm/29 in, green	5/pk	56
700-0010-12	ECG Lead Wires, DIN unshielded, snap, AAMI, 74 cm/29 in, red	5/pk	56
700-0010-13	ECG Lead Wires, DIN unshielded, snap, AAMI, 74 cm/29 in, white	5/pk	56
700-0010-14	ECG Combiner Lead Wire Set, 5L, multi-pinch, AAMI,102 cm/ 40 in	1 Set	42
700-0010-15	ECG Combiner Lead Wire Set, 5L, snap, AAMI, 61 cm/24 in	1 Set	42
700-0010-16	ECG Combiner Lead Wire Set, 5L, multi-pinch, AAMI, 61 cm/24 in	1 Set	42
700-0014-00	Spacelabs adapter cable	ea	59
700-0015-00	TruLink ABP Interface Cable, 90343/91343/96217, 60 cm/24 in	ea	60
700-0027-00	CO2 sampling line, sidestream, 91517	ea	81
700-0028-00	INVBP interface adapter cable, dual pressure	ea	80
700-0029-00	Novamatrix sensor adapter cable	ea	71
700-0030-00	SPO2 cable	ea	69, 76, 85
700-0031-00	Dual temperature adapter cable, 30 cm/12 in	ea	78
700-0292-01	INVBP integrated adapter cable, Edwards Lifesciences, dual pressure	ea	80
700-0293-01	INVBP integrated adapter cable, Edwards Lifesciences, single pressure	ea	80
700-0294-01	INVBP integrated adapter cable, Abbott Labs/Hospira/ICU Medical, dual pressure	ea	80
700-0295-01	Single pressure integrated cable	ea	80
700-0316-00	Patient cable, 3-Channel, 6-lead, 22 in	ea	96
700-0317-00	Patient cable, 3-Channel, 3-lead, Belt Clip, 41 in	ea	96
700-0318-00	Patient cable, 3-channel	ea	96
700-0319-00	Patient cable, 3-Channel, 3-lead, Neck Pouch, 26 in	ea	96
700-0320-00	Patient cable, 2-Channel	ea	96
700-0321-00	Patient cable, 2-Channel, short	ea	96
700-0322-00	Patient cable, 3-Channel, 7-lead, 23 in	ea	96
700-0323-00	Patient cable, 2-Channel, 4-lead, Belt Clip, long	ea	96
700-0324-00	Patient cable, 2-Channel, 4-lead, Neck Pouch, short	ea	96
700-0789-01	Masimo LNOP sensor adapter cable	ea	74, 87
700-0791-00	INVBP interface adapter cable, single pressure	ea	80
700-0792-00	Nellcor OxiMax sensor adapter cable	ea	76, 88
700-0906-01	Masimo SET SPO2 sensor adapter cables	ea	72, 86
700-4000-00	Universal temperature probe extension cable, YSI 400 series compatible	ea	78
701-0122-00	Temperature probe overs	200/bx	88
701-4000-00	Skin temperture sensor	20/bx	79
701-4009-00	Esophageal/stethoscope temperature sensor, 9FR	20/bx	79
701-4012-00	Esophageal/stethoscope temperature sensor, 12FR	20/bx	79
701-4018-00	Esophageal/stethoscope temperature sensor, 18FR	20/bx	79
701-4024-00	Esophageal/stethoscope temperature sensor, 24 FR	20/bx	79
701-4109-00	Esophageal/rectal temperature probe	20/bx	79
701-4112-00	Esophageal/stethoscope temperature sensor, 12 FR	20/bx	79
701-4214-00	Foley catheter w/temperature probe, 14FR	20/bx	79
701-4216-00	Foley catheter w/temperature probe, 16FR	20/bx	79
701-4218-00	Foley catheter w/temperature probe, 18FR	20/bx	79
701-4300-00	Tympanic temperature sensor, adult	20/bx	79
704-0001-00	CO2 mainstream sensor, 91517	ea	65
704-0002-00	CO2 airway adapter, mainstream, 91517, reusable, adult	ea	65
704-0003-00	CO2 airway adapter, mainstream, 91517, reusable, neonate	ea	65
704-0004-00	CO2 airway adapter, mainstream, 91517, single use, neonate	24/bx	65

Order Code	Description	UOM	Page(s)
704-0005-00	CO2 airway adapter, mainstream, 91517, single use, adult	24/bx	65
704-0006-00	Cable holding clip	5/pk	65
704-0010-00	CO2 nasal cannula, 91517, adult	10/bx	63
704-0011-00	CO2 nasal cannula, 91517, infant	10/bx	63
704-0012-00	CO2 nasal cannula, 91517, pediatric	10/bx	63
704-0014-00	CO2 nasal cannula w/O2 delivery, 91517, adult	10/bx	63
704-0015-00	CO2 nasal cannula w/O2 delivery, 91517, infant/neonatal	10/bx	63
704-0016-00	CO2 nasal cannula w/O2 delivery, 91517, pediatric	10/bx	63
704-0017-00	CO2 oral/nasal cannula, 91517, adult	10/bx	63
704-0019-00	CO2 oral/nasal cannula, 91517, pediatric	10/bx	63
704-0020-00	CO2 oral/nasal cannula w/O2 delivery, 91517, adult	10/bx	63
704-0022-00	CO2 oral/nasal cannula w/O2 delivery, 91517, pediatric	10/bx	63
704-0023-00	Airway adapter + sampling line, sidestream, 91517, adult/pediatric	10/bx	64
704-0024-00	Airway adapter + sampling line, sidestream, 91517, pediatric/infant	10/bx	64
704-0025-00	Airway adapter + sampling line w/nafion, sidestream, 91517, adult/ped	10/bx	64
704-0026-00	Airway adapter + sampling line w/nafion, sidestream, 91517, ped/infant	10/bx	64
704-0027-00	CO2 sampling line, sidestream, 91517	10/bx	64
704-0028-00	CO2 sampling line w/nafion, sidestream, 91517	10/bx	64
704-0173-00	Airway adapter, mainstream, 92517, adult/pediatric	25/bx	62
704-0174-00	Airway adapter, mainstream, 92517, infant	10/bx	62
704-0175-00	Airway T-adapter, sidestream, adult/pediatric	25/bx	62
704-0176-00	CO2 nasal cannula, adult, 2m/6.5 ft	25/cs	61
704-0177-00	CO2 nasal cannula, pediatric, 2m/6.5 ft	25/cs	61
704-0178-00	CO2 nasal cannula, infant, 2m/6.5 ft	25/cs	61
704-0179-00	CO2 nasal cannula w/O2 delivery, infant, 2m/7 ft	25/cs	61
704-0180-00	CO2 nasal cannula w/O2 delivery, pediatric, 2m/7 ft	25/cs	61
704-0181-00	CO2 nasal cannula w/O2 delivery, adult, 2m/7 ft	25/cs	61
704-0182-00	CO2 oral/nasal cannula w/O2, delivery, pediatric, 2m/7 ft	25/cs	61
704-0183-00	CO2 oral/nasal cannula w/O2 delivery, adult, 2m/7 ft	25/cs	61
704-0184-00	CO2 nasal cannula, adult, 3m/10 ft	25/cs	61
704-0185-00	CO2 nasal cannula, pediatric, 3m/10 ft	25/cs	61
704-0186-00	CO2 nasal cannula w/O2 delivery, adult, 3m/10 ft	25/cs	61
704-0187-00	CO2 oral/nasal cannula w/O2 delivery, adult, 3m/10 ft	25/cs	61
704-0189-00	CO2 oral/nasal cannula w/O2 delivery, pediatric, 3m/10 ft	25/cs	61
704-0191-00	CO2 sampling line w/T-adapter, adult, 2.5 m/6 ft	25/bx	62
704-0192-00	CO2 sampling line w/T-adapter, adult, 3 m/10 ft	25/bx	62
704-0193-00	CO2 emergency nasal cannula, adult	25/cs	62
706-0182-00	Disposable cuff cover for reusable ABP cuffs, standard	50/bx	89
706-0183-00	Disposable cuff cover for reusable ABP cuffs, large	50/bx	89
706-0189-00	Disposable pouch for evo recorder	100/pk	97
706-0191-00	ECG Lead Wire Grouper, 5 Lead, 96281	5/bg	57
706-0200-00	Disposable Pouch, 90227	50/cs	91
712-1039-00	Quick Release Metal Fittings, female bayonet	3/bg	90
714-0018-00	BP hose, single tube, adult	ea	1, 2, 3, 4, 6, 7, 9, 84
714-0018-12	BP hose, dual tube to single tube, adult	ea	9
714-0019-01	BP hose, single tube, neonatal	ea	5, 8, 11
714-0019-12	BP hose, dual tube to single tube, neonatal	ea	11
714-0020-00	BP cuff, TruLink Nylon, 1T, child, 12-19 cm, MHP	ea	6
714-0021-00	BP cuff, Spacelabs Nylon, 1T, sm adult, 17-25 cm, MHP	ea	6
714-0023-00	BP cuff, Spacelabs Nylon, 1T, adult, 23-33 cm, MHP	ea	6
714-0025-00	BP cuff, TruLink Nylon, 1T, lg adult, 31-40 cm, MHP	ea	6
714-0027-00	BP cuff, TruLink Nylon, 1T, thigh, 38-50 cm, MHP	ea	6

Order Code	Description	UOM	Page(s)
714-0038-00	BP Slip Fit Connector, 3mm ID Adult NIBP Hose	ea	9
714-0043-00	BP TruLink dual hose conversion kit, adult hose assembly	ea	10
714-0045-00	BP Slip Fit Connector, 6mm ID Adult NIBP Hose	ea	9
714-0046-00	BP Slip Fit Connector, 5mm ID Adult NIBP Hose	ea	9
714-0047-00	BP Slip Fit Connector, 4mm ID Adult NIBP Hose	ea	9
714-0052-50	BP hose fitting, plastic, neonatal, 50/pk	50/pk	11
714-0053-10	BP hose fitting, neonatal	10/pk	11
714-0054-50	BP hose fitting, plastic, infant, 50/pk	50/pk	11
714-0056-00	BP Y adapter, adult, 3 mm/1/8 in, barb	5/pk	10
714-0057-00	BP Y adapter, dual to single hose, adult	5/pk	11
714-0057-10	BP hose hose connector, non-luer	10/pk	11
714-0060-00	BP hose, single tube, adult	ea	1, 2, 3, 4, 6, 7, 9, 84
714-0061-00	BP hose, single tube, adult	ea	1, 2, 3, 4, 6, 7, 9
714-0062-00	BP hose, single tube, neonatal	ea	5, 8, 11
714-0070-00	BP cuff, Spacelabs Nylon, 1T, child long, 12-19 cm, MHP	ea	6
714-0071-00	BP cuff, Spacelabs Nylon, 1T, sm adult long, 17-25 cm, MHP	ea	6
714-0073-00	BP cuff, Spacelabs Nylon, 1T, adult long, 23-33 cm, MHP	ea	6
714-0075-00	BP cuff, Spacelabs Nylon, 1T, lg adult long, 31-40 cm, MHP	ea	6
714-1020-00	BP cuff, Spacelabs Vinyl, 1T, child, 12-19 cm, HP	5/pk	7
714-1021-00	BP cuff, Spacelabs Vinyl, 1T, small adult, 17-25 cm, HP	5/pk	7
714-1023-00	BP cuff, Spacelabs Vinyl, 1T, adult, 23-33 cm, HP	5/pk	7
714-1025-00	BP cuff, Spacelabs Vinyl, 1T, large adult, 31-40 cm, HP	5/pk	7
714-1027-00	BP cuff, Spacelabs Vinyl, 1T, thigh, 38-50 cm, HP	5/pk	7
714-1028-01	BP cuff, Spacelabs Soft Neonatal, 1T, size 1, 3-6 cm, M316	5/pk	8
714-1029-01	BP cuff, Spacelabs Soft Neonatal, 1T, size 2, 4-8 cm, M316	5/pk	8
714-1030-01	BP cuff, Spacelabs Soft Neonatal, 1T, size 3, 6-11cm, M316	5/pk	8
714-1032-01	BP cuff, Spacelabs Soft Neonatal, 1T, size 4, 7-13 cm, M316	5/pk	8
714-1033-01	BP cuff, Spacelabs Soft Neonatal, 1T, size 5, 8-15 cm, M316	5/pk	8
714-1073-00	BP cuff, Spacelabs Vinyl, 1T, adult long, 23-33 cm, HP	5/pk	7
714-1075-00	BP cuff, Spacelabs Vinyl, 1T, large adult long, 31-40 cm, HP	5/pk	7
715-1121-10	BP cuff, TruLink Nylon, 1T, child, 13-20 cm, HP	ea	1, 84
715-1122-10	BP cuff, TruLink Vinyl, 1T, child, 13-20 cm, HP	5/bg	3, 84
715-1123-10	BP cuff, TruLink Soft, 1T, child, 13-20 cm, HP	5/bg	2, 84
715-1131-10	BP cuff, TruLink Nylon, 1T, sm adult, 18-26 cm, HP	ea	1, 84
715-1132-10	BP cuff, TruLink Vinyl, 1T, sm adult, 18-26 cm, HP	5/bg	3, 84
715-1133-10	BP cuff, TruLink Soft, 1T, sm adult, 18-26 cm, HP	5/bg	2, 84
715-1135-10	BP cuff, TruLink Yellow Vinyl, 1T, sm adult, 18-26 cm, HP	5/bg	4
715-1141-10	BP cuff, TruLink Nylon, 1T, adult, 26-36 cm, HP	ea	1, 84
715-1142-10	BP cuff, TruLink Vinyl, 1T, adult, 26-35 cm, HP	5/bg	3, 84
715-1143-10	BP cuff, TruLink Soft, 1T, adult, 26-35 cm, HP	5/bg	2, 84
715-1145-10	BP cuff, TruLink Yellow Vinyl, 1T, adult, 26-35 cm, HP	5/bg	4
715-1151-10	BP cuff, TruLink Nylon, 1T, adult long, 23-33 cm, HP	ea	1, 84
715-1152-10	BP cuff, TruLink Vinyl, 1T, adult long, 23-33 cm, HP	5/bg	3, 84
715-1153-10	BP cuff, TruLink Soft, 1T, adult long, 23 - 33 cm, HP	5/bg	2, 84
715-1161-10	BP cuff, TruLink Nylon, 1T, lg adult, 32-42 cm, HP	ea	1, 84
715-1162-10	BP cuff, TruLink Vinyl, 1T, lg adult, 32-42 cm, HP	5/bg	3, 84
715-1163-10	BP cuff, TruLink Soft, 1T, lg adult, 32-42 cm, HP	5/bg	2, 84
715-1165-10	BP cuff, TruLink Yellow Vinyl, 1T, lg adult, 32-42 cm, HP	5/bg	4
715-1171-10	BP cuff, TruLink Nylon, 1T, lg adult long, 35-44 cm, HP	ea	1, 84
715-1172-10	BP cuff, TruLink Vinyl, 1T, lg adult long, 35-44 cm, HP	5/bg	3, 84
715-1173-10	BP cuff, TruLink Soft, 1T, lg adult long, 35-44 cm, HP	5/bg	2, 84

Order Code	Description	UOM	Page(s)
715-1181-10	BP cuff, TruLink Nylon, 1T, thigh, 42-50 cm, HP	ea	1, 84
715-1182-10	BP cuff, TruLink Vinyl, 1T, thigh, 42-50 cm, HP	5/bg	3, 84
715-1183-10	BP cuff, TruLink Soft, 1T, thigh, 42-50 cm, HP	5/bg	2, 84
715-2114-10	BP cuff, TruLink Non-PVC Neonatal, 1T, size 1, 3-6 cm, M316	10/bx	5
715-2124-10	BP cuff, TruLink Non-PVC Neonatal, 1T, size 2, 4-8 cm, M316	10/bx	5
715-2134-10	BP cuff, TruLink Non-PVC Neonatal, 1T, size 3, 6-11 cm, M316	10/bx	5
715-2144-10	BP cuff, TruLink Non-PVC Neonatal, 1T, size 4, 7-13 cm, M316	10/bx	5
715-2154-10	BP cuff, TruLink Non-PVC Neonatal, 1T, size 5, 8-15 cm, M316	10/bx	5
717-0071-00	ECG Electrode, AMBU	5/pk	95
717-0072-00	electrodes, adult, foam, w/snap	25/pk	100
717-0408-00	Exhaust Line	5/pk	66
717-0446-00	D-FEND Water Trap Container	5/pk	66
717-0447-00	D-FEND Water Trap Cartridge	10/pk	67
717-0448-00	D-FEND+ Water Trap Cartridge	10/pk	67
719-0001-10	BISX Pediatric Sensor	25/bx	83
719-0002-10	BISX Quarto Sensor	25/bx	83
719-0003-10	BISX Extend Sensor	25/bx	83
719-0004-01	BISX Patient Interface cable	ea	83
719-0006-20	BISx Pod	ea	83
719-0007-01	BISX cable	ea	83
803FGA-24	Unifusor 500cc w/aneroid gauge and thumbwheel valve	24/bx	37
803SGA-24	Unifusor 500cc w/aneroid gauge and stopcock valve	24/bx	37
804FGA-24	Unifusor 1000cc w/aneroid gauge and thumbwheel valve	24/bx	37
804SGA-24	Unifusor 1000cc w/aneroid gauge and stopcock valve	24/bx	37
805FGA-10	Unifusor 3000cc w/aneroid gauge and thumbwheel valve	10/bx	37
805SGA-10	Unifusor 3000cc w/aneroid gauge and stopcock valve	10/bx	37
899F-04	Unifusor II aneroid gauge and thumbwheel valve assy; semi-disposable	4/bx	38,39
899S-04	Unifusor II aneroid gauge and stopcock valve assy; semi-disposable	4/bx	38,39
903SGA-24	Unifusor Classic 500cc w/piston gauge and stopcock valve	24/bx	37
903TGA-24	Unifusor Classic 500cc w/piston gauge and thumbwheel valve	24/bx	37
904SGA-24	Unifusor Classic 1000cc w/piston gauge and stopcock valve	24/bx	37
904TGA-24	Unifusor Classic 1000cc w/piston gauge and thumbwheel valve	24/bx	37
905SGA-10	Unifusor Classic 3000cc w/piston gauge and stopcock valve	10/bx	37
905TGA-10	Unifusor Classic 3000cc w/piston gauge and thumbwheel valve	10/bx	37
91-0013-00	Cable, Tango, RS-232	ea	99
98-0002-50	Microphone pads, 5 PK	5/pk	99
98-0006-00	K-Sound microphone	ea	99
98-0061-01	Cuff, Orbit-K, small adult	ea	99
999S-04	Unifusor II piston gauge and stopcock valve assembly; semi-disposable	4/bx	38,39
999T-04	Unifusor II piston gauge and thumbwheel valve assembly; semi-disposable	4/bx	38,39
ADP-182YBPF	BP Y adapter, 182YBPF	ea	35
ADP-182YRX4	BP Y adapter, 182YRX4	ea	35
ADP-5322YBPF	BP Y adapter, 5322YBPF	ea	35
ADP-5322YRX4	BP Y adapter, 5322YRX4	ea	35
BRDT3854-05	BP cuff, SoftCheck Curve, 2T, XL-Adult, 38-54 cm, no fitting	5/bx	25
BRDT3854DM-05	BP cuff, SoftCheck Curve, 2T, XL-adult, 38-54 cm, DM	5/bx	25
BRDT3854FL-05	BP cuff, SoftCheck Curve, 2T, XL-adult, 38-54 cm, FL	5/pk	25
BRDT3854HP-05	BP cuff, SoftCheck Curve, 2T, XL-adult, 38-54 cm, HP	5/bx	25
BRDT3854MB-05	BP cuff, SoftCheck Curve, 2T, XL-adult, 38-54 cm, manual	5/bx	25
BRDT3854ML-05	BP cuff, SoftCheck Curve, 2T, XL-adult, 38-54 cm, ML	5/pk	25
BRDT3854MQ-05	BP cuff, SoftCheck Curve, 2T, XL-adult, 38-54 cm, MQ	5/bx	25
BRDT3854MQMF-05	BP cuff, SoftCheck Curve, 2T, XL-adult, 38-54 cm, MQMF	5/bx	25
BRDT3854MQML-05	BP cuff, SoftCheck Curve, 2T, XL-adult, 38-54 cm, MQML	5/pk	25

Order Code	Description	UOM	Page(s)
BRST3854-05	BP cuff, SoftCheck Curve, 1T, XL-adult, 38-54 cm, DM	5/bx	25
BRST3854DM-05	BP cuff, SoftCheck Curve, 1T, XL-adult, 38-54 cm, FL	5/bx	25
BRST3854FL-05	BP cuff, SoftCheck Curve, 1T, XL-adult, 38-54 cm, HP	5/bx	25
BRST3854HP-05	BP cuff, SoftCheck Curve, 1T, XL-adult, 38-54 cm, ML	5/bx	25
BRST3854ML-05	BP cuff, SoftCheck Curve, 1T, XL-adult, 38-54 cm, MQ	5/bx	25
BRST3854MQ-05	BP cuff, SoftCheck Curve, 1T, XL-adult, 38-54 cm, MQML	5/bx	25
BRST3854MQML-05	BP cuff, SoftCheck Curve, 1T, XL-adult, 38-54 cm, MQML	5/bx	25
BRUD3854	BP cuff, Ultracheck Curve, 2T, XL adult, 38-54 cm, no fitting	ea	18
BRUD3854DM	BP cuff, Ultracheck Curve, 2T, XL adult, 38-54 cm, DM	ea	18
BRUD3854HP	BP cuff, Ultracheck Curve, 2T, XL adult, 38-54 cm, HP	ea	18
BRUD3854MB	BP cuff, Ultracheck Curve, 2T, XL adult, 38-54 cm, manual	ea	18
BRUD3854ML	BP cuff, Ultracheck Curve, 2T, XL adult, 38-54 cm, ML	ea	18
BRUD3854MQ	BP cuff, Ultracheck Curve, 2T, XL adult, 38-54 cm, MQ	ea	18
BRUD3854MQMF	BP cuff, Ultracheck Curve, 2T, XL adult, 38-54 cm, MQMF	ea	18
BRUD3854MQML	BP cuff, Ultracheck Curve, 2T, XL adult, 38-54 cm, MQML	ea	18
BRUS3854	BP cuff, Ultracheck Curve, 1T, XL adult, 38-54 cm, no fitting	ea	18
BRUS3854DM	BP cuff, Ultracheck Curve, 1T, XL adult, 38-54 cm, DM	ea	18
BRUS3854FL	BP cuff, Ultracheck Curve, 1T, XL adult, 38-54 cm, FL	ea	18
BRUS3854HP	BP cuff, Ultracheck Curve, 1T, XL adult, 38-54 cm, HP	ea	18
BRUS3854MB	BP cuff, Ultracheck Curve, 1T, XL adult, 38-54 cm, manual	ea	18
BRUS3854MHP	BP cuff, Ultracheck Curve, 1T, XL adult, 38-54 cm, MHP	ea	18
BRUS3854ML	BP cuff, Ultracheck Curve, 1T, XL adult, 38-54 cm, ML	ea	18
BRUS3854MQ	BP cuff, Ultracheck Curve, 1T, XL adult, 38-54 cm, MQ	ea	18
CFC256MB	Flashcard LifeCard CF, 256MB	ea	98
CFC90MB	Flashcard Memory, 90MByte	ea	98
DT0814-05	BP cuff, SoftCheck Fabric, 2T infant, 8-14 cm, no fitting	5/bg	24
DT0814DM-05	BP cuff, SoftCheck Fabric, 2T, infant, 8-14 cm, DM	5/bg	22
DT0814DM-20	BP cuff, SoftCheck Fabric, 2T, infant, 8-14 cm, DM	20/bx	22
DT0814MB-05	BP cuff, SoftCheck Fabric, 2T, infant, 8-14 cm, manual	5/bg	24
DT0814MQ-05	BP cuff, SoftCheck Fabric, 2T, infant, 8-14 cm, MQ	5/bg	23
DT0814MQ-20	BP cuff, SoftCheck Fabric, 2T, infant, 8-14 cm, MQ	20/bx	23
DT0814MQMF-05	BP cuff, SoftCheck Fabric, 2T, infant, 8-14 cm, MQMF	5/bg	24
DT1320-05	BP cuff, SoftCheck Fabric, 2T, child, 13-20 cm, no fitting	5/bg	24
DT1320DM-05	BP cuff, SoftCheck Fabric, 2T, child, 13-20 cm, DM	5/bg	22
DT1320DM-20	BP cuff, SoftCheck Fabric, 2T, child, 13-20 cm, DM	20/bx	22
DT1320MB-05	BP cuff, SoftCheck Fabric, 2T, child, 13-20 cm, manual	5/bg	24
DT1320MQ-05	BP cuff, SoftCheck Fabric, 2T, child, 13-20 cm, MQ	5/bg	23
DT1320MQ-20	BP cuff, SoftCheck Fabric, 2T, child, 13-20 cm, MQ	20/bx	23
DT1320MQMF-05	BP cuff, SoftCheck Fabric, 2T, child, 13-20 cm, MQMF	5/bg	24
DT1826-05	BP cuff, SoftCheck Fabric, 2T, sm adult, 18-26 cm, no fitting	5/bg	24
DT1826DM-05	BP cuff, SoftCheck Fabric, 2T, sm adult, 18-26 cm, DM	5/bg	22
DT1826DM-20	BP cuff, SoftCheck Fabric, 2T, sm adult, 18-26 cm, DM	20/bx	22
DT1826MB-05	BP cuff, SoftCheck Fabric, 2T, sm adult, 18-26 cm, manual	5/bg	24
DT1826MQ-05	BP cuff, SoftCheck Fabric, 2T, sm adult, 18-26 cm, MQ	5/bg	23
DT1826MQ-20	BP cuff, SoftCheck Fabric, 2T, sm adult, 18-26 cm, MQ	20/bx	23
DT1826MQMF-05	BP cuff, SoftCheck Fabric, 2T, sm adult, 18-26 cm, MQMF	5/bg	24
DT2635-05	BP cuff, SoftCheck Fabric, 2T, adult, 26-35 cm, no fitting	5/bg	24
DT2635DM-05	BP cuff, SoftCheck Fabric, 2T, adult, 26-35 cm, DM	5/bg	22
DT2635DM-20	BP cuff, SoftCheck Fabric, 2T, adult, 26-35 cm, DM	20/bx	22
DT2635MB-05	BP cuff, SoftCheck Fabric, 2T, adult, 26-35 cm, manual	5/bg	24
DT2635MQ-05	BP cuff, SoftCheck Fabric, 2T, adult, 26-35 cm, MQ	5/bg	23
DT2635MQ-20	BP cuff, SoftCheck Fabric, 2T, adult, 26-35 cm, MQ	20/bx	23
DT2635MQMF-05	BP cuff, SoftCheck Fabric, 2T, adult, 26-35 cm, MQMF	5/bg	24

Order Code	Description	UOM	Page(s)
DT2938-05	BP cuff, SoftCheck Fabric, 2T, adult long, 29-38 cm, no fitting	5/bg	24
DT2938DM-05	BP cuff, SoftCheck Fabric, 2T, adult long, 29-38 cm, DM	5/bg	22
DT2938DM-20	BP cuff, SoftCheck Fabric, 2T, adult long, 29-38 cm, DM	20/bx	22
DT2938MB-05	BP cuff, SoftCheck Fabric, 2T, adult long, 29-38 cm, manual	5/bg	24
DT2938MQ-05	BP cuff, SoftCheck Fabric, 2T, adult long, 29-38 cm, MQ	5/bg	23
DT2938MQ-20	BP cuff, SoftCheck Fabric, 2T, adult long, 29-38 cm, MQ	20/bx	23
DT2938MQMF-05	BP cuff, SoftCheck Fabric, 2T, adult long, 29-38 cm, MQMF	5/bg	24
DT3242-05	BP cuff, SoftCheck Fabric, 2T, lg adult, 32-42 cm, no fitting	5/bg	24
DT3242DM-05	BP cuff, SoftCheck Fabric, 2T, lg adult, 32-42 cm, DM	5/bg	22
DT3242DM-20	BP cuff, SoftCheck Fabric, 2T, lg adult, 32-42 cm, DM	20/bx	22
DT3242MB-05	BP cuff, SoftCheck Fabric, 2T, lg adult, 32-42 cm, manual	5/bg	24
DT3242MQ-05	BP cuff, SoftCheck Fabric, 2T, lg adult, 32-42 cm, MQ	5/bg	23
DT3242MQ-20	BP cuff, SoftCheck Fabric, 2T, lg adult, 32-42 cm, MQ	20/bx	23
DT3242MQMF-05	BP cuff, SoftCheck Fabric, 2T, lg adult, 32-42 cm, MQMF	5/bg	24
DT3544-05	BP cuff, SoftCheck Fabric, 2T, lg adult long, 35-44 cm, no fitting	5/bg	24
DT3544DM-05	BP cuff, SoftCheck Fabric, 2T, lg adult long, 35-44 cm, DM	5/bg	22
DT3544DM-20	BP cuff, SoftCheck Fabric, 2T, lg adult long, 35-44 cm, DM	20/bx	22
DT3544MB-05	BP cuff, SoftCheck Fabric, 2T, lg adult long, 35-44 cm, manual	5/bg	24
DT3544MQ-05	BP cuff, SoftCheck Fabric, 2T, lg adult long, 35-44 cm, MQ	5/bg	23
DT3544MQ-20	BP cuff, SoftCheck Fabric, 2T, lg adult long, 35-44 cm, MQ	20/bx	23
DT3544MQMF-05	BP cuff, SoftCheck Fabric, 2T, lg adult long, 35-44 cm, MQMF	5/bg	24
DT4250-05	BP cuff, SoftCheck Fabric, 2T, thigh, 42-50 cm, no fitting	5/bg	24
DT4250DM-05	BP cuff, SoftCheck Fabric, 2T, thigh, 42-50 cm, DM	5/bg	22
DT4250DM-20	BP cuff, SoftCheck Fabric, 2T, thigh, 42-50 cm, DM	20/bx	22
DT4250MB-05	BP cuff, SoftCheck Fabric, 2T, thigh, 42-50 cm, manual	5/bg	24
DT4250MQ-05	BP cuff, SoftCheck Fabric, 2T, thigh, 42-50 cm, MQ	5/bg	23
DT4250MQ-20	BP cuff, SoftCheck Fabric, 2T, thigh, 42-50 cm, MQ	20/bx	23
DT4250MQMF-05	BP cuff, SoftCheck Fabric, 2T, thigh, 42-50 cm, MQMF	5/bg	24
M1050-4	Reusable Y patient circuit, adult	ea	103
MS1R-107048-A0	ECG cable, 4MM, banana connector, AHA	ea	100
MS1R-107371-A0	Precordial suction electrodes, adult, 6/set	ea	100
MS1R-107584-A0	ECG cable, 4MM, grabber connector, AHA	ea	100
MS1R-32461	Precordial suction electrodes, child, 6/set	ea	100
MS1R-32462	Limb clamp electrodes, child, 4/set	ea	100
MS1R-32463	Banana socet adapters, clip/snap, 10/ST	ea	100
NN1DT-10	BP cuff, SoftCheck Fabric Neonatal, 2T, size 1, 3-6 cm, ML	10/bx	31
NN1ST-10	BP cuff, SoftCheck Fabric Neonatal, 1T, size 1, 3-6 cm, ML	10/bx	31
NN1ST-HP-10	BP cuff, SoftCheck Fabric Neonatal, 1T, size 1, 3-6 cm, HP	10/bx	31
NN1ST-M316-10	BP cuff, SoftCheck Fabric Neonatal, 1T, size 1, 3-6 cm, M316	10/bx	31
NN2DT-10	BP cuff, SoftCheck Fabric Neonatal, 2T, size 2, 4-8 cm, ML	10/bx	31
NN2ST-10	BP cuff, SoftCheck Fabric Neonatal, 1T, size 2, 4-8 cm, ML	10/bx	31
NN2ST-HP-10	BP cuff, SoftCheck Fabric Neonatal, 1T, size 2, 4-8 cm, HP	10/bx	31
NN2ST-M316-10	BP cuff, SoftCheck Fabric Neonatal, 1T, size 2, 4-8 cm, M316	10/bx	31
NN3DT-10	BP cuff, SoftCheck Fabric Neonatal, 2T, size 3, 6-11 cm, ML	10/bx	31
NN3ST-10	BP cuff, SoftCheck Fabric Neonatal, 1T, size 3, 6-11 cm, ML	10/bx	31
NN3ST-HP-10	BP cuff, SoftCheck Fabric Neonatal, 1T, size 3, 6-11 cm, HP	10/bx	31
NN3ST-M316-10	BP cuff, SoftCheck Fabric Neonatal, 1T, size 3, 6-11 cm, M316	10/bx	31
NN4DT-10	BP cuff, SoftCheck Fabric Neonatal, 2T, size 4, 7-13 cm, ML	10/bx	31
NN4ST-10	BP cuff, SoftCheck Fabric Neonatal, 1T, size 4, 7-13 cm, ML	10/bx	31
NN4ST-HP-10	BP cuff, SoftCheck Fabric Neonatal, 1T, size 4, 7-13 cm, HP	10/bx	31
NN4ST-M316-10	BP cuff, SoftCheck Fabric Neonatal, 1T, size 4, 7-13 cm, M316	10/bx	31
NN5DT-10	BP cuff, SoftCheck Fabric Neonatal, 2T, size 5, 8-15 cm, ML	10/bx	31
NN5ST-10	BP cuff, SoftCheck Fabric Neonatal, 1T, size 5, 8-15 cm, ML	10/bx	31

Order Code	Description	UOM	Page(s)
NN5ST-HP-10	BP cuff, SoftCheck Fabric Neonatal, 1T, size 5, 8-15 cm, HP	10/bx	31
NN5ST-M316-10	BP cuff, SoftCheck Fabric Neonatal, 1T, size 5, 8-15 cm, M316	10/bx	31
PK0001	Fabric waist belt	ea	97
PNN1DT-10	BP cuff, SoftCheck Free Neonatal, 2T, size 1, 3-6 cm, ML	10/bx	34
PNN1ST-10	BP cuff, SoftCheck Free Neonatal, 1T, size 1, 3-6 cm, ML	10/bx	34
PNN1ST-HP-10	BP cuff, SoftCheck Free Neonatal, 1T, size 1, 3-6 cm, HP	10/bx	34
PNN1ST-M316-10	BP cuff, SoftCheck Free Neonatal, 1T, size 1, 3-6 cm, M316	10/bx	34
PNN2DT-10	BP cuff, SoftCheck Free Neonatal, 2T, size 2, 4-8 cm, ML	10/bx	34
PNN2ST-10	BP cuff, SoftCheck Free Neonatal, 1T, size 2, 4-8 cm, ML	10/bx	34
PNN2ST-HP-10	BP cuff, SoftCheck Free Neonatal, 1T, size 2, 4-8 cm, HP	10/bx	34
PNN2ST-M316-10	BP cuff, SoftCheck Free Neonatal, 1T, size 2, 4-8 cm, M316	10/bx	34
PNN3DT-10	BP cuff, SoftCheck Free Neonatal, 2T, size 3, 6-11 cm, ML	10/bx	34
PNN3ST-10	BP cuff, SoftCheck Free Neonatal, 1T, size 3, 6-11 cm, ML	10/bx	34
PNN3ST-HP-10	BP cuff, SoftCheck Free Neonatal, 1T, size 3, 6-11 cm, HP	10/bx	34
PNN3ST-M316-10	BP cuff, SoftCheck Free Neonatal, 1T, size 3, 6-11 cm, M316	10/bx	34
PNN4DT-10	BP cuff, SoftCheck Free Neonatal, 2T, size 4, 7-13 cm, ML	10/bx	34
PNN4ST-10	BP cuff, SoftCheck Free Neonatal, 1T, size 4, 7-13 cm, ML	10/bx	34
PNN4ST-HP-10	BP cuff, SoftCheck Free Neonatal, 1T, size 4, 7-13 cm, HP	10/bx	34
PNN4ST-M316-10	BP cuff, SoftCheck Free Neonatal, 1T, size 4, 7-13 cm, M316	10/bx	34
PNN5DT-10	BP cuff, SoftCheck Free Neonatal, 2T, size 5, 8-15 cm, ML	10/bx	34
PNN5ST-10	BP cuff, SoftCheck Free Neonatal, 1T, size 5, 8-15 cm, ML	10/bx	34
PNN5ST-HP-10	BP cuff, SoftCheck Free Neonatal, 1T, size 5, 8-15 cm, HP	10/bx	34
PNN5ST-M316-10	BP cuff, SoftCheck Free Neonatal, 1T, size 5, 8-15 cm, M316	10/bx	34
RT-41	Skintact electrode tab, (100/PK)	100/pk	100
ST0814-05	BP cuff, SoftCheck Fabric, 1T, infant, 8-14 cm, no fitting	5/bg	21
ST0814DM-05	BP cuff, SoftCheck Fabric, 1T, infant, 8-14 cm, DM	5/bg	20
ST0814FL-05	BP cuff, SoftCheck Fabric, 1T, infant, 8-14 cm, FL	5/bg	21
ST0814HP-05	BP cuff, SoftCheck Fabric, 1T, infant, 8-14 cm, HP	5/bg	19
ST0814HP-20	BP cuff, SoftCheck Fabric, 1T, infant, 8-14 cm, HP	20/bx	19
ST0814ML-05	BP cuff, SoftCheck Fabric, 1T, infant, 8-14 cm, ML	5/bg	19
ST0814ML-20	BP cuff, SoftCheck Fabric, 1T, infant, 8-14 cm, ML	20/bx	19
ST0814MQ-05	BP cuff, SoftCheck Fabric, 1T, infant, 8-14 cm, MQ	5/bg	20
ST0814MQ-20	BP cuff, SoftCheck Fabric, 1T, infant, 8-14 cm, MQ	20/bx	20
ST0814MQML-05	BP cuff, SoftCheck Fabric, 1T, infant, 8-14 cm, MQML	5/bg	21
ST1320-05	BP cuff, SoftCheck Fabric, 1T, child, 13-20 cm, no fitting	5/bg	21
ST1320DM-05	BP cuff, SoftCheck Fabric, 1T, child, 13-20 cm, DM	5/bg	20
ST1320FL-05	BP cuff, SoftCheck Fabric, 1T, child, 13-20 cm, FL	5/bg	21
ST1320HP-05	BP cuff, SoftCheck Fabric, 1T, child, 13-20 cm, HP	5/bg	19
ST1320HP-20	BP cuff, SoftCheck Fabric, 1T, child, 13-20 cm, HP	20/bx	19
ST1320ML-05	BP cuff, SoftCheck Fabric, 1T, child, 13-20 cm, ML	5/bg	19
ST1320ML-20	BP cuff, SoftCheck Fabric, 1T, child, 13-20 cm, ML	20/bx	19
ST1320MQ-05	BP cuff, SoftCheck Fabric, 1T, child, 13-20 cm, MQ	5/bg	20
ST1320MQ-20	BP cuff, SoftCheck Fabric, 1T, child, 13-20 cm, MQ	20/bx	20
ST1320MQML-05	BP cuff, SoftCheck Fabric, 1T, child, 13-20 cm, MQML	5/bg	21
ST1826-05	BP cuff, SoftCheck Fabric, 1T, sm adult, 18-26 cm, no fitting	5/bg	21
ST1826DM-05	BP cuff, SoftCheck Fabric, 1T, sm adult, 18-26 cm, DM	5/bg	20
ST1826FL-05	BP cuff, SoftCheck Fabric, 1T, sm adult, 18-26 cm, FL	5/bg	21
ST1826HP-05	BP cuff, SoftCheck Fabric, 1T, sm adult, 18-26 cm, HP	5/bg	19
ST1826HP-20	BP cuff, SoftCheck Fabric, 1T, sm adult, 18-26 cm, HP	20/bx	19
ST1826ML-05	BP cuff, SoftCheck Fabric, 1T, sm adult, 18-26 cm, ML	5/bg	19
ST1826ML-20	BP cuff, SoftCheck Fabric, 1T, sm adult, 18-26 cm, ML	20/bx	19
ST1826MQ-05	BP cuff, SoftCheck Fabric, 1T, sm adult, 18-26 cm, MQ	5/bg	20
ST1826MQ-20	BP cuff, SoftCheck Fabric, 1T, sm adult, 18-26 cm, MQ	20/bx	20

Order Code	Description	UOM	Page(s)
ST1826MQML-05	BP cuff, SoftCheck Fabric, 1T, sm adult, 18-26 cm, MQML	5/bg	21
ST2111	2.0 L Rebreathing bag, 22mm neck	ea	103
ST2635-05	BP cuff, SoftCheck Fabric, 1T, adult, 25-35 cm, no fitting	5/bg	21
ST2635DM-05	BP cuff, SoftCheck Fabric, 1T, adult, 25-35 cm, DM	5/bg	20
ST2635FL-05	BP cuff, SoftCheck Fabric, 1T, adult, 25-35 cm, FL	5/bg	21
ST2635HP-05	BP cuff, SoftCheck Fabric, 1T, adult, 25-35 cm, HP	5/bg	19
ST2635HP-20	BP cuff, SoftCheck Fabric, 1T, adult, 25-35 cm, HP	20/bx	19
ST2635ML-05	BP cuff, SoftCheck Fabric, 1T, adult, 25-35 cm, ML	5/bg	19
ST2635ML-20	BP cuff, SoftCheck Fabric, 1T, adult, 25-35 cm, ML	20/bx	19
ST2635MQ-05	BP cuff, SoftCheck Fabric, 1T, adult, 25-35 cm, MQ	5/bg	20
ST2635MQ-20	BP cuff, SoftCheck Fabric, 1T, adult, 25-35 cm, MQ	20/bx	20
ST2635MQML-05	BP cuff, SoftCheck Fabric, 1T, adult, 25-35 cm, MQML	5/bg	21
ST2670	Malleable face mask, size 4	ea	103
ST2938-05	BP cuff, SoftCheck Fabric, 1T, adult long, 29-38 cm, no fitting	5/bg	21
ST2938DM-05	BP cuff, SoftCheck Fabric, 1T, adult long, 29-38 cm, DM	5/bg	21
ST2938FL-05	BP cuff, SoftCheck Fabric, 1T, adult long, 29-38 cm, FL	5/bg	21
ST2938HP-05	BP cuff, SoftCheck Fabric, 1T, adult long, 29-38 cm, HP	5/bg	19
ST2938HP-20	BP cuff, SoftCheck Fabric, 1T, adult long, 29-38 cm, HP	20/bx	19
ST2938ML-05	BP cuff, SoftCheck Fabric, 1T, adult long, 29-38 cm, ML	5/bg	19
ST2938ML-20	BP cuff, SoftCheck Fabric, 1T, adult long, 29-38 cm, ML	20/bx	19
ST2938MQ-05	BP cuff, SoftCheck Fabric, 1T, adult long, 29-38 cm, MQ	5/bg	20
ST2938MQ-20	BP cuff, SoftCheck Fabric, 1T, adult long, 29-38 cm, MQ	20/bx	20
ST2938MQML-05	BP cuff, SoftCheck Fabric, 1T, adult long, 29-38 cm, MQML	5/bg	21
ST3242-05	BP cuff, SoftCheck Fabric, 1T, lg adult, 32-42 cm, no fitting	5/bg	21
ST3242DM-05	BP cuff, SoftCheck Fabric, 1T, lg adult, 32-42 cm, DM	5/bg	20
ST3242FL-05	BP cuff, SoftCheck Fabric, 1T, lg adult, 32-42 cm, FL	5/bg	21
ST3242HP-05	BP cuff, SoftCheck Fabric, 1T, lg adult, 32-42 cm, HP	5/bg	19
ST3242HP-20	BP cuff, SoftCheck Fabric, 1T, lg adult, 32-42 cm, HP	20/bx	19
ST3242ML-05	BP cuff, SoftCheck Fabric, 1T, lg adult, 32-42 cm, ML	5/bg	19
ST3242ML-20	BP cuff, SoftCheck Fabric, 1T, lg adult, 32-42 cm, ML	20/bx	19
ST3242MQ-05	BP cuff, SoftCheck Fabric, 1T, lg adult, 32-42 cm, MQ	5/bg	20
ST3242MQ-20	BP cuff, SoftCheck Fabric, 1T, lg adult, 32-42 cm, MQ	20/bx	20
ST3242MQML-05	BP cuff, SoftCheck Fabric, 1T, lg adult, 32-42 cm, MQML	5/bg	21
ST3544-05	BP cuff, SoftCheck Fabric, 1T, lg adult long, 35-44 cm, no fitting	5/bg	21
ST3544DM-05	BP cuff, SoftCheck Fabric, 1T, lg adult long, 35-44 cm, DM	5/bg	20
ST3544FL-05	BP cuff, SoftCheck Fabric, 1T, lg adult long, 35-44 cm, FL	5/bg	21
ST3544HP-05	BP cuff, SoftCheck Fabric, 1T, lg adult long, 35-44 cm, HP	5/bg	19
ST3544HP-20	BP cuff, SoftCheck Fabric, 1T, lg adult long, 35-44 cm, HP	20/bx	19
ST3544ML-05	BP cuff, SoftCheck Fabric, 1T, lg adult long, 35-44 cm, ML	5/bg	19
ST3544ML-20	BP cuff, SoftCheck Fabric, 1T, lg adult long, 35-44 cm, ML	20/bx	19
ST3544MQ-05	BP cuff, SoftCheck Fabric, 1T, lg adult long, 35-44 cm, MQ	5/bg	20
ST3544MQ-20	BP cuff, SoftCheck Fabric, 1T, lg adult long, 35-44 cm, MQ	20/bx	20
ST3544MQML-05	BP cuff, SoftCheck Fabric, 1T, lg adult long, 35-44 cm, MQML	5/bg	21
ST4250-05	BP cuff, SoftCheck Fabric, 1T, thigh, 42-50 cm, no fitting	5/bg	21
ST4250DM-05	BP cuff, SoftCheck Fabric, 1T, thigh, 42-50 cm, DM	5/bg	20
ST4250FL-05	BP cuff, SoftCheck Fabric, 1T, thigh, 42-50 cm, FL	5/bg	21
ST4250HP-05	BP cuff, SoftCheck Fabric, 1T, thigh, 42-50 cm, HP	5/bg	19
ST4250HP-20	BP cuff, SoftCheck Fabric, 1T, thigh, 42-50 cm, HP	20/bx	19
ST4250ML-05	BP cuff, SoftCheck Fabric, 1T, thigh, 42-50 cm, ML	5/bg	19
ST4250ML-20	BP cuff, SoftCheck Fabric, 1T, thigh, 42-50 cm, ML	20/bx	19
ST4250MQ-05	BP cuff, SoftCheck Fabric, 1T, thigh, 42-50 cm, MQ	5/bg	20
ST4250MQ-20	BP cuff, SoftCheck Fabric, 1T, thigh, 42-50 cm, MQ	20/bx	20
ST4250MQML-05	BP cuff, SoftCheck Fabric, 1T, thigh, 42-50 cm, MQML	5/bg	21

Order Code	Description	UOM	Page(s)
ST7557	Magill breathing circuit w/face mask	ea	103
UD0814-05	BP cuff, UltraCheck Nylon, 2T, infant, 8-14 cm, no fitting	5/bg	17
UD0814DM-05	BP cuff, UltraCheck Nylon, 2T, infant, 8-14 cm, DM	5/bg	16
UD0814MB-05	BP cuff, UltraCheck Nylon, 2T, infant, 8-14 cm, manual	5/bg	17
UD0814MQ-05	BP cuff, UltraCheck Nylon, 2T, infant, 8-14 cm, MQ	5/bg	16
UD0814MQMF-05	BP cuff, UltraCheck Nylon, 2T, infant, 8-14 cm, MQMF	5/bg	16
UD1320-05	BP cuff, UltraCheck Nylon, 2T, child, 13-20 cm, no fitting	5/bg	17
UD1320DM-05	BP cuff, UltraCheck Nylon, 2T, child, 13-20 cm, DM	5/bg	16
UD1320MB-05	BP cuff, UltraCheck Nylon, 2T, child, 13-20 cm, manual	5/bg	17
UD1320MQ-05	BP cuff, UltraCheck Nylon, 2T, child, 13-20 cm, MQ	5/bg	16
UD1320MQMF-05	BP cuff, UltraCheck Nylon, 2T, child, 13-20 cm, MQMF	5/bg	16
UD1826-05	BP cuff, UltraCheck Nylon, 2T, sm adult, 18-26 cm, no fitting	5/bg	17
UD1826DM-05	BP cuff, UltraCheck Nylon, 2T, sm adult, 18-26 cm, DM	5/bg	16
UD1826MB-05	BP cuff, UltraCheck Nylon, 2T, sm adult, 18-26 cm, manual	5/bg	17
UD1826MQ-05	BP cuff, UltraCheck Nylon, 2T, sm adult, 18-26 cm, MQ	5/bg	16
UD1826MQMF-05	BP cuff, UltraCheck Nylon, 2T, sm adult, 18-26 cm, MQMF	5/bg	16
UD2635-05	BP cuff, UltraCheck Nylon, 2T, adult, 26-35 cm, no fitting	5/bg	17
UD2635DM-05	BP cuff, UltraCheck Nylon, 2T, adult, 26-35 cm, 2T, DM	5/bg	16
UD2635MB-05	BP cuff, UltraCheck Nylon, 2T, adult, 26-35 cm, 2T, manual	5/bg	17
UD2635MQ-05	BP cuff, UltraCheck Nylon, 2T, adult, 26-35 cm, 2T, MQ	5/bg	16
UD2635MQMF-05	BP cuff, UltraCheck Nylon, 2T, adult, 26-35 cm, 2T, MQMF	5/bg	16
UD2938-05	BP cuff, UltraCheck Nylon, 2T, adult long, 29-39 cm, 2T, no fitting	5/bg	17
UD2938DM-05	BP cuff, UltraCheck Nylon, 2T, adult long, 29-39 cm, 2T, DM	5/bg	16
UD2938MB-05	BP cuff, UltraCheck Nylon, 2T, adult long, 29-39 cm, 2T, manual	5/bg	17
UD2938MQ-05	BP cuff, UltraCheck Nylon, 2T, adult long, 29-39 cm, 2T, MQ	5/bg	16
UD2938MQMF-05	BP cuff, UltraCheck Nylon, 2T, adult long, 29-39 cm, 2T, MQMF	5/bg	16
UD3242-05	BP cuff, UltraCheck Nylon, 2T, lg adult, 32-42 cm, no fitting	5/bg	17
UD3242DM-05	BP cuff, UltraCheck Nylon, 2T, lg adult, 32-42 cm, DM	5/bg	16
UD3242MB-05	BP cuff, UltraCheck Nylon, 2T, lg adult, 32-42 cm, manual	5/bg	17
UD3242MQ-05	BP cuff, UltraCheck Nylon, 2T, lg adult, 32-42 cm, MQ	5/bg	16
UD3242MQMF-05	BP cuff, UltraCheck Nylon, 2T, lg adult, 32-42 cm, MQMF	5/bg	16
UD3544-05	BP cuff, UltraCheck Nylon, 2T, lg adult long, 35-44 cm, no fitting	5/bg	17
UD3544DM-05	BP cuff, UltraCheck Nylon, 2T, lg adult long, 35-44 cm, DM	5/bg	16
UD3544MB-05	BP cuff, UltraCheck Nylon, 2T, lg adult long, 35-44 cm, manual	5/bg	17
UD3544MQ-05	BP cuff, UltraCheck Nylon, 2T, lg adult long, 35-44 cm, MQ	5/bg	16
UD3544MQMF-05	BP cuff, UltraCheck Nylon, 2T, lg adult long, 35-44 cm, MQMF	5/bg	16
UD4250-05	BP cuff, UltraCheck Nylon, 2T, thigh, 42-50 cm, no fitting	5/bg	17
UD4250DM-05	BP cuff, UltraCheck Nylon, 2T, thigh, 42-50 cm, DM	5/bg	16
UD4250MB-05	BP cuff, UltraCheck Nylon, 2T, thigh, 42-50 cm, manual	5/bg	17
UD4250MQ-05	BP cuff, UltraCheck Nylon, 2T, thigh, 42-50 cm, MQ	5/bg	16
UD4250MQMF-05	BP cuff, UltraCheck Nylon, 2T, thigh, 42-50 cm, MQMF	5/bg	16
US0814-05	BP cuff, UltraCheck Nylon, 1T, infant, 8-14 cm, no fitting	5/bg	15
US0814DM-05	BP cuff, UltraCheck Nylon, 1T, infant, 8-14 cm, DM	5/bg	13
US0814FL-05	BP cuff, UltraCheck Nylon, 1T, infant, 8-14 cm, FL	5/bg	14
US0814HP	BP cuff, UltraCheck Nylon, 1T, infant, 8-14 cm, HP	ea	12
US0814HP-05	BP cuff, UltraCheck Nylon, 1T, infant, 8-14 cm, HP	5/bg	12
US0814MB	BP cuff, UltraCheck Nylon, 1T, infant, 8-14 cm, manual	ea	15
US0814MHP-05	BP cuff, UltraCheck Nylon, 1T, infant, 8-14 cm, MHP	5/bg	13
US0814ML-05	BP cuff, UltraCheck Nylon, 1T, infant, 8-14 cm, ML	5/bg	13
US0814MQ-05	BP cuff, UltraCheck Nylon, 1T, infant, 8-14 cm, MQ	5/bg	14
US0814MQML-05	BP cuff, UltraCheck Nylon, 1T, infant, 8-14 cm, MQML	5/bg	14
US1320-05	BP cuff, UltraCheck Nylon, 1T, child, 13-20 cm, no fitting	5/bg	15
US1320DM-05	BP cuff, UltraCheck Nylon, 1T, child, 13-20 cm, DM	5/bg	13

Order Code	Description	UOM	Page(s)
US1320FL-05	BP cuff, UltraCheck Nylon, 1T, child, 13-20 cm, FL	5/bg	14
US1320HP	BP cuff, UltraCheck Nylon, 1T, child, 13-20 cm, HP	ea	12
US1320HP-05	BP cuff, UltraCheck Nylon, 1T, child, 13-20 cm, HP	5/bg	12
US1320MB	BP cuff, UltraCheck Nylon, 1T, child, 13-20 cm, manual with sphyg	ea	15
US1320MHP-05	BP cuff, UltraCheck Nylon, 1T, child, 13-20 cm, MHP	5/bg	13
US1320ML-05	BP cuff, UltraCheck Nylon, 1T, child, 13-20 cm, ML	5/bg	13
US1320MQ-05	BP cuff, UltraCheck Nylon, 1T, child, 13-20 cm, MQ	5/bg	14
US1320MQML-05	BP cuff, UltraCheck Nylon, 1T, child, 13-20 cm, MQML	5/bg	14
US1826-05	BP cuff, UltraCheck Nylon, 1T, sm adult, 18-26 cm, 1T, no fitting	5/bg	15
US1826DM-05	BP cuff, UltraCheck Nylon, 1T, sm adult, 18-26 cm, 1T, DM	5/bg	13
US1826FL-05	BP cuff, UltraCheck Nylon, 1T, sm adult, 18-26 cm, 1T, FL	5/bg	14
US1826HP	BP cuff, UltraCheck Nylon, 1T, sm adult, 18-26 cm, 1T, HP	ea	12
US1826HP-05	BP cuff, UltraCheck Nylon, 1T, sm adult, 18-26 cm, 1T, HP	5/bg	12
US1826MB	BP cuff, UltraCheck Nylon, 1T, sm adult, 18-26 cm, 1T, manual	ea	15
US1826MHP-05	BP cuff, UltraCheck Nylon, 1T, sm adult, 18-26 cm, 1T, MHP	5/bg	13
US1826ML-05	BP cuff, UltraCheck Nylon, 1T, sm adult, 18-26 cm, 1T, ML	5/bg	13
US1826MQ-05	BP cuff, UltraCheck Nylon, 1T, sm adult, 18-26 cm, 1T, MQ	5/bg	14
US1826MQML-05	BP cuff, UltraCheck Nylon, 1T, sm adult, 18-26 cm, 1T, MQML	5/bg	14
US2635-05	BP cuff, UltraCheck Nylon, 1T, adult, 26-35 cm, no fitting	5/bg	15
US2635DM-05	BP cuff, UltraCheck Nylon, 1T, adult, 26-35 cm, DM	5/bg	13
US2635FL-05	BP cuff, UltraCheck Nylon, 1T, adult, 26-35 cm, FL	5/bg	14
US2635HP	BP cuff, UltraCheck Nylon, 1T, adult, 26-35 cm, HP	ea	12
US2635HP-05	BP cuff, UltraCheck Nylon, 1T, adult, 26-35 cm, HP	5/bg	12
US2635MB	BP cuff, UltraCheck Nylon, 1T, adult, 26-35 cm, manual with sphyg	ea	15
US2635MHP-05	BP cuff, UltraCheck Nylon, 1T, adult, 26-35 cm, MHP	5/bg	13
US2635ML-05	BP cuff, UltraCheck Nylon, 1T, adult, 26-35 cm, ML	5/bg	13
US2635MQ-05	BP cuff, UltraCheck Nylon, 1T, adult, 26-35 cm, MQ	5/bg	14
US2635MQML-05	BP cuff, UltraCheck Nylon, 1T, adult, 26-35 cm, MQML	5/bg	14
US2938-05	BP cuff, UltraCheck Nylon, 1T, adult long, 29-38 cm, no fitting	5/bg	15
US2938DM-05	BP cuff, UltraCheck Nylon, 1T, adult long, 29-38 cm, DM	5/bg	13
US2938FL-05	BP cuff, UltraCheck Nylon, 1T, adult long, 29-38 cm, FL	5/bg	14
US2938HP	BP cuff, UltraCheck Nylon, 1T, adult long, 29-38 cm, HP	ea	12
US2938HP-05	BP cuff, UltraCheck Nylon, 1T, adult long, 29-38 cm, HP	5/bg	12
US2938MB	BP cuff, UltraCheck Nylon, 1T, adult long, 29-38 cm, manual with sphyg	ea	15
US2938MHP-05	BP cuff, UltraCheck Nylon, 1T, adult long, 29-38 cm, 1T, MHP	5/bg	13
US2938ML-05	BP cuff, UltraCheck Nylon, 1T, adult long, 29-38 cm, ML	5/bg	13
US2938MQ-05	BP cuff, UltraCheck Nylon, 1T, adult long, 29-38 cm, MQ	5/bg	14
US2938MQML-05	BP cuff, UltraCheck Nylon, 1T, adult long, 29-38 cm, MQML	5/bg	14
US3242-05	BP cuff, UltraCheck Nylon, 1T, lg adult, 32-42 cm, no fitting	5/bg	15
US3242DM-05	BP cuff, UltraCheck Nylon, 1T, lg adult, 32-42 cm, DM	5/bg	13
US3242FL-05	BP cuff, UltraCheck Nylon, 1T, lg adult, 32-42 cm, FL	5/bg	14
US3242HP	BP cuff, UltraCheck Nylon, 1T, lg adult, 32-42 cm, HP	ea	12
US3242HP-05	BP cuff, UltraCheck Nylon, 1T, lg adult, 32-42 cm, HP	5/bg	12
US3242MB	BP cuff, UltraCheck Nylon, 1T, lg adult, 32-42 cm, manual with sphyg	ea	15
US3242MHP-05	BP cuff, UltraCheck Nylon, 1T, lg adult, 32-42 cm, MHP	5/bg	13
US3242ML-05	BP cuff, UltraCheck Nylon, 1T, lg adult, 32-42 cm, ML	5/bg	13
US3242MQ-05	BP cuff, UltraCheck Nylon, 1T, lg adult, 32-42 cm, MQ	5/bg	14
US3242MQML-05	BP cuff, UltraCheck Nylon, 1T, lg adult, 32-42 cm, MQML	5/bg	14
US3544-05	BP cuff, UltraCheck Nylon, 1T, lg adult long, 29-39 cm, no fitting	5/bg	15
US3544DM-05	BP cuff, UltraCheck Nylon, 1T, sm adult, 18-26 cm, 1T, DM	5/bg	13
US3544FL-05	BP cuff, UltraCheck Nylon, 1T, lg adult long, 29-39 cm, FL	5/bg	14
US3544HP	BP cuff, UltraCheck Nylon, 1T, lg adult long, 29-39 cm, HP	ea	12
US3544HP-05	BP cuff, UltraCheck Nylon, 1T, lg adult long, 29-39 cm, HP	5/bg	12

Order Code	Description	UOM	Page(s)
US3544MB	BP cuff, UltraCheck Nylon, 1T, lg adult long, 29-39 cm, manual	ea	15
US3544MHP-05	BP cuff, UltraCheck Nylon, 1T, lg adult long, 35-44 cm, MHP	5/bg	13
US3544ML-05	BP cuff, UltraCheck Nylon, 1T, lg adult long, 29-39 cm, ML	5/bg	13
US3544MQ-05	BP cuff, UltraCheck Nylon, 1T, lg adult long, 29-39 cm, MQ	5/bg	14
US3544MQML-05	BP cuff, UltraCheck Nylon, 1T, lg adult long, 29-39 cm, MQML	5/bg	14
US4250-05	BP cuff, UltraCheck Nylon, 1T, thigh, 42-50 cm, no fitting	5/bg	15
US4250DM-05	BP cuff, UltraCheck Nylon, 1T, thigh, 42-50 cm, DM	5/bg	13
US4250FL-05	BP cuff, UltraCheck Nylon, 1T, thigh, 42-50 cm, FL	5/bg	14
US4250HP	BP cuff, UltraCheck Nylon, 1T, thigh, 42-50 cm, HP	ea	12
US4250HP-05	BP cuff, UltraCheck Nylon, 1T, thigh, 42-50 cm, HP	5/bg	12
US4250MB	BP cuff, UltraCheck Nylon, 1T, thigh, 42-50 cm, manual	ea	15
US4250MHP-05	BP cuff, UltraCheck Nylon, 1T, thigh, 42-50 cm, MHP	5/bg	13
US4250ML-05	BP cuff, UltraCheck Nylon, 1T, thigh, 42-50 cm, ML	5/bg	13
US4250MQ-05	BP cuff, UltraCheck Nylon, 1T, thigh, 42-50 cm, MQ	5/bg	14
US4250MQML-05	BP cuff, UltraCheck Nylon, 1T, thigh, 42-50 cm, MQML	5/bg	14
VDT0814-20	BP cuff, SoftCheck Vinyl, 2T, infant, 8-14 cm, no fitting	20/bx	29
VDT0814DM-05	BP cuff, SoftCheck Vinyl, 2T, infant, 8-14 cm, DM	5/bg	29
VDT1320-20	BP cuff, SoftCheck Vinyl, 2T, child, 13-20 cm, no fitting	20/bx	29
VDT1320DM-05	BP cuff, SoftCheck Vinyl, 2T, child, 13-20 cm, DM	5/bg	29
VDT1826-05	BP cuff, SoftCheck Vinyl, 2T, sm adult, 18-26 cm, no fitting	5/bg	29
VDT1826DM-05	BP cuff, SoftCheck Vinyl, 2T, sm adult, 18-26 cm, DM	5/bg	29
VDT2635-20	BP cuff, SoftCheck Vinyl, 2T, adult, 26-35 cm, no fitting	20/bx	29
VDT2635DM-05	BP cuff, SoftCheck Vinyl, 2T, adult, 26-35 cm, DM	5/bg	29
VDT2938-20	BP cuff, SoftCheck Vinyl, 2T, adult long, 29-38 cm, no fitting	20/bx	29
VDT2938DM-05	BP cuff, SoftCheck Vinyl, 2T, adult long, 29-38 cm, DM	5/bg	29
VDT3242-20	BP cuff, SoftCheck Vinyl, 2T, lg adult, 32-42 cm, no fitting	20/bx	29
VDT3242DM-05	BP cuff, SoftCheck Vinyl, 2T, lg adult, 32-42 cm, DM	5/bg	29
VDT4250-20	BP cuff, SoftCheck Vinyl, 2T, thigh, 42-50 cm, no fitting	20/bx	29
VDT4250DM-05	BP cuff, SoftCheck Vinyl, 2T, thigh, 42-50 cm, DM	5/bg	29
VNN1DT-10	BP cuff, SoftCheck Vinyl Neonatal, 2T, size 1, 3-6 cm, ML	10/bx	33
VNN1ST-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 1, 3-6 cm, ML	10/bx	32
VNN1ST-HP-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 1, 3-6 cm, HP	10/bx	32
VNN1ST-M316-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 1, 3-6 cm, M316	10/bx	32
VNN1ST-MQ-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 1, 3-6 cm, MQ	10/bx	32
VNN1ST-SQ-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 1, 3-6 cm, SQ	10/bx	32
VNN2DT-10	BP cuff, SoftCheck Vinyl Neonatal, 2T, size 2, 4-8 cm, ML	10/bx	33
VNN2ST-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 2, 4-8 cm, ML	10/bx	32
VNN2ST-HP-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 2, 4-8 cm, HP	10/bx	32
VNN2ST-M316-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 2, 4-8 cm, M316	10/bx	32
VNN2ST-MQ-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 2, 4-8 cm, MQ	10/bx	32
VNN2ST-SQ-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 2, 4-8 cm, SQ	10/bx	32
VNN3DT-10	BP cuff, SoftCheck Vinyl Neonatal, 2T, size 3, 6-11 cm, ML	10/bx	33
VNN3ST-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 3, 6-11 cm, ML	10/bx	32
VNN3ST-HP-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 3, 6-11 cm, HP	10/bx	32
VNN3ST-M316-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 3, 6-11 cm, M316	10/bx	32
VNN3ST-MQ-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 3, 6-11 cm, MQ	10/bx	32
VNN3ST-SQ-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 3, 6-11 cm, SQ	10/bx	32
VNN4DT-10	BP cuff, SoftCheck Vinyl Neonatal, 2T, size 4, 7-13 cm, ML	10/bx	33
VNN4ST-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 4, 7-13 cm, ML	10/bx	32
VNN4ST-HP-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 4, 7-13 cm, HP	10/bx	32
VNN4ST-M316-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 4, 7-13 cm, M316	10/bx	32
VNN4ST-MQ-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 4, 7-13 cm, MQ	10/bx	32
VNN4ST-SQ-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 4, 7-13 cm, SQ	10/bx	32

Order Code	Description	UOM	Page(s)
VNN5DT-10	BP cuff, SoftCheck Vinyl Neonatal, 2T, size 5, 8-15 cm, ML	10/bx	33
VNN5ST-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 5, 8-15 cm, ML	10/bx	32
VNN5ST-HP-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 5, 8-15 cm, HP	10/bx	32
VNN5ST-M316-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 5, 8-15 cm, M316	10/bx	32
VNN5ST-MQ-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 5, 8-15 cm, MQ	10/bx	32
VNN5ST-SQ-10	BP cuff, SoftCheck Vinyl Neonatal, 2T, size 5, 8-15 cm, SQ	10/bx	32
VST0814HP-05	BP cuff, SoftCheck Vinyl, 1T, infant, 8-14 cm, HP	5/bg	26
VST0814HP-20	BP cuff, SoftCheck Vinyl, 1T, infant, 8-14 cm, HP	20/bx	26
VST0814ML-05	BP cuff, SoftCheck Vinyl, 1T, infant, 8-14 cm, ML	5/bg	27
VST0814ML-20	BP cuff, SoftCheck Vinyl, 1T, infant, 8-14 cm, ML	20/bx	27
VST0814MQ-05	BP cuff, SoftCheck Vinyl, 1T, infant, 8-14 cm, MQ	5/bg	28
VST0814MQ-20	BP cuff, SoftCheck Vinyl, 1T, infant, 8-14 cm, MQ	20/bx	28
VST0814MQML-20	BP cuff, SoftCheck Vinyl, 1T, infant, 8-14 cm, MQML	20/bx	28
VST1320HP-05	BP cuff, SoftCheck Vinyl, 1T, child, 13-20 cm, HP	5/bg	26
VST1320HP-20	BP cuff, SoftCheck Vinyl, 1T, child, 13-20 cm, HP	20/bx	26
VST1320ML-05	BP cuff, SoftCheck Vinyl, 1T, child, 13-20 cm, ML	5/bg	27
VST1320ML-20	BP cuff, SoftCheck Vinyl, 1T, child, 13-20 cm, ML	20/bx	27
VST1320MQ-05	BP cuff, SoftCheck Vinyl, 1T, child, 13-20 cm, MQ	5/bg	28
VST1320MQ-20	BP cuff, SoftCheck Vinyl, 1T, child, 13-20 cm, MQ	20/bx	28
VST1320MQML-20	BP cuff, SoftCheck Vinyl, 1T, child, 13-20 cm, MQML	20/bx	28
VST1826HP-05	BP cuff, SoftCheck Vinyl, 1T, sm adult, 18-26 cm, HP	5/bg	26
VST1826HP-20	BP cuff, SoftCheck Vinyl, 1T, sm adult, 18-26 cm, HP	20/bx	26
VST1826ML-05	BP cuff, SoftCheck Vinyl, 1T, sm adult, 18-26 cm, ML	5/bg	27
VST1826ML-20	BP cuff, SoftCheck Vinyl, 1T, sm adult, 18-26 cm, ML	20/bx	27
VST1826MQ-05	BP cuff, SoftCheck Vinyl, 1T, sm adult, 18-26 cm, MQ	5/bg	28
VST1826MQML-20	BP cuff, SoftCheck Vinyl, 1T, sm adult, 18-26 cm, MQML	20/bx	28
VST2635HP-05	BP cuff, SoftCheck Vinyl, 1T, adult, 26-35 cm, HP	5/bg	26
VST2635HP-20	BP cuff, SoftCheck Vinyl, 1T, adult, 26-35 cm, HP	20/bx	26
VST2635ML-05	BP cuff, SoftCheck Vinyl, 1T, adult, 26-35 cm, ML	5/bg	27
VST2635ML-20	BP cuff, SoftCheck Vinyl, 1T, adult, 26-35 cm, ML	20/bx	27
VST2635MQ-05	BP cuff, SoftCheck Vinyl, 1T, adult, 26-35 cm, MQ	5/bg	28
VST2635MQ-20	BP cuff, SoftCheck Vinyl, 1T, adult, 26-35 cm, MQ	20/bx	28
VST2635MQML-20	BP cuff, SoftCheck Vinyl, 1T, adult, 25-35 cm, MQML	20/bx	28
VST2938HP-05	BP cuff, SoftCheck Vinyl, 1T, adult long, 29-38 cm, HP	5/bg	26
VST2938HP-20	BP cuff, SoftCheck Vinyl, 1T, adult long, 29-38 cm, HP	20/bx	26
VST2938ML-05	BP cuff, SoftCheck Vinyl, 1T, adult long, 29-38 cm, ML	5/bg	27
VST2938ML-20	BP cuff, SoftCheck Vinyl, 1T, adult long, 29-38 cm, ML	20/bx	27
VST2938MQ-05	BP cuff, SoftCheck Vinyl, 1T, adult long, 29-38 cm, MQ	5/bg	28
VST2938MQ-20	BP cuff, SoftCheck Vinyl, 1T, adult long, 29-38 cm, MQ	20/bx	28
VST2938MQML-20	BP cuff, SoftCheck Vinyl, 1T, adult long, 29-38 cm, MQML	20/bx	28
VST3242HP-05	BP cuff, SoftCheck Vinyl, 1T, lg adult, 32-42 cm, HP	5/bg	26
VST3242HP-20	BP cuff, SoftCheck Vinyl, 1T, lg adult, 32-42 cm, HP	20/bx	26
VST3242ML-05	BP cuff, SoftCheck Vinyl, 1T, lg adult, 32-42 cm, ML	5/bg	27
VST3242ML-20	BP cuff, SoftCheck Vinyl, 1T, lg adult, 32-42 cm, ML	20/bx	27
VST3242MQ-05	BP cuff, SoftCheck Vinyl, 1T, lg adult, 32-42 cm, MQ	5/bg	28
VST3242MQ-20	BP cuff, SoftCheck Vinyl, 1T, lg adult, 32-42 cm, MQ	20/bx	28
VST3242MQML-20	BP cuff, SoftCheck Vinyl, 1T, lg adult, 32-42 cm, MQML	20/bx	28
VST3544HP-05	BP cuff, SoftCheck Vinyl, 1T, lg adult long, 35-44 cm, HP	5/bg	26
VST3544HP-20	BP cuff, SoftCheck Vinyl, 1T, lg adult long, 35-44 cm, HP	20/bx	26
VST3544ML-05	BP cuff, SoftCheck Vinyl, 1T, lg adult long, 35-44 cm, ML	5/bg	27
VST3544ML-20	BP cuff, SoftCheck Vinyl, 1T, lg adult long, 35-44 cm, ML	20/bx	27
VST3544MQML-20	BP cuff, SoftCheck Vinyl, 1T, lg adult long, 35-44 cm, MQML	20/bx	28
VST4250HP-05	BP cuff, SoftCheck Vinyl, 1T, thigh, 42-50 cm, HP	5/bg	26

Order Code	Description	UOM	Page(s)
VST4250HP-20	BP cuff, SoftCheck Vinyl, 1T, thigh, 42-50 cm, HP	20/bx	26
VST4250ML-05	BP cuff, SoftCheck Vinyl, 1T, thigh, 42-50 cm, ML	5/bg	27
VST4250ML-20	BP cuff, SoftCheck Vinyl, 1T, thigh, 42-50 cm, ML	20/bx	27
VST4250MQ-05	BP cuff, SoftCheck Vinyl, 1T, thigh, 42-50 cm, MQ	5/bg	28
VST4250MQ-20	BP cuff, SoftCheck Vinyl, 1T, thigh, 42-50 cm, MQ	20/bx	28
VST4250MQML-20	BP cuff, SoftCheck Vinyl, 1T, thigh, 42-50 cm, MQML	20/bx	28
YVST1826HP-05	BP cuff, SoftCheck Yellow Vinyl, 1T, sm adult, 18-26 cm, HP	5/bg	30
YVST2635HP-05	BP cuff, SoftCheck Yellow Vinyl, 1T, adult, 26-35 cm, HP	5/bg	30
YVST3242HP-05	BP cuff, SoftCheck Yellow Vinyl, 1T, lg adult, 32-42 cm, HP	5/bg	30
ZFT-BPF230	BP fitting, BPF230	ea	35
ZFT-BPF240	BP fitting, BP240	ea	35
ZFT-CS901	BP fitting, CS901	ea	36
ZFT-DM-20	BP fitting, DM	20/bx	36
ZFT-DM2YRX1	BP Y adapter, DM2YRX1	ea	35
ZFT-DM532-20	BP fitting, DM532	20/bx	36
ZFT-DMML-20	BP fitting, DMML	20/bx	36
ZFT-FL532-20	BP fitting, FL532	20/bx	36
ZFT-FLL240-20	BP fitting, FLL240	20/bx	36
ZFT-HP-20	BP fitting, HP	20/bx	35
ZFT-HP532-20	BP fitting, HP532	20/bx	35
ZFT-HPDMAD-05	BP Y adapter, HPDMAD	5/bg	36
ZFT-ML2YFL	BP Y adapter, ML2YFL	ea	35
ZFT-ML532-20	BP fitting, ML532	20/bx	36
ZFT-MLL532-20	Male Locking Luer	20/bx	36
ZFT-MQ-20	BP fitting, MQ	20/bx	36
ZFT-MQ2YRX1	BP Y adapter, MQ2YRX1	ea	35
ZFT-MQMF2YHP	BP Y adapter, MQMF2YHP	ea	36
ZFT-MQMF2YRX1	BP Y adapter, MQMF2YRX1	ea	35
ZFT-MQML-20	BP fitting, MQML	20/bx	36
ZFT-MQML2YRX1	BP Y adapter, MQML2YRX1	ea	35
ZFT-MQML5-20	BP fitting, MQML5	20/bx	36
ZFT-QC8030	BP fitting, QC8030	ea	36
ZFT-QC8220	BP fitting, QC8220	ea	35
ZFT-QC8230	BP fitting, QC8230	ea	35
ZFT-RX1DM	BP adapter, RX1DM	ea	35
ZFT-RX1FL	BP adapter, RX1FL	ea	35
ZFT-RX1ML	5 in Tube w/Metal Female Bayonet and Male Luer	ea	35
ZFT-RX1MQ	BP adapter, RX1MQ	ea	35
ZFT-RX1MQML	BP adapter, RX1MQML	ea	35
ZGA-BPGBHB	BP cuff Inflator w/gauge, BPGHB	ea	36
ZGA-BPHBRX1	BP cuff Inflator w/gauge, BPHBA1241	ea	36
ZHB-4FBPADPRX1	BP manual bulb adapter, 4FBPADPRX1	ea	36
ZHB-BPCOILADPRX1	BP manual bulb adapter, BPCOILADPRX1	ea	35
ZHB-HBTHM	Hand bulb, THM	ea	36