



## SIDESTREAM SELECTION GUIDE FOR 91517 CAPNOGRAPHY MODULES

IS THE PATIENT INTUBATED?

YES

NO

**SHORT TERM**  
Up to 12 hours  
(patient non-humidified)

**LONG TERM**  
Up to 120 hours  
(patient humidified)

**Key Applications:**  
Critical Care, Operating Room, Recovery Room

**Key Applications:**  
Critical Care, Operating Room



704-0023-00  
Adult/Pediatric  
(ET Tube >4.0mm)  
704-0024-00  
Pediatric/Infant  
(ET Tube <4.0mm)



704-0027-00  
Sample line with Luer  
(Recommended for Heat  
Moisture Exchange Filters)



704-0025-00  
Adult/Pediatric  
(ET Tube >4.0mm;  
with Nafion)



704-0026-00  
Pediatric/Infant  
(ET Tube <4.0mm;  
with Nafion)



704-0028-00  
Sample line with Luer  
and Nafion  
(Recommended for Heat  
Moisture Exchange Filters)

DOES THE PATIENT NEED OXYGEN?

YES

NO

Monitoring O<sub>2</sub>/CO<sub>2</sub>  
Nasal  
Up to 12 hours

Monitoring O<sub>2</sub>/CO<sub>2</sub>  
Oral/Nasal  
Up to 12 hours

Nasal  
Up to 12 hours

Oral/Nasal  
Up to 12 hours

**Key Applications:**  
ED, EMS, Transport, Conscious Sedation, Critical Care

**Key Applications:**  
Sleep Lab, ED, EMS, Transport,  
Conscious Sedation, Critical Care



704-0014-00 Adult  
704-0016-00 Pediatric  
704-0015-00 Infant



704-0020-00 Adult  
704-0022-00 Pediatric



704-0010-00 Adult  
704-0012-00 Pediatric  
704-0011-00 Infant



704-0017-00 Adult  
704-0019-00 Pediatric