



CAPNOGRAPHY SELECTION GUIDE

FOR 92516, 92517, 92518
CAPNOGRAPHY & MULTIGAS MODULES



IS THE PATIENT INTUBATED?

YES

NO

SIDESTREAM

MAINSTREAM
(92517 ONLY)

Key Applications:

Emergency Dept, Critical Care, Operating Room, Recovery Room

Key Applications:
Critical Care, NICU



103-0234-00

Multi-Patient Use Nomoline Adapter



010-1980-00

Mainstream Sensor



166-7085-00

6.5ft Sampling Line

166-7121-00

10ft Sampling Line



704-0175-00

Adult/Ped T-Adapter



704-0173-00

Airway Adapter, Adult/Ped



704-0174-00

Airway Adapter, Infant

DOES THE PATIENT NEED OXYGEN?

YES

NO

Monitoring O₂/CO₂
Nasal Only

Monitoring O₂/CO₂
Oral/Nasal

CO₂ Nasal
Cannula

Key Applications:

Emergency Dept, Transport, Procedural Sedation,
Critical Care, Recovery, General Floors

Key Applications:

Emergency Dept, Transport, Procedural Sedation,
Critical Care, General Floors



103-0234-00

Multi-Patient Use Nomoline Adapter



103-0234-00

Multi-Patient Use Nomoline Adapter

704-0181-00 Adult
704-0180-00 Pediatric
704-0179-00 Infant



704-0183-00 Adult

704-0182-00 Pediatric



704-0176-00 Adult

704-0177-00 Pediatric

704-0178-00 Infant

Spacelabs Healthcare reserves the right to make changes to the products at any time without notice. All brands and product names are trademarks of their respective owners. Product is CE Marked. Product may not be available for sale in all locations. Check with your local Spacelabs representative.