




10630 Wiles Road, Coral Springs, FL 33076 • toll free: 800.358.2468 • office: 954.345.9800 • fax: 954.345.0507 • http://www.advantagemed.com

## PHILIPS® QUICK REFERENCE GUIDE

**MOST COMMONLY PURCHASED PHILIPS® ACCESSORIES SHOWN. FOR COMPLETE LISTING OF PRODUCTS CALL OR VISIT OUR WEBSITE.**

### YOKEMATE LWS® (UNIVERSAL LEAD WIRE SYSTEM)

### CMS PATIENT CABLES



#### YOKEMATE LWS® ADAPTERS

- 3 Lead - YM-813DIN
- 5 Lead - YM-815DIN



#### TELEMETRY YOKEMATE LWS®

- 3 / 5 Lead - YM-815DIN-TM

Allows for AMC&E's Direct Connect IntelliVue lead wires to be used with 2601A transmitters.



#### UNIVERSAL CABLES\*

- 3 Lead - CB-73385R (M1669A)
- 5 Lead - CB-73585R (M1668A)
- 6 Lead - CB-73685R (M1667A)
- 10 Lead (5+5) - CB-735585R (M1663A)
- 10 Lead (6+4) - CB-736485R (M1665A)



- 3 Lead - CB-81385R (M1500A)
- 5 Lead - CB-81585R (M1520A)

\* Universal cables work with AMC&E's Philips® IntelliVue direct connect lead wires, but are not compatible with original manufacturer's Philips® IntelliVue lead wires.


### COMMON PHILIPS® DIRECT CONNECT LEAD WIRES (HARDWIRED AND TELEMETRY)

**PHILIPS® OEM LEAD WIRES ARE 40" OR 62" – AMC&E LEAD WIRES ARE 29" OR 51"**



SNAP, PINCH, & MINI PINCH\*

\* Mini Pinch Available by Special Order Only



HARDWIRED


TELEMETRY

#### PHILIPS® CMS

SNAP		
<b>3 Lead</b>		
▶ LW-2300029/3A	29 inch	M1605A
▶ LW-2300051/3A	51 inch	
▶ LW-2300029/3AT	29 inch [Telemetry]	M2590A
<b>5 Lead</b>		
▶ LW-2300029/5A	29 inch	
▶ LW-23000MX/5A	Mixed Length	M1625A
▶ LW-2300051/5A	51 inch	
▶ LW-2300029/5AT	29 inch [Telemetry]	M2592A

PINCH		
<b>3 Lead</b>		
▶ LW-3300029/3A	29 inch	M1603A
▶ LW-3300051/3A	51 inch	
▶ LW-3300029/3AT	29 inch [Telemetry]	M2591A
<b>5 Lead</b>		
▶ LW-3300029/5A	29 inch	
▶ LW-33000MX/5A	Mixed Length	M1623A
▶ LW-3300051/5A	51 inch	
▶ LW-3300029/5AT	29 inch [Telemetry]	M2593A


#### PHILIPS® INTELLIVUE



REUSABLE

SNAP		
<b>3 Lead</b>		
▶ LW-2090S29/3A	29 inch	M1673A
▶ LW-2090S51/3A	51 inch	
<b>5 Lead</b>		
▶ LW-2090S29/5A	29 inch	
▶ LW-2090SMX/5A	Mixed Length	M1644A
▶ LW-2090S51/5A	51 inch	
▶ LW-2090S29/5V	29 inch [V-Leads]	M1602A
▶ LW-2090S29/5AT	29 inch [Telemetry]	9898031 52071

PINCH		
<b>3 Lead</b>		
▶ LW-3090S29/3A	29 inch	M1671A
▶ LW-3090S51/3A	51 inch	
<b>5 Lead</b>		
▶ LW-3090S29/5A	29 inch	
▶ LW-3090SMX/5A	Mixed Length	M1968A
▶ LW-3090S51/5A	51 inch	
▶ LW-3090S29/5V	29 inch [V-Leads]	M1976A
▶ LW-3090S29/5AT	29 inch [Telemetry]	9898031 52051



DISPOSABLE

SNAP		
<b>3 Lead</b>		
▶ LW-209DS50/3A	50 inch	No Equivalent
<b>5 Lead</b>		
▶ LW-209DS50/5A	50 inch	No Equivalent
▶ LW-209DS50/5AT	50 inch [V-Leads]	No Equivalent
▶ LW-209DS50/5V	50 inch [Telemetry]	No Equivalent

PINCH		
<b>3 Lead</b>		
▶ LW-309DS50/3A	50 inch	No Equivalent
<b>5 Lead</b>		
▶ LW-309DS50/5A	50 inch	No Equivalent
▶ LW-309DS50/5AT	50 inch [V-Leads]	No Equivalent
▶ LW-309DS50/5V	50 inch [Telemetry]	No Equivalent

■ AMC&E Part Number

■ Philips® OEM Part Number



10630 Wiles Road, Coral Springs, FL 33076 • toll free: 800.358.2468 • office: 954.345.9800 • fax: 954.345.0507 • http://www.advantagemed.com

## PHILIPS® QUICK REFERENCE GUIDE (CONT.)

**MOST COMMONLY PURCHASED PHILIPS® ACCESSORIES SHOWN. FOR COMPLETE LISTING OF PRODUCTS CALL OR VISIT OUR WEBSITE.**

### CUFFS

For NeoNatal Cuffs Call AMC&E



\* AMC&E cuff sizes differ slightly from Philips OEM. Please refer to sizes of each.

Reusable		
▶ CU-US0814HP	M4552A	1 Cuff
▶ CU-US1320HP	M4553A	1 Cuff
▶ CU-US1826HP	M4554A	1 Cuff
▶ CU-US2635HP	M4555A	1 Cuff
▶ CU-US2345HP	No Equivalent	1 Cuff
▶ CU-US4560HP	M4557A	1 Cuff
▶ CU-US3860HP	M4559A	1 Cuff

- Infant (8-14 cm) -
- Child (13-20 cm) -
- Small Adult (18-26 cm) -
- Adult (26-35 cm) -
- Long Adult (23-45 cm) -
- Extra Long Adult (45-60 cm) -
- Large Adult / Thigh (38-60 cm) -

Disposable		
▶ CU-ST0814HP-20	M1874A	Box of 20 Cuffs
▶ CU-ST1320HP-20	M1875A	Box of 20 Cuffs
▶ CU-ST1826HP-20	M1876A	Box of 20 Cuffs
▶ CU-ST2635HP-20	M1877A	Box of 20 Cuffs
▶ CU-ST2345HP-20	No Equivalent	Box of 20 Cuffs
▶ CU-ST4560HP-10	M1878A	Box of 10 Cuffs
▶ CU-ST3860HP-10	M1879A	Box of 10 Cuffs

### HOSES



▶ HOS16158-10 - M1599B




▶ HO-S16368-10 - No Equivalent

### SPO<sub>2</sub>


Philips® and Nellcor® (Non-OxiMax) Probes

#### SOFT TIP



▶ PR-A900-1006 - M1190A


#### LARGE FINGER CLIP



▶ PR-A520-1006 - No Equivalent




▶ PR-A900-1006V - M1191A, M1191B



▶ PR-A520-1006V - M1196A

#### Extension Cable




▶ CB-A400-1006VM - M1941A


#### Adapter Cable



▶ CB-A400-1006V - M1940A



▶ PR-A900-1011N - M1191T



▶ PR-A520-1011N - M1196T

#### Adapter Cables



▶ CB-A400-1006VN - M1943A



▶ CB-A400-1006A - M1900B