


SUPPLIES & ACCESSORIES ELECTROCARDIOGRAPHY



The ECG Combiner System can standardize Spacelabs monitors and Ultraview & Aria Telemetry transmitters to use one set of lead wires.

Additional cables, adapters and lead wires are also available to standardize most other manufacturer monitors and telemetry transmitters, allowing your hospital to use one set of lead wires per patient throughout the continuum of care.



Reusable

Disposable

Multi-Pinch

Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Combiner Lead Wire Sets - requires ECG Combiner Cable or Adapter

DIN Shielded, Multi-Pinch, Disposable, AAMI



Multi-Pinch

700-0006-30	3 Lead	Multi-Pinch, Disposable	127cm /50 in	25 Sets
700-0006-31	4 Lead (V-Lead)	Multi-Pinch, Disposable	127cm /50 in	25 Sets
700-0006-32	5 Lead	Multi-Pinch, Disposable	127cm /50 in	25 Sets
700-0006-33	5 Lead (V-Lead)	Multi-Pinch, Disposable	127cm /50 in	25 Sets
700-0006-34	6 Lead	Multi-Pinch, Disposable	127cm /50 in	25 Sets

DIN Shielded, Snap, Disposable, AAMI



Snap

700-0007-30	3 Lead	Snap, Disposable	127cm /50 in	25 Sets
700-0007-31	4 Lead (V-Lead)	Snap, Disposable	127cm /50 in	25 Sets
700-0007-32	5 Lead	Snap, Disposable	127cm /50 in	25 Sets
700-0007-33	5 Lead (V-Lead)	Snap, Disposable	127cm /50 in	25 Sets
700-0007-34	6 Lead	Snap, Disposable	127cm /50 in	25 Sets

DIN Shielded, Multi-Pinch, Reusable, AAMI



Multi-Pinch

700-0006-35	3 Lead	Multi-Pinch, Reusable	74cm /29 in	1 Set
700-0006-36	4 Lead (V-Lead)	Multi-Pinch, Reusable	74cm /29 in	1 Set
700-0006-37	5 Lead	Multi-Pinch, Reusable	74cm /29 in	1 Set
700-0006-38	5 Lead (V-Lead)	Multi-Pinch, Reusable	74cm /29 in	1 Set
700-0006-39	6 Lead	Multi-Pinch, Reusable	74cm /29 in	1 Set
700-0010-00	3 Lead	Multi-Pinch, Reusable	130cm/51 in	1 Set
700-0010-01	5 Lead	Multi-Pinch, Reusable	130cm/51 in	1 Set
700-0010-02	5 Lead (V-Lead)	Multi-Pinch, Reusable	130cm/51 in	1 Set

DIN Shielded, Snap, Reusable, AAMI




Snap

700-0007-35	3 Lead	Snap, Reusable	74cm /29 in	1 Set
700-0007-36	4 Lead (V-Lead)	Snap, Reusable	74cm /29 in	1 Set
700-0007-37	5 Lead	Snap, Reusable	74cm /29 in	1 Set
700-0007-38	5 Lead (V-Lead)	Snap, Reusable	74cm /29 in	1 Set
700-0007-39	6 Lead	Snap, Reusable	74cm /29 in	1 Set

The ECG Combiner System can standardize Spacelabs monitors and Ultraview & Aria Telemetry transmitters to use one set of lead wires.

Additional cables, adapters and lead wires are also available to standardize most other manufacturer monitors and telemetry transmitters, allowing your hospital to use one set of lead wires per patient throughout the continuum of care.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Combiner Radiolucent Lead Wire Sets - requires ECG Combiner Cable or Adapter

DIN Shielded, Multi-Pinch, Reusable, AAMI



700-0010-03	3 Lead	Radiolucent, Multi-Pinch, Reusable	74cm /29 in	1 Set
700-0010-04	5 Lead	Radiolucent, Multi-Pinch, Reusable	74cm /29 in	1 Set
700-0010-05	3 Lead	Radiolucent, Multi-Pinch, Reusable	122cm/48 in	1 Set
700-0010-06	5 Lead	Radiolucent, Multi-Pinch, Reusable	122cm/48 in	1 Set
700-0010-07	3 Lead	Radiolucent, Multi-Pinch, Reusable	183cm/72 in	1 Set
700-0010-08	5 Lead	Radiolucent, Multi-Pinch, Reusable	183cm/72 in	1 Set

ECG Combiner Lead Wire Sets with Spacelabs Connection - plugs directly into Command Modules 90496/91496

Shielded, Spacelabs Connection, Multi-Pinch, Reusable, AAMI



700-0010-16	5 Lead	Spacelabs Connection, Multi-Pinch	60 cm/24 in	1 Set
700-0010-14	5 Lead	Spacelabs Connection, Multi-Pinch	102 cm/40 in	1 Set



Spacelabs Connection

Shielded, Spacelabs Connection, Snap, Reusable, AAMI



700-0010-15	5 Lead	Spacelabs Connection, Snap	60 cm/24 in	1 Set
-------------	--------	----------------------------	-------------	-------

The ECG Combiner System can standardize Spacelabs monitors and Ultraview & Aria Telemetry transmitters to use one set of lead wires.

Additional cables, adapters and lead wires are also available to standardize most other manufacturer monitors and telemetry transmitters, allowing your hospital to use one set of lead wires per patient throughout the continuum of care.



Reusable



Disposable



Multi-Pinch



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Combiner Cables

Spacelabs Connect Hardwired Cables for monitors using Command Module (90496/91496)



700-0008-55	3 Lead	AAMI	305 cm/10 ft	1 Cable
700-0008-57	5 Lead	AAMI	305 cm/10 ft	1 Cable
700-0008-59	6 Lead	AAMI	305 cm/10 ft	1 Cable
700-0008-56	10 Lead (5+5)	AAMI	305 cm/10 ft	1 Cable
700-0008-58	10 Lead (6+4)	AAMI	305 cm/10 ft	1 Cable

ECG Combiner Adapters

Spacelabs Monitors, AAMI



700-0008-78	3 Lead	Command Modules (90496/ 91496)		1 Adapter
700-0008-79	5 Lead	Command Modules (90496/ 91496)		1 Adapter


Spacelabs Telemetry, AAMI



700-0008-80	5 Lead	Ultraview Telemetry (9034X/ 9134X)		1 Adapter
700-0008-88	5 Lead	Aria Telemetry (96281)		1 Adapter

The ECG Combiner System can standardize Spacelabs monitors and Ultraview & Aria Telemetry transmitters to use one set of lead wires.

Additional cables, adapters and lead wires are also available to standardize most other manufacturer monitors and telemetry transmitters, allowing your hospital to use one set of lead wires per patient throughout the continuum of care.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Combiner Cables for Standardization

Other manufacturers



700-0008-30	3 Lead	Datascope	305 cm/10 ft	1 Cable
700-0008-31	5 Lead	Datascope	305 cm/10 ft	1 Cable
700-0008-33	3 Lead	Datex-Ohmeda	305 cm/10 ft	1 Cable
700-0008-34	5 Lead	Datex-Ohmeda	305 cm/10 ft	1 Cable
700-0008-35	6 lead	Datex-Ohmeda	305 cm/10 ft	1 Cable
700-0008-36	3 Lead	GE	305 cm/10 ft	1 Cable
700-0008-37	5 Lead	GE	305 cm/10 ft	1 Cable
700-0008-38	10 Lead (5+5)	GE	305 cm/10 ft	1 Cable
700-0008-39	6 Lead	GE	305 cm/10 ft	1 Cable
700-0008-40	10 Lead (6+4)	GE	305 cm/10 ft	1 Cable
700-0008-47	3 Lead	Philips	305 cm/10 ft	1 Cable

The ECG Combiner System can standardize Spacelabs monitors and Ultraview & Aria Telemetry transmitters to use one set of lead wires.

Additional cables, adapters and lead wires are also available to standardize most other manufacturer monitors and telemetry transmitters, allowing your hospital to use one set of lead wires per patient throughout the continuum of care.



Reusable



Disposable



Multi-Pinch



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Combiner Cables for Standardization



Other manufacturers					
700-0008-48	10 Lead (5+5)	Philips	305 cm/10 ft	1 Cable	
700-0008-49	5 Lead	Philips	305 cm/10 ft	1 Cable	
700-0008-50	10 Lead (6+4)	Philips	305 cm/10 ft	1 Cable	
700-0008-51	6 Lead	Philips	305 cm/10 ft	1 Cable	
700-0008-41	3 Lead	Medtronic LifePak	305 cm/10 ft	1 Cable	
700-0008-42	5 Lead	Medtronic LifePak	305 cm/10 ft	1 Cable	
700-0008-43	6 lead	Medtronic LifePak	305 cm/10 ft	1 Cable	
700-0008-44	3 Lead	Nihon Kohden	305 cm/10 ft	1 Cable	
700-0008-45	5 Lead	Nihon Kohden	305 cm/10 ft	1 Cable	
700-0008-46	6 lead	Nihon Kohden	305 cm/10 ft	1 Cable	
700-0008-52	3 Lead	Zoll, 90° angle	305 cm/10 ft	1 Cable	
700-0008-53	5 Lead	Zoll, 90° angle	305 cm/10 ft	1 Cable	
700-0008-54	6 lead	Zoll, 90° angle	305 cm/10 ft	1 Cable	

The ECG Combiner System can standardize Spacelabs monitors and Ultraview & Aria Telemetry transmitters to use one set of lead wires.

Additional cables, adapters and lead wires are also available to standardize most other manufacturer monitors and telemetry transmitters, allowing your hospital to use one set of lead wires per patient throughout the continuum of care.



Reusable



Disposable



Multi-Pinch



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Combiner Adapters for Standardization

Other manufacturers

	700-0008-60	5 Lead	Datascope Telemetry	1 Adapter
	700-0008-61	3 Lead	Datascope	1 Adapter
	700-0008-62	5 Lead	Datascope	1 Adapter
	700-0008-63	3 Lead	Datex-Ohmeda	1 Adapter
	700-0008-64	5 Lead	Datex-Ohmeda	1 Adapter
	700-0008-65	5 Lead	GE Apex Pro FH Telemetry	1 Adapter
	700-0008-66	6 Lead	GE Apex Pro FH Telemetry	1 Adapter

The ECG Combiner System can standardize Spacelabs monitors and Ultraview & Aria Telemetry transmitters to use one set of lead wires.

Additional cables, adapters and lead wires are also available to standardize most other manufacturer monitors and telemetry transmitters, allowing your hospital to use one set of lead wires per patient throughout the continuum of care.



Reusable



Disposable



Multi-Pinch



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Combiner Adapters for Standardization

Other manufacturers

	700-0008-67	3 Lead	GE	1 Adapter
	700-0008-68	5 Lead	GE	1 Adapter
	700-0008-69	6 Lead	GE Monitors	1 Adapter
	700-0008-70	3 Lead	Nihon Kohden	1 Adapter
	700-0008-71	6 Lead	Nihon Kohden	1 Adapter
	700-0008-72	3/5 Lead	Philips Telemetry	1 Adapter
	700-0008-73	3 Lead	Philips	1 Adapter
	700-0008-74	5 Lead	Philips	1 Adapter
	700-0008-75	3 Lead	Siemens Drager	1 Adapter
	700-0008-76	5 Lead	Siemens Drager	1 Adapter

Shielded ECG lead wires and cables for use with Spacelabs Command Modules 90496 & 91496 and mCare Monitor 91220 (3- and 5- lead wire sets only).



Shielded

Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Lead Wire Sets, Shielded, Reusable, AAMI



Snap

Snap Sets

700-0007-00	3 Lead, Shielded	61 cm/24 in	1 Set
700-0007-02	3 Lead, Shielded	102 cm/40 in	1 Set
700-0007-08	5 Lead, Shielded	61 cm/24 in	1 Set
700-0007-10	5 Lead, Shielded	102 cm/40 in	1 Set
700-0007-16*	10 Lead, limb/chest, Shielded	61 cm/24 in	1 Set
700-0007-18*	10 Lead, limb/chest, Shielded	102 cm/40 in	1 Set

*10 Lead ECG Lead Wire Sets not compatible with mCare Monitor 91220



Pinch

Pinch Sets

700-0006-00	3 Lead, Shielded	61 cm/24 in	1 Set
700-0006-02	3 Lead, Shielded	102 cm/40 in	1 Set
700-0006-08	5 Lead, Shielded	61 cm/24 in	1 Set
700-0006-10	5 Lead, Shielded	102 cm/40 in	1 Set
700-0006-16	10 Lead, limb/chest, Shielded	61 cm/24 in	1 Set
700-0006-18	10 Lead, limb/chest, Shielded	102 cm/40 in	1 Set

*10 Lead ECG Lead Wire Sets not compatible with mCare Monitor 91220

Shielded ECG lead wires and cables for use with Spacelabs Command Modules 90496 & 91496 and mCare Monitor 91220 (3- and 5- lead wire sets only).



Shielded

Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Lead Wires, Individual Colors, Shielded, Reusable, AAMI



Snap

24" Snap Individual Color

012-0498-00	Red	61 cm/24 in	3/Pack
012-0498-04	Black	61 cm/24 in	3/Pack
012-0498-34	Black	61 cm/24 in	5/Pack
012-0498-05	Brown	61 cm/24 in	3/Pack
012-0498-03	White	61 cm/24 in	3/Pack
012-0498-01	Green	61 cm/24 in	3/Pack
012-0498-11	Yellow/Brown	61 cm/24 in	3/Pack
012-0498-12	Green/Brown	61 cm/24 in	3/Pack
012-0498-13	Blue/Brown	61 cm/24 in	3/Pack
012-0498-14	Orange/Brown	61 cm/24 in	3/Pack
012-0498-15	Violet/Brown	61 cm/24 in	3/Pack



Snap

40" Snap Individual Color

012-0498-06	Red	102 cm/40 in	3/Pack
012-0498-09	Black	102 cm/40 in	3/Pack
012-0498-39	Black	102 cm/40 in	5/Pack
012-0498-27	Brown	102 cm/40 in	3/Pack
012-0498-08	White	102 cm/40 in	3/Pack
012-0498-07	Green	102 cm/40 in	3/Pack

Shielded ECG lead wires and cables for use with Spacelabs Command Modules 90496 & 91496 and mCare Monitor 91220 (3- and 5- lead wire sets only).



Shielded

Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Cables, Shielded, AAMI



With 17 PIN Adapter

700-0008-08	3 Lead, Shielded	305 cm/10 ft	1 Cable
700-0008-06	5 Lead, Shielded	305 cm/10 ft	1 Cable
700-0008-00	10 Lead, Shielded	305 cm/10 ft	1 Cable

ECG Accessories



Electrodes

369062-001	ECG electrode, tape, single use	Adult	100/Box
690-0262-00	ECG electrode, transparent tape, general/stress	Adult	30/Box

Shielded ECG lead wires and cables for use with Spacelabs Command Modules 90496 & 91496 and mCare Monitor 91220 (3- and 5- lead wire sets only).



Shielded

Order Code	Description	Length	Qty
------------	-------------	--------	-----

Neonatal ECG Supplies & Accessories



DIN Unshielded, Neonatal ECG Cable

700-0008-12	3 Lead		1 Cable
-------------	--------	--	---------



DIN Unshielded, Neonatal Electrodes


685-0037-00	Pre-wired neonatal electrodes, single use Pre-gelled. Silver/Silver Chloride	3-Lead	100 Sets/Box
-------------	---	--------	--------------



Electrode

369063-001	ECG neonatal electrode, single use, Individually packaged. Pre-gelled. Silver/Silver Chloride.	Ped/Neo	100/Box
------------	--	---------	---------

DIN Unshielded ECG lead wires and cables for use with Spacelabs Multiparameter Module 90470 ; available in 3- and 5-lead configurations.



DIN Unshielded

Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Lead Wire Sets, DIN Unshielded, Reusable, AAMI



Snap


Snap Sets			
012-0283-00	3 Lead	46 cm/18 in	1 Set
012-0283-11	3 Lead	61 cm/24 in	1 Set
012-0283-02	3 Lead	102 cm/40 in	1 Set
012-0285-00	5 Lead	46 cm/18 in	1 Set
012-0285-11	5 Lead	61 cm/24 in	1 Set
012-0285-02	5 Lead	102 cm/40 in	1 Set



Pinch

Pinch Sets			
012-0284-00	3 Lead	46 cm/18 in	1 Set
012-0284-11	3 Lead	61 cm/24 in	1 Set
012-0284-02	3 Lead	102 cm/40 in	1 Set
012-0286-00	5 Lead	46 cm/18 in	1 Set
012-0286-11	5 Lead	61 cm/24 in	1 Set
012-0286-02	5 Lead	102 cm/40 in	1 Set

DIN Unshielded ECG lead wires and cables for use with Spacelabs Multiparameter Module 90470 ; available in 3- and 5-lead configurations.



DIN Unshielded

Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Lead Wires, Individual Colors, DIN Unshielded, Reusable, AAMI




Snap

Snap Individual Color

012-0295-00	Red	46 cm/18 in	3/Pack
012-0295-10	Black	46 cm/18 in	3/Pack
012-0295-20	Brown	46 cm/18 in	3/Pack
012-0295-30	White	46 cm/18 in	3/Pack
012-0295-40	Green	46 cm/18 in	3/Pack
012-0295-06	Red	61 cm/24 in	3/Pack
012-0295-26	Black	61 cm/24 in	3/Pack
012-0295-46	Brown	61 cm/24 in	3/Pack
012-0295-16	White	61 cm/24 in	3/Pack
012-0295-36	Green	61 cm/24 in	3/Pack
012-0295-02	Red	102 cm/40 in	3/Pack
012-0295-22	Black	102 cm/40 in	3/Pack
012-0295-42	Brown	102 cm/40 in	3/Pack
012-0295-12	White	102 cm/40 in	3/Pack
012-0295-32	Green	102 cm/40 in	3/Pack

DIN Unshielded ECG lead wires and cables for use with Spacelabs Multiparameter Module 90470 ; available in 3- and 5- lead configurations.



DIN Unshielded

Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Lead Wires, Individual Colors, DIN Unshielded, Reusable, AAMI



Pinch

Pinch Individual Color

012-0296-00	Red	46 cm/18 in	3/Pack
012-0296-20	Black	46 cm/18 in	3/Pack
012-0286-00	Brown	46 cm/18 in	3/Pack
012-0296-10	White	46 cm/18 in	3/Pack
012-0296-04	Red	61 cm/24 in	3/Pack
012-0296-24	Black	61 cm/24 in	3/Pack
012-0296-44	Brown	61 cm/24 in	3/Pack
012-0296-14	White	61 cm/24 in	3/Pack
012-0296-34	Green	61 cm/24 in	3/Pack
012-0296-02	Red	102 cm/40 in	3/Pack
012-0296-22	Black	102 cm/40 in	3/Pack
012-0296-42	Brown	102 cm/40 in	3/Pack
012-0296-12	White	102 cm/40 in	3/Pack
012-0296-32	Green	102 cm/40 in	3/Pack

DIN Unshielded ECG lead wires and cables for use with Spacelabs Multiparameter Module 90470 ; available in 3- and 5-lead configurations.



DIN Unshielded

Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Adapter and Cable, AAMI



Required for use with Multiparameter Module 90470

103-0059-00	17-pin to 6-pin adapter, allows the ECG cables designed for Command Modules 90496/91496 to work on the Multiparameter Module 90470		1 Each
-------------	--	--	--------

