

SUPPLIES & ACCESSORIES

INVASIVE PRESSURE



INVASIVE PRESSURE CABLES



For use with Spacelabs Command Modules (90496 & 91496) configured with the invasive pressure option. Option B supports two pressures; Option C supports four pressures.



(9X496-B) (9X496-C)
Command Modules

Order Code	Description	Size	Qty
------------	-------------	------	-----

Integrated Adapter Cables (connects directly to transducer)

Edward Lifesciences



700-0293-01	Single pressure integrated cable	396 cm/13 ft	1 Each
-------------	----------------------------------	--------------	--------



700-0292-01	Dual pressure integrated cable	396 cm/13 ft	1 Each
-------------	--------------------------------	--------------	--------

Abbott Laboratories/Hospira/ICU Medical



700-0295-01	Single pressure integrated cable	396 cm/13 ft	1 Each
-------------	----------------------------------	--------------	--------



700-0294-01	Dual pressure integrated cable	396 cm/13 ft	1 Each
-------------	--------------------------------	--------------	--------

Interface Cable & Adapters

Interface Cable



690-0021-00	Compatible with Datex-Ohmeda TC-VTK or Becton Dickinson model 684085 pressure transducers. Requires adapter cable 700-0091-00 or 700-0028-00	366 cm/12 ft	1 Each
-------------	--	--------------	--------

Interface Adapter- Transducer interface cable required from transducer company




700-0791-00	Single pressure cable	122 cm/4 ft	1 Each
-------------	-----------------------	-------------	--------



700-0028-00	Dual pressure cable	122 cm/4 ft	1 Each
-------------	---------------------	-------------	--------

Diaphragm Dome



690-0018-00	Diaphragm dome, single use Rotating luer lock fittings Compatible with Datex-Ohmeda TA101OND		12/Box
-------------	--	--	--------

