

At Statcorp Medical, we care about the quality of patient care that you provide. That's why for over 20 years, we have manufactured the highest quality Unifusors, to help you achieve the best possible clinical outcomes.

Convenient, durable, leak-proof, certified accurate

Designed for rapid IV infusion, invasive pressure monitoring, and autotransfusion procedures; the Unifusor line is the clear choice in pressure infusors.

Features

- Latex free
- Made in the USA
- 100% QA tested during manufacturing process

Multi-purpose

- 3 models:
 - Unifusor
 - Unifusor Classic
 - Unifusor Plus
- 3 sizes:
 - 500cc
 - 1000cc
 - 3000cc





Unifusor features

- IV pole loop
- Fluid bag hook
- Durable hospital white nylon cuff
- Reinforced translucent sleeve
 - Can read through the sleeve
 - Type of fluid
 - Fluid levels
- No-crimp tubing
- Palm-fitted handbulb
- Available with either a traditional thumbwheel-type air control or a leak-proof stopcock valve
- 100% leak tested during manufacturing

Unifusor Classic features

- IV pole loop
- Fluid bag hook
- Durable hospital white nylon cuff
- Reinforced translucent sleeve
 - Can read through the sleeve
 - Type of fluid
 - Fluid levels
- Fluid levels
- No-crimp tubing
- Palm-fitted handbulb
- Available with either a traditional thumbwheel-type air control or a leak-proof stopcock valve
- Piston gauge with:
 - 360° gauge visibility
 - Bleed valve prevents over-inflation
- 100% leak tested during manufacturing

Unifusor Plus Reusable features

- IV pole loop
- Fluid bag hook
- Ultra-durable blue nylon cuff
- Reinforced translucent sleeve
 - Can read through the sleeve
 - Type of fluid
 - Fluid levels
- No-crimp tubing
- Palm-fitted handbulb / gauge assembly
 - Easy-to-use trigger-type airway valve
 - Large, easy-to-read gauge
- 100% leak tested during manufacturing

Product Numbers			
	500cc	1000cc	3000cc
Unifusor with aneroid gauge & thumbwheel valve	803FGA 24/cs	804FGA 24/cs	805FGA 10/cs
Unifusor with aneroid gauge & stopcock valve	803SGA 24/cs	804SGA 24/cs	805SGA 10/cs
Unifusor Classic with piston gauge & thumbwheel valve	903TGA 24/cs	904TGA 24/cs	905TGA 10/cs
Unifusor Classic with piston gauge & stopcock valve	903SGA 24/cs	904SGA 24/cs	905SGA 10/cs
Unifusor Plus with trigger / gauge assembly	1103 1 ea 1103-05 5/bx	1104 1 ea 1104-05 5/bx	1105 1 ea 1105-05 5/bx
Unifusor Plus Bag only	1103X 1 ea 1103X-05 5/bx	1104X 1 ea 1104X-05 5/bx	1105X 1 ea 1105X-05 5/bx
Trigger / Gauge assembly only	1199, 1 ea	1199-05, 5/bx	

Aneroid gauge



Stopcock valve



Piston gauge



Thumbwheel valve



Trigger gauge



Product is  Marked



ISO 13485
 Certificate Number FM 40458

Statcorp Medical
 an OSI Optoelectronics company
 14476 Duval Place West, Suite 303
 Jacksonville, FL 32218
 800-992-0014, 904-786-5113
 www.statcorpmedical.com

Dealer Imprint