

SPACELABS SOFT CUFF DISPOSABLE BLOOD PRESSURE CUFFS

Spacelabs brand cuffs provide a comprehensive solution for your NIBP needs. The Soft Cuff Disposable Blood Pressure Cuffs are available in both regular and long sizes.

- Soft Cuff fabric
- Long size cuffs for large adults
- Long cuffs feature additional length allowing ¼ extra turn for greater security and elimination of “pop-offs”



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

Spacelabs Soft Cuff Disposable BP Cuffs

Single Tube

For use with Command Modules 90496/91496, DM3 Monitor 91330 and mCare Monitor 91220



Bayonet (HP)

714-5020-00	Single Tube, HP	Child	12 - 19 cm	5/Pack
714-5021-00	Single Tube, HP	Small Adult	17 - 25 cm	5/Pack
714-5023-00	Single Tube, HP	Adult	23 - 33 cm	5/Pack
714-5032-00	Single Tube, HP	Adult Long	23 - 33 cm	5/Pack
714-5025-00	Single Tube, HP	Large Adult	31 - 40 cm	5/Pack
714-5033-00	Single Tube, HP	Large Adult Long	31 - 40 cm	5/Pack
714-5027-00	Single Tube, HP	Thigh	35 - 44 cm	5/Pack

Single Tube Hose



714-0018-00	90496/91496/91330	Adult	275 cm/9 ft	1 Each
714-0060-00	90496/91496/91330	Adult	366 cm/12 ft	1 Each
714-0061-00	91220	Adult	275 cm/9 ft	1 Each

Double Tube

For use with Command Modules 90496/91496



Bayonet (HP)

714-3020-00	Double Tube, HP	Child	12 - 19 cm	5/Pack
714-3021-00	Double Tube, HP	Small Adult	17 - 25 cm	5/Pack
714-3023-00	Double Tube, HP	Adult	23 - 33 cm	5/Pack
714-3032-00	Double Tube, HP	Adult Long	23 - 33 cm	5/Pack
714-3025-00	Double Tube, HP	Large Adult	31 - 40 cm	5/Pack
714-3033-00	Double Tube, HP	Large Adult Long	31 - 40 cm	5/Pack
714-3027-00	Double Tube, HP	Thigh	38 - 50 cm	5/Pack

Double Tube to Single Tube Hose



714-0018-12	90496/91496	Adult	275 cm/9 ft	1 Each
-------------	-------------	-------	-------------	--------