

TRULINK® NYLON REUSABLE BLOOD PRESSURE CUFFS

TruLink Nylon Reusable Blood Pressure Cuffs are available in seven adult/pediatric sizes for accurate measurements.

- Meets AHA / AAMI standards
- Color coded for accurate sizing
- Manufactured without PVC/DEHP
- Provides a fluid barrier
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

TruLink Nylon Reusable BP Cuffs

Single Tube

For use with Command Modules 90496/91496, DM3 Monitor 91330 and mCare Monitor 91220



Bayonet (HP)

715-1121-10	Single tube, HP	Child	13 - 20 cm	1 Each
715-1131-10	Single tube, HP	Small Adult	18 - 26 cm	1 Each
715-1141-10	Single tube, HP	Adult	26 - 35 cm	1 Each
715-1151-10	Single tube, HP	Adult Long	29 - 38 cm	1 Each
715-1161-10	Single tube, HP	Large Adult	32 - 42 cm	1 Each
715-1171-10	Single tube, HP	Large Adult Long	35 - 44 cm	1 Each
715-1181-10	Single tube, HP	Thigh Cuff	42 - 50 cm	1 Each

Single Tube Hose



714-0018-00	90496/91496/91330	Adult	275 cm/9 ft	1 Each
714-0060-00	90496/91496/91330	Adult	366 cm/12 ft	1 Each
714-0061-00	91220	Adult	275 cm/9 ft	1 Each

Double Tube

For use with Command Modules 90496/91496



Bayonet (HP)

715-1221-10	Double tube, HP	Child	13 - 20 cm	1 Each
715-1231-10	Double tube, HP	Small Adult	18 - 26 cm	1 Each
715-1241-10	Double tube, HP	Adult	26 - 35 cm	1 Each
715-1251-10	Double tube, HP	Adult Long	29 - 38 cm	1 Each
715-1261-10	Double tube, HP	Large Adult	32 - 42 cm	1 Each
715-1271-10	Double tube, HP	Large Adult Long	35 - 44 cm	1 Each
715-1281-10	Double tube, HP	Thigh Cuff	42 - 50 cm	1 Each

Double Tube to Single Tube Hose



714-0018-12	90496/91496	Adult	275 cm/9 ft	1 Each
-------------	-------------	-------	-------------	--------