

## TRULINK® SOFT DISPOSABLE BLOOD PRESSURE CUFFS

TruLink Soft Disposable Blood Pressure Cuffs are available in seven adult/pediatric sizes for accurate measurements.

- Soft fabric helps prevent perspiration under cuff
- Meets AHA / AAMI standards
- Color coded for accurate sizing
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

### TruLink Soft Disposable BP Cuffs

#### Single Tube

*For use with Command Modules 90496/91496, DM3 Monitor 91330 and mCare Monitor 91220*



Bayonet (HP)

715-1123-10	Single tube, HP	Child	13 - 20 cm	5/Box
715-1133-10	Single tube, HP	Small Adult	18 - 26 cm	5/Box
715-1143-10	Single tube, HP	Adult	26 - 35 cm	5/Box
715-1153-10	Single tube, HP	Adult Long	29 - 38 cm	5/Box
715-1163-10	Single tube, HP	Large Adult	32 - 42 cm	5/Box
715-1173-10	Single tube, HP	Large Adult Long	35 - 44 cm	5/Box
715-1183-10	Single tube, HP	Thigh Cuff	42 - 50 cm	5/Box

#### Single Tube Hose



714-0018-00	90496/91496/91330	Adult	275 cm/9 ft	1 Each
714-0060-00	90496/91496/91330	Adult	366 cm/12 ft	1 Each
714-0061-00	91220	Adult	275 cm/9 ft	1 Each

#### Double Tube

*For use with Command Modules 90496/91496*



Bayonet (HP)

715-1223-10	Double tube, HP	Child	13 - 20 cm	5/Box
715-1233-10	Double tube, HP	Small Adult	18 - 26 cm	5/Box
715-1243-10	Double tube, HP	Adult	26 - 35 cm	5/Box
715-1253-10	Double tube, HP	Adult Long	29 - 38 cm	5/Box
715-1263-10	Double tube, HP	Large Adult	32 - 42 cm	5/Box
715-1273-10	Double tube, HP	Large Adult Long	35 - 44 cm	5/Box
715-1283-10	Double tube, HP	Thigh Cuff	42 - 50 cm	5/Box

#### Double Tube to Single Tube Hose



714-0018-12	90496/91496	Adult	275 cm/9 ft	1 Each
-------------	-------------	-------	-------------	--------